



With the highest capacity of Health o meter® Professional's remote display scales, the 752KL is a great supplement to accommodate ambulatory bariatric patients. The remote display and portable base make it convenient to move the scale anywhere it is needed.

Advantages

- ✓ **EMR Capable**
Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions
- ✓ **Remote Display**
The remote display provides versatility with the ability to be mounted on a wall or placed on a table
- ✓ **600 lb Capacity Floor Scale**
The 752KL is the market's highest capacity floor scale, providing more versatile use
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off, Parent/Child Weighing Function) included at no additional cost
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI



Remote display can be mounted on a wall or placed on a table



Optional carrying case (order # 64771)



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
14 1/4" x 14 1/4" x 2 5/8"
362 mm x 362 mm x 67 mm
- Product Footprint (w x h x d):
14 1/2" x 14 1/2" x 2 3/4"
368 mm x 368 mm x 70 mm
- Product Weight: 14 lb / 6 kg
- Display: 1" / 25 mm LCD
- Connectivity: USB, optional Pelstar® wireless technology
- Power Source: 100-240V adapter (ADPT50) included, 6 C-cell batteries (not included)
- Optional Carrying Case: 64771
- Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function
- 2 year limited warranty

Model	Description
752KG	Heavy duty remote display digital scale, kilograms only
SS-752KL	ScaleSurance 2 year extended warranty for 752 models*
ADPT50	100-240V adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America
64771	Optional carrying case for 752KL & 752KG
PELSTARONE-P	Wireless thermal printer