


THE 752KL

HEAVY DUTY REMOTE DISPLAY DIGITAL SCALE

The 752KL is a remote display digital scale with a high capacity that supports up to 660 pounds and provides the optimum in quality, reliability, and ease of use. The remote display provides versatility with the ability to be mounted on a wall or placed on a table. The portable base and optional carrying case make it convenient to move the scale anywhere it is needed. Designed to meet the needs of doctor's offices, clinics, and ambulatory surgery centers, the scale's large non-skid low profile platform requires a minimal step-up providing a safe weighing environment.

The scale's large LCD screen provides a user-friendly interface that includes BMI, hold/release, zero, LB/KG, tare, and more. The 752KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising it's not surprising that Health o meter® Professional Scales is the industry leading name in medical scales.

120+ LOCAL
FIELD SPECIALISTS

100%
CONNECTIVITY

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®

P r o f e s s i o n a l

Health o meter® Professional 752KL

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

660 lb / 300 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via USB, optional Pelstar® wireless technology

Display

1" / 25 mm LCD Screen

Power Source

100-240V adapter (ADPT50) included, 6 C-cell batteries not included

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Hold/release
Auto zero
Auto off
Parent/child weighing function

Platform Size (w x d x h)

14 1/4" x 14 1/4" x 2 5/8"
(362 mm x 362 mm x 67 mm)

Product Footprint - Base (w x d x h)

14 1/2" x 14 1/2" x 2 3/4"
(368 mm x 368 mm x 70 mm)

Product Weight

14 lb / 6 kg

Portability

64771 - optional carrying case

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
752KL	Heavy duty remote display digital scale, LB & KG
752KG	Heavy duty remote display digital scale, KG only
SS-752KL	ScaleSurance® 2 year extended warranty for 752 series*
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America
64771	Optional Carrying Case for 752KL & 752KG

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscales.com/ConnectivitySolutions



64771
Optional Carrying Case

To learn more, call **1-800-815-6615**
or visit **homscales.com/752KL**



*Not available in all countries.

Health o meter®
Professional