



The 844KL offers a high capacity paired with the highest resolution in the bath scale line. A good option for patients who need greater precision in their weight measurement.

Advantages

✓ **440 lb / 200 kg Capacity Floor Scale**

The high capacity of this floor scale provides more versatile use

✓ **1 1/4" LCD Display**

The large display shows the user an easy-to-read weight measurement

✓ **EMR Connectivity**

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions

✓ **Low Profile Platform**

Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience

Pack of 2



844KL Display



Specifications

- Capacity: 440 lb / 200 kg / 31 st
- Resolution: 0.05 kg / 0.1 lb / 1/8 lb
- Platform Size (w x d x h):
12 5/8" x 12 5/8" x 1 7/8" /
321 mm x 321 mm x 48 mm
- Product Footprint (w x d x h):
12 5/8" x 12 5/8" x 1 7/8" /
321 mm x 321 mm x 48 mm
- Product Weight: 4 lb / 2 kg
- Display: 1 1/4" / 32 mm LCD
- Power Source: (2) 3-Volt Lithium batteries included
- Connectivity: USB
- Optional Carrying Case item # 64771
- Functions: LB / KG / ST Switch, Auto Zero, 5 Second Weight Hold, Auto Off
- 1 Year Limited Warranty

Model	Description
844KLS	Digital platform scale, pounds & kilograms
64771	Digital platform scale, kilograms only



Carrying Case for 844KL
Item# 64771