

THE PTS-1000KL

PATIENT TRANSFER SCALE

The PTS-1000KL is a transfer board with a built-in medical scale that provides a weight measurement when transferring a patient. All departments can benefit from the Patient Transfer Scale, particularly Stroke Units, the Emergency Room, ICU, and Radiology.

Designed by Clinicians for Clinicians the Patient Transfer Scale weighs critical or immobile patients instantly and accurately helping shorten door-to-needle times. It offers a superior alternative to other weighing equipment, such as in-bed scales, that are often not calibrated or inaccurate.

The scale is light and portable and can be utilized as part of an existing patient transfer process when a weight is needed, and will not disrupt the patient flow process.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that Health o meter® Professional Scales is the industry-leading name in medical scales.

120+ LOCAL
FIELD SPECIALISTS

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE



PTS-1000KL

Health o meter®
Professional

Health o meter® Professional PTS-1000KL

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

550 lb / 250 kg

Resolution

1 lb / 0.5 kg

Display

3" x 1" / 76 mm x 25 mm Backlight LCD

Power Source

3.7V rechargeable battery

Functions

LB/KG conversion

Zero

Hold/release

Auto off

Weighing Mat Size (w x d)

71" x 21"

(1805mm x 533mm)

Product Footprint (w x d x h)

71" x 27 1/2" x 1"

(1805mm x 700mm x 30mm)

Product Weight

25 lb / 11.4 kg

Measuring Tape

Height guide shown on mat 100cm – 180cm
shown at 5 cm


Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
PTS-1000KL	Patient Transfer Scale, LB & KG
SS-PTS-1000KL*	ScaleSurance® 2 year extended warranty for Patient Transfer Scale

Always use a Slide Sheet
with the Patient Transfer Scale!



Check out the Single Patient Use
and Reuseable Slide Sheets.

To learn more, call **1-800-815-6615**
or visit **homscapes.com/PTS-1000KL**



*Not available in all countries.



Health o meter®
Professional