Optional Height Rod for 2101KL & 2101KG

Easy to Use

Smooth sliding headpiece and telescoping rod allows the user to quickly and easily read the height measurement.

Durable Construction

Aluminum height rod withstands frequent use. Rigorous laboratory tests showed that the Health o meter[®] Professional STROD lasts 2 ½ times longer than similar models.

Specifications

- Measuring Range: 2"- 90" / 6 cm 230 cm
- Graduation: 1/16" / 1mm
- 2 Year Warranty

Large Measurement Range The 2"- 90" / 6 cm - 230 cm measuring range provides height readings for a wide range of patients.



Manufactured by: Pelstar, LLC. 9500 West 55th St. McCook, IL 60525-7110 USA Phone: 1-800-815-6615 Fax: 1-800-654-7330

 σ

O

S

S

Ф

Ω

STROD_F&B_20200211