



Health o meter[®]

P r o f e s s i o n a l

**WEIGH
EASIER[®]**

Product Catalog

why choose us?

Health o meter®

P r o f e s s i o n a l

weigh**easier**®

What sets Health o meter® Professional Scales apart is that everything we do - from product design to technical support - is done with our customers in mind. Each decision and each response is formed by first asking *What will be weigh**easier** for our customer?* – whether the customer is a physician, healthcare system, distributor, or sales team.

#1 Medical Scale Brand in North America

Health o meter® has been a trusted brand since 1919. Today, Health o meter® Professional Scales is the U.S. market leader and the fastest growing brand in Medical Scales.* We offer:

- High quality products that serve the entire continuum of medical care
- The industry's largest network of trained and qualified field representatives
- Direct toll-free number for Technical Support with customer advocates based in the USA
- Fast shipments and industry-leading fill rates
- Standardization programs
- No-cost product evaluations

100% Focus on Healthcare

Because all we do is support the healthcare market, we understand the impact of healthcare reform, Value-Based Purchasing, and other industry-specific initiatives. We understand that providers are overwhelmed with patient care and are focused on operating efficiently and with an eye on costs. We concern ourselves with patient safety, supply chain efficiency, and technology standardization because our customers are concerned about those things. At Health o meter® Professional Scales, we design our products to help patients and the providers and institutions that care for them.

Products Designed in the USA

Health o meter® Professional's US based engineers and product specialists design high-value, high-quality devices for users around the globe.

Rated #1 by Medical Distributors & Healthcare Professionals!

Health o meter® Professional products are manufactured, designed and owned by Pelstar LLC. Health o meter® is a registered trademark of Sunbeam, Inc. used under license. ©2021 Pelstar®, LLC – PRODUCT CAT-ONLINE 2021 0121

*As reported by Global Industry Analysts, Inc., Patient Scales market survey and Global Healthcare Exchange, (GHX) sales data. All Health o meter® Professional products are registered with the US Food & Drug Administration.

table of contents

02

ScaleSurance® Extended Warranty

03

Connectivity Solutions

04-51

Products

- ▲ Stand-On
- ▲ Body Composition Scale
- ▲ Pediatric
- ▲ Wheelchair / Chair / Multi-Use
- ▲ Specialty
- ▲ Professional Home Care
- ▲ Height Measuring
- ▲ Accessories

52-66

Product Specifications

67

Glossary

Symbols

68

Global Support

ISO 9001:2015
ISO 13485:2016
CERTIFIED



Extended Warranty Program

Health o meter® Professional Scales offers an extended warranty that protects a facility's scale investment. This program extends the standard Limited Warranty that is available with select Health o meter® Professional products.

Benefits of an Extended Warranty:

- Peace of mind against the financial hardship of unexpected equipment replacement costs
- Quick product replacement – including free outbound shipping
- Hassle-free coverage provided by the manufacturer, not a third party warranty company
- Immediate, live interaction with product specialists for troubleshooting, diagnostics, or EMR connectivity assistance

A fraction of the cost of buying a new scale!

Additional Details:

- Can be purchased with a new scale or for a facility's existing scale before the current warranty has expired.
- Available on Health o meter® Professional digital scales and digital height rods that have a 2 year limited warranty and G6 Body Comp Scales

For more information visit
www.homscales.com/scalesurance

Not available in all countries

To purchase an extended warranty or for more information, contact your medical supply distributor today!

Connectivity Solutions

Health o meter® Professional Scales offers connectivity solutions that transmit weight, height, and Body Mass Index (BMI) measurements to vital signs devices, PCs and EMR systems. Automatic data transfer reduces transcription errors and is the fastest, easiest way to capture patient information. Device connectivity is an important tool to improve a healthcare facility's workflow efficiency, patient safety and EMR compliance efforts.

Our scales are designed to interface with any EMR, and our list of connectivity partners is always growing. We will work with any EMR company at the request of our customers to facilitate direct interfaces with their systems. Among the systems we are compatible with are:



The following solutions are available for connecting Health o meter® Professional digital scales to a vital signs device, PC or EMR system.

PELSTARONE

External wireless solution to EMR capable Health o meter® Professional digital scale. There is no need to purchase a new scale. Allows virtually all Health o meter® Professional digital scales to interface directly with Welch Allyn Connex® Vital Signs Monitors and Spot Vital Signs® LXi monitors, Allscripts TouchWorks® and Professional™ systems, and Midmark® IQ Manager® software.

L-HOMWA -1

Wired connection for Welch Allyn® Spot Vital Signs® LXi Monitor. (Adapter and USB cable connect the scale to the Welch Allyn® Spot Vital Signs® LXi Monitor.)

C-HOMWA-1

Wired connection for Welch Allyn® Connex Devices. (Adapter and USB cable connect the scale to Welch Allyn® Connex® Devices.)



Products bearing this symbol can connect to EMR systems via these Connectivity Solutions.



Stand-On



★ TOP SELLER
500KL

EYE-LEVEL DIGITAL SCALE

The 500KL is the #1 selling physician scale in North America. With a built-in height rod, the scale enhances workflow by measuring weight and height and calculating Body Mass Index without moving the patient. Optional handlebars provide stability to keep both patients and caregivers safe.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Built-in Height Rod: The easy-to-use height rod provides height measurements from 30" - 84 ¼" / 76 cm - 214 cm.

Calculates BMI: Scale functions include the ability to calculate the patient's BMI.



Optional 500HB Handlebar Accessory can be mounted to your new or existing 500KL, 500KG, 500KL-BT, 500KG-BT, 501KL and 501KG.

Capacity	Resolution	Platform Size
----------	------------	---------------

550 lb / 250 kg	0.2 lb / 0.1 kg	13 ¾" (w) x 16 ½" (d) x 2 ⅜" (h) 349 mm (w) x 419 mm (d) x 60 mm (h)
-----------------	-----------------	---

Height Rod	Height Rod Resolution	Display
------------	-----------------------	---------

30" - 84 ¼" 76 cm - 214 cm	⅛" / 1 mm	1" / 25 mm LCD
-------------------------------	-----------	----------------

Options & Accessories:

Model	Description
500KG	Scale in kilograms only
500KLHB	Scale with handlebars included
500KLAD	500KL with adapter included
500KGAD	500KG with adapter included
501KL	TAA Compliant eye-level scale
501KG	TAA Compliant eye-level scale, KG only
500KL-BT	500KL with built-in Pelstar® wireless technology
500KG-BT	500KG with built-in Pelstar® wireless technology
500HB	Optional handlebar accessory for 500 series scales
SS-500KL	2 year extended warranty
ADPT31	Adapter for 500 series scales
PELSTARONE-P	Wireless thermal printer

Power Source: 6 AA batteries included, 100-240V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Auto Zero, Auto Off, Reweigh, Recall, Everlock®, Parent / Child Weighing Function

Wheels: 2 **Connectivity:** USB, optional Pelstar® wireless technology

499KL

WAIST-HIGH DIGITAL SCALE



The 499KL is a compact version of our most popular scale, the 500KL. The waist-high display is easy to read and the scale is an attractive, robust option for use in medical offices. With an optional height rod, the 499KL is an economical choice for almost any environment.

Capacity: 500 lb / 220 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

13 3/4" (w) x 16 1/2" (d) x 2 3/8" (h)

349 mm (w) x 419 mm (d) x 60 mm (h)

Display: 1" / 25 mm LCD

Connectivity: USB, optional Pelstar® wireless technology

Wheels: 2

Power Source: 6 AA batteries included, 120V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Options & Accessories:

Model	Description
499KG	Scale in kilograms only
499KLHR	499KL with height rod included
499KGHR	499KG with height rod included
499KLAD	499KL with adapter included
499KGAD	499KG with adapter included
499KLHRAD	499KL with height rod and adapter included
499KGHRAD	499KG with height rod and adapter included
SS-499KL	2 year extended warranty
499KLROD	Height rod for 499 series scales
ADPT31	Adapter for 499 series scales
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer



599KL

HEAVY DUTY WAIST-HIGH DIGITAL SCALE



Great for doctor offices and clinics with minimal floor space, the 599KL has a high capacity with waist-high display that can swivel 350°. For providers wishing to measure height, an optional height rod is available.

Capacity: 600 lb / 272 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

14 1/4" (w) x 14 1/4" (d) x 2 5/8" (h)

362 mm (w) x 362 mm (d) x 67 mm (h)

Display: 1" / 25 mm LCD, 350° Swivel

Connectivity: USB, optional Pelstar® wireless technology

Wheels: 2

Power Source: 100-240V adapter included, 6 C batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Options & Accessories:

Model	Description
599KG	Scale in kilograms only
599KLHR	599KL with height rod included
599KGHR	599KG with height rod included
SS-599KL	2 year extended warranty
499KLROD	Height rod for 599 series scales
ADPT50	Adapter included, item # for replacement or spare
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer





★ TOP SELLER

597KL

HEAVY DUTY EYE-LEVEL DIGITAL SCALE

With a high capacity and mechanical height rod, this scale properly weighs and measures the heaviest ambulatory patients. The 350° swivel head makes it easy for clinicians to view measurements and interact with the scale.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Built-in Height Rod: The easy-to-use height rod provides height measurements 23 5/8" - 84" / 60 cm - 213 cm.

Calculates BMI: Scale functions include the ability to calculate the patient's BMI.



Capacity	Resolution	Platform Size
----------	------------	---------------

600 lb / 272 kg	0.2 lb / 0.1 kg	14 1/4" (w) x 14 1/4" (d) x 2 5/8" (h) 362 mm (w) x 362 mm (d) x 67 mm (h)
-----------------	-----------------	---

Height Rod	Height Rod Graduation	Display
------------	-----------------------	---------

23 5/8" - 84" 60 cm - 213 cm	1/8" / 1 mm	1" / 25 mm LCD 350° Swivel
---------------------------------	-------------	-------------------------------

Options & Accessories:

Model	Description
597KG	Scale in kilograms only
SS-597KL	2 year extended warranty
ADPT50	Adapter included, item # for replacement or spare
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer

Power Source: 100-240V adapter included, 6 C batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out , BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Wheels: 2 **Connectivity:** USB, optional Pelstar® wireless technology

500KLF





★ TOP SELLER

600KL & ELEVATE

HEAVY DUTY EYE-LEVEL DIGITAL SCALE WITH DIGITAL HEIGHT ROD

Improve workflow with this top-of-the-line physician office scale. The built-in digital height rod eliminates the need to manually enter height for computing BMI. The 350° swivel head allows for line-of-sight viewing and patient privacy.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Automatically Calculates Weight, Height and BMI: Improve efficiency and accuracy of vital sign measurements.

Built-in Digital Height Rod: The easy-to-use, state of the art digital height rod provides measurements from 23 5/8" - 88" / 60 cm - 222 cm.

Patented Design: The 600KL's patented design provides a user-friendly and intuitive interface for measuring patient weight, height and BMI with a convenient rotating display head.

ELEVATE models include a connectivity kit that allows the scale to interface directly with Welch Allyn® vital signs monitors.



Capacity	Resolution	Platform Size
600 lb / 272 kg	0.2 lb / 0.1 kg	14 1/4" (w) x 14 1/4" (d) x 2 5/8" (h) 362 mm (w) x 362 mm (d) x 67 mm (h)
Height Rod	Height Rod Resolution	Display
23 5/8" - 88" 60 cm - 222 cm	0.1" / 0.1 cm	1" / 25 mm LCD 350° Swivel

Options & Accessories:

Model	Description
ELEVATE-C	Scale with Connex connectivity kit
ELEVATE-KG	ELEVATE-C in kilograms only
600KL-BT	600KL with built-in Pelstar® wireless technology & HOMWA connectivity kit
600KG-BT	600KG with built-in Pelstar® wireless technology & HOMWA connectivity kit
ELEVATE-BT	ELEVATE-C with built-in Pelstar® wireless technology & HOMWA connectivity kit
ELEVATE-KG-BT	ELEVATE-KG with built-in Pelstar® wireless technology & HOMWA connectivity kit
ADPT50	Adapter included, item # for replacement or spare
SS-600KL	2 year extended warranty
SS-ELEVATE-C	2 year extended warranty
SS-600KL-BT	2 year extended warranty
SS-ELEVATE-BT	2 year extended warranty
PELSTARONE-P	Wireless thermal printer

Power Source: 100-240V adapter included, 6 C batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out , BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Wheels: 2 **Connectivity:** USB, optional Pelstar® wireless technology

498KL

REMOTE DISPLAY DIGITAL SCALE



The 498KL is a versatile choice for almost any clinical environment. Its portability and remote display head allow the scale to be moved from one location to another.

Capacity: 500 lb / 220 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

13 3/4" (w) x 16 1/2" (d) x 2 5/8" (h)

349 mm (w) x 419 mm (d) x 60 mm (h)

Display: 1" / 25 mm LCD

Connectivity: USB, optional Pelstar® wireless technology

Power Source: 6 AA batteries included, 120V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Options & Accessories:

Model	Description
498KLAD	498KL with adapter included
SS-498KL	2 year extended warranty
ADPT31	Adapter for 498KL
64771	Carrying case for 498KL
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer



Optional Height Rods, see page 46

752KL

HEAVY DUTY REMOTE DISPLAY DIGITAL SCALE



With the highest capacity of Health o meter® Professional's remote display scales, the 752KL is a great supplement to accommodate ambulatory bariatric patients. The remote display and portable base make it convenient to move the scale anywhere it is needed.

Capacity: 600 lb / 272 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

14 1/4" (w) x 14 1/4" (d) x 2 5/8" (h)

362 mm (w) x 362 mm (d) x 67 mm (h)

Display: 1" / 25 mm LCD

Connectivity: USB, optional Pelstar® wireless technology

Power Source: 100-240V adapter included, 6 C-cell batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Options & Accessories:

Model	Description
752KG	Scale in kilograms only
SS-752KL	2 year extended warranty
ADPT50	Adapter included, item # for replacement or spare
64771	Carrying case for 752KL
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer



Optional Height Rods, see page 46



★ TOP SELLER

349KLX

DIGITAL FLOOR SCALE

The 349KLX is a top selling scale that can be used in a host of locations, from MD offices to clinics to home care. Users can place the remote display head where it is easiest for them to read the weight. The scale's small footprint makes it a great choice for weight loss clubs and home use, too.

Interfaces with Welch Allyn® Spot Vital Signs® LXi Monitor:

Users can connect the scale to display weight measurement on the Welch Allyn® monitor.

Low Profile Platform: Low platform height offers patient a comfortable and stable weighing experience.

Remote Display: The remote display provides versatility with the ability to be mounted on a wall or placed on a table.

Capacity	Resolution	Platform Size
----------	------------	---------------

400 lb / 180 kg	0.2 lb / 0.1 kg	12 ½" (w) x 12" (d) x 1 ⅞" (h) 318 mm (w) x 305 mm (d) x 48 mm (h)
-----------------	-----------------	---

Display

1" / 25 mm LCD

Options & Accessories:

Model	Description
349KLXAD	Scale with adapter included
SS-349KLX	2 year extended warranty
ADPT40	Adapter for 349KLX
64771	Carrying case for 349KLX

Optional Height Rods, see page 46

Power Source: 6 AA batteries included, 120V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, Zero, Tare, Auto Zero, Auto Off, Parent / Child Weighing Function

Connectivity: Serial Port

3001KLE-AM



★ TOP SELLER

3001KL-AM & 3001KL-AMX



ANTIMICROBIAL DIGITAL PLATFORM SCALE

Incredibly robust, these platform scales feature feet and wheels specially designed to withstand frequent use and regular movement. Antimicrobial coating on the handrails, transport handle, and keypad help reduce the spread of harmful microbes. The 3001KL-AMX offers extended handrails that provide extra support and help patients maintain balance when stepping onto the platform.

Infection Prevention: Silver ionic Antimicrobial Powder Coating suppresses the growth of bacteria and helps protect against cross-contamination.

Fall Prevention: “Live” handrails are advanced technology that allow the patient to hold on for stability while the scale accurately calculates weight.

Low Profile Platform: Requires only minimal step up, offering the patient a comfortable and stable weighing experience.

Durability: Convenient transport handle allows the user to easily move the scale to where it is needed in the facility. Rugged wheels allow for smooth and easy portability.



3001KL-AMX
ANTIMICROBIAL DIGITAL
PLATFORM SCALE WITH
EXTENDED HANDRAILS



Capacity	Resolution	Platform Size
----------	------------	---------------

1,000 lb / 454 kg	0.2 lb / 0.1 kg	22" (w) x 15 3/4" (d) x 2 5/8" (h) 559 mm (w) x 400 mm (d) x 67 mm (h)
-------------------	-----------------	---

Display

1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen

Options & Accessories:

Model	Description
3001KG-AM	3001KL-AM, assembled scale in kilograms only
3001KL-AMUA	3001KL-AM, ships unassembled
3001KG-AMUA	3001KG-AM, ships unassembled
3001KL-AMHR	3001KL-AM with height rod
3001KG-AMHR	3001KG-AM with height rod
3001KL-AMX	Assembled scale with extended handrails
3001KG-AMX	3001KL-AMX in kilograms only
SS-3001KL-AM	2 year extended warranty
SS-3001KL-AMX	2 year extended warranty
3001HR	Height rod for 3001KL-AM & 3001KG-AM
ADPT30	Adapter included, item # for replacement or spare
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer

Power Source: 100-240V adapter ADPT30 (included), 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Parent / Child Weighing Function, Everlock®

Wheels: 2 **Connectivity:** USB, optional Pelstar® wireless technology

1100KL

DIGITAL PLATFORM SCALE



A comprehensive platform scale that supports up to 1000 pounds offers “live” handrails and a low profile platform to provide a stable and comfortable weighing experience for patients of any size. The keypad features Microban® antimicrobial protection. Ideal for any type of facility including primary care, rehab or hospital.

Capacity: 1,000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

22" (w) x 15 3/4" (d) x 2 1/2" (h)

559 mm (w) x 400 mm (d) x 64 mm (h)

Display: 1 1/2" / 38 mm High-Contrast

Color TFT -LCD Screen

Connectivity: USB, optional Pelstar® wireless technology

Options & Accessories:

Wheels: 2

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Parent / Child Weighing Function, Everlock®

Model	Description
1100KG	Scale in kilograms only
1100KL-BT	1100KL with built-in Pelstar® wireless technology
1100KG-BT	1100KG with built-in Pelstar® wireless technology
1100KLHR	1100KL with height rod included
1100KGHR	1100KG with height rod included
1100KL-EHR	1100KL with digital height rod included
SS-1100KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
PROPLUSROD	Height rod for 1100 series of scales
PELSTARONE-P	Wireless thermal printer



2101KL

DIGITAL PLATFORM SCALE



This large bariatric scale features “live”, long handrails to help frail and balance-challenged people get onto the scale and hold on while taking an accurate weight. The keypad features Microban® antimicrobial protection. Suitable for use in MD offices, long term care, hospitals or anywhere there is adequate space for the device.

Capacity: 1,000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

26" (w) x 22" (d) x 2 3/8" (h)

660 mm (w) x 559 mm (d) x 60 mm (h)

Display: 1 1/2" / 38 mm High-Contrast

Color TFT-LCD Screen

Connectivity: USB, optional Pelstar® wireless technology

Options & Accessories:

Wheels: 2

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Parent / Child Weighing Function, Everlock®

Model	Description
2101KG	Scale in kilograms only
2101KL-BT	2101KL with built-in Pelstar® wireless technology
2101KG-BT	2101KG with built-in Pelstar® wireless technology
2101KLHR	2101KL with height rod included
2101KGHR	2101KG with height rod included
SS-2101KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
STROD	Height rod for 2101KL
PELSTARONE-P	Wireless thermal printer





★ TOP SELLER

402KL

MECHANICAL BEAM SCALE

The classic physician scale is built with a robust steel frame and powder coated column to prevent rust. Includes a convenient height rod so clinicians can measure height and weight all in one place.

Height Rod Included: The easy-to-use height rod provides height measurements from 23 5/8" - 84" / 60 cm - 213 cm. Headpiece can be folded down when not in use.

Exclusive 10 Year Limited Warranty: Health o meter® Professional's industry-leading warranty protects against defects of materials or workmanship.

Capacity	Graduation	Platform Size
390 lb / 180 kg (Counterweights raise to 490 lb / 210 kg)	1/4 lb / 100 g	10 1/2" (w) x 14" (d) x 3 1/4" (h) 267 mm (w) x 356 mm (d) x 83 mm (h)
Height Rod	Height Rod Graduation	
23 5/8" - 84" 60 cm - 213 cm	1/8" x 1 mm 60 cm - 213 cm	

Options & Accessories:

Model	Description
402KLCW	402KL with counterweights included
402KLWH	402KL with wheels included
402KLWHCW	402KL with wheels and counterweights included
402KG50	402KL with 50 gram graduation (special order)
402LB	Scale in pounds only (400 lb capacity)
402LBWH	Scale in LB only with wheels included
402KG	Scale in kilograms only (special order)
201HR-400	Optional telescopic metal height rod for 400 series
55090	Counterweights for 402KL
55000	Wheel brackets for 402KL, 402KG & 402LB

Functions: Measures weight in pounds and kilograms. Measures height in inches and centimeters.

400KL

MECHANICAL BEAM SCALE



A traditional mechanical scale with an industry-leading 10 year warranty. Features a unique coating inside the column to prevent rust when used in gyms, locker rooms and humid environments.

Capacity: 390 lb / 180 kg
(Counterweights raise capacity to 490 lb / 210 kg)

Graduation: ¼ lb / 100 g

Platform Size:

10 ½" (w) x 14" (d) x 3 ¼" (h)

267 mm (w) x 356 mm (d) x 83 mm (h)

Options & Accessories:

Functions: Measures weight in pounds and kilograms.

Model	Description
400KLCW	Scale with counterweights included
400KLWH	Scale with wheels included
400KLWHCW	Scale with wheels and counterweights included
201HR-400	Optional telescopic height rod for 400KL
402KLROD	Optional metal height rod for 400KL
55090	Counterweights for 400KL
55000	Wheel brackets for 400KL

450KL

MECHANICAL BEAM SCALE WITH ROTATING POISE BARS



The 450KL is the top of the line mechanical scale and a market leader in doctor offices for decades. Special features include long-lasting engraved numerals on rotating (LB or KG) poise bars and included counterweights that increase capacity to 500 lb / 200 kg.

Capacity: 500 lb / 200 kg (with included counter weights)

Graduation: ¼ lb / 100 g

Height Rod:

23 ⅝" - 84" / 60 cm - 213 cm

Height Rod Graduation: ⅛" x 1 mm

Platform Size:

10 ½" (w) x 14" (d) x 3 ¼" (h)

Functions: Measures weight in pounds and kilograms. Measures height in inches and centimeters.

Model	Description
450KLWH	Scale with wheels included
201HR-400	Optional telescopic height rod for 400KL
55000	Wheel brackets for 450KL

3102KL-AM

HEAVY-DUTY ANTIMICROBIAL PLATFORM SCALE



“Live” handrails allow the patient to hold onto handrails for stability while the scale accurately calculates weight. Low profile platform requires only minimal step up, offering the patient a comfortable and stable weighing experience. Heavy-duty construction increases stability for unsteady patients and reduces the chance of tipping over. Antimicrobial keypad, handrails and transport handle help reduce patient and staff exposure to microbes.

Capacity: 1000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

19 ¼" (w) x 21 ¾" (d) x 1 ¾" (h)

489 mm (w) x 552 mm (d) x 44 mm (h)

Display: ¾" / 19 mm LCD

Connectivity: USB, optional Pelstar® wireless technology

Options and Accessories:

Model	Description
SS-3102KL-AM	2 year extended warranty
ADPT50	Adapter included, item # for replacement or spare
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer

Wheels: 2 rugged, scuff-proof

Power Source: 6 D batteries and 100-240V adapter included

Functions: LB / KG Conversion, LB / KG Lock Out, Hold, Recall, Reweigh, Zero, Auto Zero, Auto Off, Audible Signal, Low Battery Indicator



Don't forget that our wheelchair scales can also be used as stand-on scales. See pages 31-34



Stand-On Scales



Body Composition



BCS-G6

BODY COMPOSITION SCALE

Affordable, accurate and attractive. The perfect body composition scale for practices desiring consistent results and accuracy regardless of when it's used. With no pre-test restrictions and proprietary Ohms Modulation Logic technology, G6 scales eliminate the Body Fat Bounce Effect and can be used with patients at any time of day.

15 Second Measurement Time: On-screen display delivers rapid access to weight, body fat %, body water %, lean mass, BMR & BMI. A full illustrated printout can be generated in under a minute.

Large Platform: The spacious platform allows patients to find a comfortable and secure stance.

Few Pretest Restrictions: Analysis can be performed at any time of day and under almost any circumstances, resulting in the greatest flexibility of any body composition analyzer.

High Capacity: The scale features an industry-leading 1,000 lb capacity that provides the flexibility and accuracy to weigh any patient.



Capacity	Resolution	Platform Size
1,000 lb / 453.5 kg	0.1 lb / 0.1 kg	14" (w) x 16" (d) x 2 7/8" (h) 356 mm x 406 mm x 48 mm

Display

6-window display (1/2" digits) including weight, % body fat, % body water, % muscle mass, BMI, REE

Optional printer, sold separately, not available through Health o meter® Professional Scales

Power Source: 110-220V power adapter (included)

Functions: 15 Second Measurement Time, Body Fat, Body Water, Fat Free Mass, BMI, Basal Metabolic Rate, Visceral Fat, Measures weight in pounds and kilograms

More About the G6:

The G6 body composition analyzer is a valuable tool that can help caregivers address obesity and the serious medical conditions that often occur in overweight and obese populations. The detailed body composition analysis provided by the device and the Patient Report Templates that it generates are excellent tools for educating patients about their individual risks and also can be used as reinforcement for patients embarking on a weight loss or health improvement regimen. Engaged patients can help practices improve outcomes and support efforts to achieve results consistent with MACRA and other regulatory, quality and reimbursement initiatives.

What sets the G6 apart from other body composition scales is the proprietary Ohms Modulation Logic that eliminates the need for extensive and intrusive pre-test restrictions, while providing accurate and repeatable results. Unlike other scales, the G6 doesn't require patients to prepare for the test by not eating, drinking or exercising prior to use. This makes the G6 a product you can use at any time with patients...and no down-time means that you can be more productive.

The **easiest** to use, most **cost-effective** body composition analyzer available...

OK To Test



The G6 is the only body composition analyzer that eliminates the “body fat bounce effect” – the frustrating outcome that can show increased body fat despite lost weight. Our proprietary Ohms Modulation Logic (OML) algorithm takes the guesswork out of body composition analysis. OML makes it possible to get accurate results without subjecting the patient or medical practice to a long list of unrealistic pre-test restrictions. The G6 – with OML – is the right choice for practices that want to get the most value from their body composition analyzer.

G6 Continued

Users can choose from 4 graphic patient report templates including Adult / Bariatric, Adult Spanish, Pediatric and Fitness versions. These personalized printouts include detailed information about a patient's body composition and educational content regarding the health risks associated with being overweight. Providing patients with a personalized risk assessment is an effective way to engage them in their care and motivate them to take control of their health.



For a current list of compatible printers visit www.homscales.com/innovations/body-comp-scale



Waiting room brochures (BCS-BROCHURE) are effective marketing tools for practices that charge patients for the G6 analysis. The flyers come in a handy stand and provide a way to inform patients that a Body Composition Analysis is available and can be performed during their visit. The brochure highlights the measurements provided by the analysis and what they can learn from the results.

Model	Description
BCS-G6-ADULT	The whole-body adult only scale is ideal for practices that exclusively see patients 18 and over
BCS-G6-DUO	The whole-body adult/pediatric version is a good choice for providers who see adults and children
BCS-G6-LIMBS	The segmented version provides adults with body composition information about their torso and each of their limbs
BCS-G61-UPPER	The upper-body-only scale provides a partial body composition analysis for adults who do not want to remove shoes and socks
BCS-G6-ADULT-SP	The whole body, adult-only scale is ideal for practices that exclusively see adults and who want a Spanish language report
SS-BCS-G6-ADULT	2 year warranty for all versions of the G6 scale
IPO-ADU	Patient Report Template - Adult (5-page report, 25 reports/package)
IPO-PED	Patient Report Template - Pediatric (7-page report, 25 reports/package)
IPO-ATH	Patient Report Template - Athletic (3-page report, 50 reports/package)
IPO-ADU-SP	Spanish Language Patient Report Template - Adult (5-page report, 25 reports/package)
BCS-BROCHURE	Waiting room patient brochure (50/package)



Pediatric

★ TOP SELLER

553KL



553KL
with 553HR



553KL
with 553EHR



553KL with
2210CART

DIGITAL PEDIATRIC TRAY SCALE

The 553KL is Health o meter® Professional's most popular infant scale and perfect for newborns, infants and toddlers. The scale operates on batteries or AC power so you can use it anywhere. Practices benefit from one-button operation that makes it easy to use while ensuring that the baby remains safely on the scale.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Measuring Tape Included: Easy-to-read measuring tape allows user to measure both weight and height on the same device.

Motion Detection: Scale weighs accurately despite movement of the child on the scale, and locks onto the measurement.



Capacity	Resolution	Tray Size
44 lb / 20 kg	0.5 oz / 10 g	24 7/8" (w) x 13" (d) x 1 3/8" (h) 632 mm (w) x 330 mm (d) x 35 mm (h)
Measuring Tape	Display	
0" - 24" 0 cm - 58 cm	1" / 25 mm LCD	

Options & Accessories:

Model	Description
553KLCT	Scale with rolling cart included
553KG	Scale in kilograms only
553KL-HR	553KL with attachable mechanical height rod
553KG-HR	553KG with attachable mechanical height rod
553KL-EHR	553KL with attachable digital height rod
553KG-EHR	553KG with attachable digital height rod
SS-553KL	2 year extended warranty
553HR	Attachable mechanical height rod for 553 series scales
553EHR	Attachable digital height rod for 553 series scales
2210CART	Rolling cart for 553KL & 553KG
553CASE	Carrying case for 553KL & 553KG
ADPT30	Adapter for 553KL & 553KG
PELSTARONE	External wireless adapter
PELSTARONE-P	Wireless thermal printer

Power Source: 6 AA batteries included, 100-240V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lockout, Zero, Tare, Hold / Release, Auto Off

Connectivity: USB, optional Pelstar® wireless technology

522KL

DIGITAL PEDIATRIC TRAY SCALE



This easy-to-use scale is perfect for busy pediatric offices or clinics where infants are seen regularly. The tray has a convenient measuring tape to determine the baby's length during the weighing session.

Capacity: 50 lb / 23 kg

Resolution: 0 - 20 lb / 0.2 oz; 20 - 50 lb / 0.5 oz;
0 - 9 kg / 5 g; 9 - 23 kg / 10 g

Tray Size:

24 1/8" (w) x 14 1/2" (d) x 2 5/8" (h)

613 mm (w) x 368 mm (d) x 67 mm (h)

Measuring Tape: 0" - 23" / 0 cm - 58 cm

Display: 1" / 25 mm LCD

Connectivity: USB, optional Pelstar® wireless technology

Power Source: 6 AA batteries included, 120V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, Hold / Release, Zero, Tare, Auto Zero, Auto Off, Measures length in inches and centimeters, Motion Detection

Options & Accessories:

Model	Description
522KG	Scale in kilograms only
522KL-HR	522KL with attachable mechanical height rod
522KG-HR	522KG with attachable mechanical height rod
522KL-EHR	522KL with attachable digital height rod
522KG-EHR	522KG with attachable digital height rod
SS-522KL	2 year extended warranty
522HR	Attachable mechanical height rod for 522 series scales
522EHR	Attachable digital height rod for 522 series scales
ADPT31	Adapter for 522 series scales
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer



524KL

DIGITAL PEDIATRIC SEAT SCALE



The deep infant seat on the 524KL reduces the movement of young children to speed up the weighing process. Weight is captured easily while the seat supports the child's head and body.

Capacity: 50 lb / 23 kg

Resolution: 0 - 20 lb / 0.2 oz; 20 - 50 lb / 0.5 oz;
0 - 9 kg / 5 g; 9 - 23 kg / 10 g

Seat Size:

24 1/4" (w) x 14 1/4" (d) x 18 1/2" (h)

616 mm (w) x 362 mm (d) x 470 mm (h)

Display: 1" / 25 mm LCD

Connectivity: USB, optional Pelstar® wireless technology

Power Source: 6 AA batteries included, 120V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Motion Detection

Options & Accessories:

Model	Description
524KG	Scale in kilograms only
SS-524KL	2 year extended warranty
ADPT31	Adapter for 524KL & 524KG
PELSTARONE	External wireless kit
PELSTARONE-P	Wireless thermal printer



549KL

DIGITAL PEDIATRIC TRAY SCALE



Lightweight and battery-powered with basic functionality, this scale is a good choice for practices with occasional infant patients or for visiting nurses.

Capacity: 44 lb / 20 kg

Resolution: 0.5 oz / 10 g

Tray Size:

19 ¼" (w) x 13" (d) x 2 ⅝" (h)

489 mm (w) x 330 mm (d) x 67 mm (h)

Display: ⅝" / 15 mm LCD

Options & Accessories:

Power Source: 9V battery included

Functions: LB / KG Conversion, LB / KG Lock Out, Zero, Tare, Auto Zero, Auto Off

Model	Description
549KG	Scale in kilograms only
64769	Carrying case for 549KL
SS-549KL	2 year extended warranty

550KL

INFANT TO TODDLER, 2 IN 1 DIGITAL SCALE



For practices looking to weigh both infants and toddlers, this dual-purpose scale can be used with the infant weighing tray or toddler step-on platform. An easy-to-use scale that is lightweight and portable with simple functionality and a digital display.

Capacity: 44 lb / 20 kg

Resolution: 0.5 oz / 10 g

Tray Size:

23 ½" (w) x 8" (d) x 2" (h)

597 mm (w) x 203 mm (d) x 51 mm (h)

Options & Accessories:

Platform Size:

12 ½" (w) x 8 ⅝" (d) x 1 ¼" (h)

318 mm (w) x 219 mm (d) x 32 mm (h)

Display: ½" / 13 mm digit height LCD

Power Source: 9V battery included

Functions: Auto Off, Hold, Tare, Unit

Model	Description
SS-550KL	2 year extended warranty
553CASE	Carrying case for 550KL
2210CART	Rolling cart for 550KL
eBABYHR	Digital height rod for 550KL (purchase with 550EHRBRACKET to mount on scale)
550EHRBRACKET	550KL mounting bracket for eBABYHR (purchase with eBABYHR to mount on scale)

8320KL

DIGITAL PORTABLE PEDIATRIC TRAY SCALE



8320KL Closed

Excellent choice for visiting nurses who need a lightweight pediatric scale. The folding scale is easy to use and compact to carry in its optional bag.

Capacity: 33 lb / 15 kg

Resolution: 0.5 oz / 10 g

Tray Size:

23 ½" (w) x 11" (d) x 7/8" (h)

597 mm (w) x 279 mm (d) x 22 mm (h)

Display: ¾" / 19 mm

Options & Accessories:

Power Source: 2 AAA batteries included

Functions: LB / KG Conversion, Zero, Tare, Auto Zero, Auto Off

Model	Description
SS-8320KL	2 year extended warranty
64771	Carrying case for 8320KL

2210KL

HIGH RESOLUTION DIGITAL NEONATAL / PEDIATRIC TRAY SCALE



The 2210KL neonatal scale was designed for use in the NICU, delivery room, nursery and pediatric practices where a high degree of accuracy is critical to ensure that babies are gaining weight. A convenient cart is available to easily move the scale from location to location.

Capacity: 45 lb / 20 kg / 20,000 g

Resolution: 1 g (to 10,000 g);
5 g (10,001 to 20,000 g); 0.1 oz (0 - 45 lb)

Tray Size:
24" (w) x 15 3/8" (d) x 3" (h)
610 mm (w) x 391 mm (d) x 76 mm (h)

Display: 7/8" / 22 mm LCD

Options & Accessories:

Power Source:

6 C batteries included

Functions: LB / GM Conversion,
LB / KG Conversion,
LB / KG / GM Lock Out, Zero, Tare,
Reweigh / Recall, Auto Zero, Auto Off,
Motion Detection

Model	Description
2210KLCT	Scale with rolling cart included
SS-2210KL	2 year extended warranty
2210CART	Rolling cart for 2210KL



2210KL4

HIGH RESOLUTION DIGITAL NEONATAL / PEDIATRIC SCALE WITH 4-SIDED TRAY



Weighing newborns, including premature infants, requires a critical degree of accuracy, and the 2210KL4 offers the highest resolution in the industry. The 4-sided weighing tray, which offers an extra degree of assurance that the baby is protected, can be removed for cleaning and features an integral measuring tape to monitor the infant's growth.

Capacity: 45 lb / 20 kg / 20,000 g

Resolution: 1 g (to 10,000 g);
5 g (10,001 to 20,000 g); 0.1 oz (0-45 lb)

Tray Size:
27" (w) x 15 3/8" (d) x 3" (h)
686 mm (w) x 391 mm (d) x 76 mm (h)

Display: 7/8" / 22 mm LCD

Options & Accessories:

Power Source: 6 C batteries included

Functions: LB / GM Conversion,
LB / KG Conversion,
LB / KG / GM Lock Out, Zero,
Tare, Reweigh / Recall, Auto Zero,
Auto Off, Motion Detection

Model	Description
2210KL4CT	Scale with rolling cart included
SS-2210KL4	2 year extended warranty
2210CART	Rolling cart for 2210KL4



1522KL

MECHANICAL PEDIATRIC BEAM SCALE



This high capacity pediatric tray scale is perfect for small infants or young children who can sit independently. The mechanical poise bar does not require electricity, allowing the scale to be used anywhere.

Capacity: 130 lb / 65 kg

Graduation: 1 oz / 20 g

Tray Size:

24 1/8" (w) x 14 1/2" (d) x 3 1/4" (h)

613 mm (w) x 368 mm (d) x 83 mm (h)

Measuring Tape: 0" - 23" / 0 cm - 58 cm

Functions: Weighs in pounds and kilograms, measures length in inches and centimeters.

1524KL

MECHANICAL PEDIATRIC SEAT SCALE



The built-in chair on the 1524KL is convenient for use with infants who are unable to sit; it can also help small children keep still for weight measurements. The mechanical poise bar does not require electricity, allowing the scale to be used anywhere.

Capacity: 130 lb / 65 kg

Graduation: 1 oz / 20 g

Seat Size:

24 1/4" (w) x 14 1/4" (d) x 18 1/2" (h)

622 mm (w) x 362 mm (d) x 470 mm (h)

Functions: Measures weight in pounds and kilograms.

386S-01 / 386KGS-01

MECHANICAL TRAY SCALE
POUNDS ONLY / KILOGRAMS ONLY



These economical mechanical tray scales are simple to use and include a built-in measuring tape. The sturdy body and convenient handle make it a great choice for visiting nurses.

Capacity:

386S-01: 50 lb

386KGS-01: 25 kg

Graduation:

386S-01: ¼ lb

386KGS-01: 100 g

Tray Size: 19 ¾" (w) x 12 ¾" (d) x 3 ¾" (h)

492 mm (w) x 314 mm (d) x 95 mm (h)

Dial Size: 3" / 76 mm

Measuring Tape: 0" - 18" / 0 cm - 46 cm

Functions:

386S-01: Measures weight in pounds.

386KGS-01: Measures weight in kilograms.

BABYSLINGSCALE

HANGING PEDIATRIC SCALE



This scale is an economical, portable option for use with young children that are able to hold their heads up. The large dial makes it easy to read.

Capacity: 25 kg

Functions: Measures weight in kilograms.

Graduation: 100 g

Dial Size: 5 ½" / 140 mm

Options and Accessories:

Model	Description
BABYSLING	Additional sling / trouser

3400KL / 3401KL

DIGITAL WET DIAPER / LAP SPONGE / ORGAN SCALE



The solid, industrial design of the 3400KL includes a stable, stainless steel base and removable, stainless steel weighing pan (3401KL comes with a flat tray) that can be autoclaved for sterilization. Powered by a long-life rechargeable battery and backlit for any lighting level, this portable scale can be used in any environment, from OR to ED to NICU and beyond.

Capacity: 176 oz / 5,000 g

Resolution: 0.05 oz / 1 g

Base Size:

11 ½" (w) x 9 ⅞" (d) x 5 ¼" (h)

292 mm (w) x 250 mm (d) x 133 mm (h)

Pan Size (3400KL/3400KG):

10" (w) x 9 ⅞" (d) x 2 ½" (h)

254 mm (w) x 250 mm (d) x 62 mm (h)

Tray Size (3401KL):

7 ¾" (w) x 9 ¾" (d) x ⅝" (h)

197 mm (w) x 248 mm (d) x 16 mm (h)

Display: ⅝" / 16 mm LCD

Options & Accessories:

Power Source: 6V rechargeable battery included, 100 - 240V AC battery charger included

Functions: Zero, Tare, Auto Zero, Auto Off, g / oz Conversion, Low Battery Indicator, Backlit Display

Model	Description
3400KG	Scale with pan, kilograms only
3400KG-C	Scale with pan, kilograms only, Canada only
SS-3400KL	2 year extended warranty
SS-3400KG	2 year extended warranty
SS-3400KG-C	2 year extended warranty
SS-3401KL	2 year extended warranty
3400PAN	3 pack of stainless steel pans
3400TRAY	3 pack of stainless steel trays
3400ADPT	Battery charger for 3400KL, 3400KG, 3400KG-C & 3401KL



YG500R

MECHANICAL DIAPER / SMALL PLATFORM SCALE



Perfect for facilities that occasionally need to weigh small items, this cost-effective scale is a great choice.

Capacity: 500 g

Graduation: 2 g

Platform Size:

6 ½" (w) x 4 ¾" (d)

152 mm (w) x 121 mm (d)

Dial Size: 6" / 152 mm

Functions: Measures weight in grams.



**Wheelchair /
Chair / Multi-Use**

★ TOP SELLER

2600KL / 2610KL



DIGITAL WHEELCHAIR RAMP SCALE /
DIGITAL WHEELCHAIR DUAL RAMP SCALE

Our most popular wheelchair scale is also the most versatile, and is a good choice for any facility. The platform folds for convenient storage and the casters (the only scale in the market that has them!) make transport a breeze. For easy on/off access the 2610KL offers dual ramps.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Casters and Folding Ramps for Easy Portability and Storage:

4 wheels, including 2 casters, afford ease of movement and the folding ramps make it simple to store.

Wide Platform for Wheelchair or Bariatric Use:

Wide platform accommodates wheelchairs of all sizes and provides a comfortable weighing experience for bariatric patients.

Calculates BMI: Scale functions include the ability to calculate patient's BMI.



Capacity	Resolution	Platform Size
----------	------------	---------------

1,000 lb / 454 kg	0.2 lb / 0.1 kg	36" (w) x 32 1/4" (d) x 2 1/2" (h) 914 mm (w) x 819 mm (d) x 64 mm (h)
-------------------	-----------------	---

Display

1 1/2" / 38 mm High-Contrast Color TFT - LCD Screen

Wheels: 4,

including 2 casters

Connectivity:

USB, optional Pelstar® wireless technology

Options & Accessories:

Model	Description
2600KG	Scale in kilograms only
2610KG	Scale in kilograms only
2600KL-BT	2600KL with built-in Pelstar® wireless technology
2610KL-BT	2610KL with built-in Pelstar® wireless technology
2600KG-BT	2600KG with built-in Pelstar® wireless technology
2610KG-BT	2610KG with built-in Pelstar® wireless technology
SS-2600KL	2 year extended warranty
SS-2610KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
B2600RAMP	Secondary ramp for 2600KL & 2600KG
PELSTARONE-P	Wireless thermal printer

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock® Audible / Mute Sound Option

2500KL

DIGITAL WHEELCHAIR RAMP SCALE



This single-ramp wheelchair scale is perfect for use with standard size wheelchairs. To weigh patients in a seated position, a chair accessory is also available, item # CHAIRACC.

Capacity: 1,000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

29 3/8" (w) x 28 1/4" (d) x 2 3/8" (h)

746 mm (w) x 718 mm (d) x 60 mm (h)

Display: 1 1/2" / 38 mm High-Contrast Color TFT - LCD Screen

Connectivity: USB, optional Pelstar® wireless technology

Wheels: 2

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Everlock®

Options and Accessories:

Model	Description
2500KG	Scale in kilograms only
2500KL-BT	2500KL with built-in Pelstar® wireless technology
2500KG-BT	2500KG with built-in Pelstar® wireless technology
SS-2500KL	2 year extended warranty
CHAIRACC	Optional fold away seat
ADPT30	Adapter included, item # for replacement or spare
B2500RAMP	Secondary ramp for 2500KL and 2500KG
PELSTARONE-P	Wireless thermal printer



2500CKL

DIGITAL WHEELCHAIR RAMP SCALE WITH FOLD AWAY SEAT



Facilities looking for a scale that can be used with wheelchairs and takes advantage of a built-in chair, should consider the 2500CKL, which features both a ramp and a folding patient seat.

Platform Capacity: 1000 lb / 454 kg

Seat Capacity: 419 lb / 190 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

29 3/8" (w) x 28 1/4" (d) x 2 3/8" (h)

746 mm (w) x 718 mm (d) x 60 mm (h)

Seat Size: 18" (w) x 15" (d) x 18 1/2" (h)

457 mm (w) x 381 mm (d) x 470 mm (h)

Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen

Connectivity: USB, optional Pelstar® wireless technology

Wheels: 2

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Everlock®

Options & Accessories:

Model	Description
2500CKG	Scale in kilograms only
SS-2500CKL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
B2500RAMP	Secondary ramp for 2500CKL
PELSTARONE-P	Wireless thermal printer



2650KL

DIGITAL WHEELCHAIR DUAL RAMP SCALE



Facilities frequently weighing patients in wheelchairs will appreciate the 2650KL, which offers a remote display head and a medium-sized platform and folds for convenient storage. The dual ramps make it easy to wheel patients on and off the platform. Like all wheelchair scales, the 2650KL can be used for walk-on patients, as well.

Capacity: 1,000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

31 ½" (w) x 31 ½" (d) x 2" (h)

800 mm (w) x 800 mm (d) x 51 mm (h)

Display: ¾" / 19 mm LCD,

Remote Display

Connectivity: USB, optional Pelstar® wireless technology

Options & Accessories:

Wheels: 2

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option

Model	Description
2650KG	Scale in kilograms only
2650KL-BT	2650KL with built-in Pelstar® wireless technology
2650KG-BT	2650KG with built-in Pelstar® wireless technology
SS-2650KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
PELSTARONE-P	Wireless thermal printer



2700KL

DIGITAL WHEELCHAIR DUAL RAMP SCALE WITH EXTRA LARGE PLATFORM



Featuring the largest platform in the wheelchair scale collection, the 2700KL can easily accommodate even the widest wheelchair. It is a good choice in facilities where the scale will remain in one location.

Capacity: 1,000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

43" (w) x 42" (d) x 2 ½" (h)

1092 mm (w) x 1067 mm (d) x 64 mm (h)

Display: ¾" / 19 mm LCD,

Remote Display

Connectivity: USB, optional Pelstar® wireless technology

Options & Accessories:

Wheels: 2

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option

Model	Description
2700KG	Scale in kilograms only
2700KL-BT	2700KL with built-in Pelstar® wireless technology
2700KG-BT	2700KG with built-in Pelstar® wireless technology
SS-2700KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
PELSTARONE-P	Wireless thermal printer



2400KL

DIGITAL PORTABLE WHEELCHAIR SCALE



The 2400KL is the perfect scale for locations where occasional wheelchair-bound patients need to be weighed. The scale folds easily and requires minimal room for storage. The scale includes a handle for portability.

Capacity: 800 lb / 363 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size: Expands and adjusts from 25 3/4" - 50 1/8" / 654 mm - 1273 mm

Display: 3/4" / 19 mm LCD, Remote Display

Connectivity: USB, optional Pelstar® wireless technology

Options & Accessories:

Wheels: 2 wheels, 2 casters

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option


Model	Description
SS-2400KL ADPT30 PELSTARONE-P	2 year extended warranty Adapter included, item # for replacement or spare Wireless thermal printer



Wheelchair scales / Multi-use scales

2842KL

LARGE PLATFORM DIGITAL SCALE



With its large platform and stainless steel base, the 2842KL can serve in multiple capacities. It can be used as a basic stand-on scale with an extra-large platform, and it can also be used in shipping/receiving or in laundry areas.

Capacity: 600 lb / 270 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size: 22 1/4" (w) x 42" (d) x 2 1/2" (h)

565 mm (w) x 1067 mm (d) x 64 mm (h)

Display: 7/8" / 22 mm LCD

Power Source: 120V adapter included, 6 C-cell batteries (not included)

Functions: LB / KG Conversion, Zero, Tare, Motion Detection, Reweigh

Options & Accessories:

Model	Description
SS-2842KL ADPT31	2 year extended warranty Adapter included, item # for replacement or spare



2595KL

DIGITAL CHAIR SCALE



This best-selling chair scale has a comfortable seat so patients can sit safely while being weighed. The moveable armrests and footrests allow for easy access on and off the scale. A great scale for use in Long-Term Care centers, rehab facilities, and senior living homes.

Capacity: 600 lb / 270 kg

Resolution: 0.2 lb / 0.1 kg

Seat Size:

18 ¼" (w) x 15 ⅝" (d) x 16 ½" (h)

464 mm (w) x 397 mm (d) x 419 mm (h)

Display: 1 ⅞" / 29 mm LCD

Connectivity: USB, optional Pelstar wireless technology

Options and Accessories:

Model	Description
SS-2595KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
PELSTARONE-P	Wireless thermal printer

Wheels: 2

Power Source: 100-240V adapter included, 6 C-cell batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, Zero, Tare, Auto Off



594KL

DIGITAL CHAIR SCALE



The 594KL has moveable arm and foot rests to allow for easy patient access on and off this economical chair scale. Two wheels make portability easy for bringing the scale to the patient's location.

Capacity: 600 lb / 270 kg

Resolution: 0.2 lb / 0.1 kg

Seat Size:

18" (w) x 18 ½" (d)"

457 mm (w) x 470 mm (d)

Display: 1" / 25 mm LCD

Options and Accessories:

Model	Description
SS-594KL	2 year extended warranty

Wheels: 4 casters

Power Source: 100-240V adapter included, 6 AA batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BSA, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off



445KL

MECHANICAL CHAIR SCALE



Nursing homes, rehab hospitals and facilities that treat elderly or unstable patients may find this basic chair scale an appropriate choice for weighing their patients.

Capacity: 440 lb / 200 kg

Graduation: ¼ lb / 100 g

Seat Size:

18 ½" (w) x 17 ¼" (d) x 15 ¾" (h)

470 mm (w) x 438 mm (d) x 400 mm (h)

Wheels: 2

Functions: Measures weight in pounds and kilograms.





Specialty

★ TOP SELLER

3400KL / 3401KL



3401KL

DIGITAL WET DIAPER / LAP SPONGE / ORGAN SCALE

The solid, industrial design of the 3400KL includes a stable, stainless steel base and removable, stainless-steel weighing pan (3401KL comes with a flat tray) that can be autoclaved for sterilization. Powered by a long-life rechargeable battery and featuring a backlit display for all lighting levels, this portable scale can be used in environments, from OR to ED to NICU and beyond.

Stainless Steel Base and Weighing Surface: The stainless steel design is easy to clean and ideal for sterile environments.

Rechargeable Battery: Long battery life due to a 6V rechargeable battery that supports continuous use for days.

Backlit Display: Easy to read from any angle and in any light condition.



Capacity	Resolution	Base Size
176 oz / 5,000 g	0.05 oz / 1 g	11 1/2" (w) x 9 7/8" (d) x 5 1/4" (h) 292 mm (w) x 250 mm (d) x 133 mm (h)
Pan Size	Display	
10" (w) x 9 7/8" (d) x 2 1/2" (h) 254 mm (w) x 250 mm (d) x 64 mm (h)	5/8" / 16 mm LCD	
Tray Size		
7 3/4" (w) x 9 3/4" (d) x 5/8" (h) 197 mm (w) x 248 mm (d) x 16 mm (h)		

Power Source: 6V rechargeable battery included, 100 - 240V AC battery charger included

Functions: Zero, Tare, Auto Zero, Auto Off, g / oz Conversion, Low Battery Indicator, Backlit Display

Options & Accessories:

Model	Description
3400KG	KG-only scale
3400KG-C	KG-only scale, Canada only
SS-3400KL	2 year extended warranty
SS-3400KG	2 year extended warranty
SS-3400KG-C	2 year extended warranty
SS-3401KL	2 year extended warranty
3400PAN	3 pack of stainless steel pans
3400TRAY	3 pack of stainless steel trays
3400ADPT	Battery charger for 3400KL, 3400KG, 3400KG-C & 3401KL

2000KL

DIGITAL IN-BED / STRETCHER SCALE



This scale can weigh immobile patients up to 400 lb / 181 kg. It has a durable steel structure, fabric sling, and wheels for portability.

Capacity: 400 lb / 181 kg

Resolution: 0.2 lb / 0.1 kg

Stretcher Size:

70 1/2" (w) x 30" (d)

1791 mm (w) x 762 mm (d)

Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen

Connectivity: USB, optional Pelstar® wireless technology

Options and Accessories:

Wheels: 4 (1 with safety lock)

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option, Hoyer lift

Model	Description
2000KG	Scale in kilograms only
2000KL-BT	2000KL with built-in Pelstar® wireless technology
2000KG-BT	2000KG with built-in Pelstar® wireless technology
SS-2000KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
PELSTARONE-P	Wireless thermal printer



YG500R

MECHANICAL DIAPER / SMALL PLATFORM SCALE



Perfect for facilities that occasionally need to weigh small items, this cost-effective scale is a great choice.

Capacity: 500 g

Graduation: 2 g

Platform Size:

6 1/2" (w) x 4 3/4" (d)

152 mm (w) x 121 mm (d)

Dial Size: 6" / 152 mm

Functions: Mechanically measures weight in grams.



Professional Home Care

800KL

DIGITAL FLOOR SCALE (PACK OF 4)



Our popular digital bath-style scale is great for use in the home with its simple design and its ability to match any décor. A great option for visiting nurses and health assessment professionals.

Capacity: 397 lb / 180 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

12 ¼" (w) x 12 ¼" (d) x 1 ⅛" (h)

311 mm (w) x 311 mm (d) x 29 mm (h)

Display: 1 ½" / 38 mm LCD

Options & Accessories:

Power Source: (2) 3V lithium battery included

Functions: LB / KG Switch, Auto Zero, 5 Second Weight Hold, Auto Off

Model	Description
800KLS	Digital floor scale (single)
64771	Carrying case for 800KL

880KL

HEAVY-DUTY DIGITAL FLOOR SCALE (PACK OF 4)



With its non-slip surface, the 880KL is the best choice for large or unstable patients. The robust design offers the largest platform available and an extra-large digital display. The display's size makes it easier to read than any of our other bath scales and is great for use in a clinic setting or for home health.

Capacity: 500 lb / 220 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

17" (w) x 13" (d) x 1 ¼" (h)

432 mm (w) x 330 mm (d) x 33 mm (h)

Display: 2 ½" / 64 mm LCD

Options & Accessories:

Power Source: (2) 3V lithium batteries included

Functions: LB / KG / ST Switch, Auto Off, 10 Second Weight Hold, Zero

Model	Description
880KLS	Heavy-duty digital floor scale (single)
64769	Carrying case for 880KL



All Health o meter® Professional scales are manufactured in FDA-registered factories.

844KL

DIGITAL FLOOR SCALE (PACK OF 2)



The 844KL offers a high capacity paired with the highest resolution in the bath scale line. A good option for patients who need greater precision in their weight measurement.

Capacity: 440 lb / 200 kg

Resolution: 0.1 lb / 0.05 kg

Platform Size:

12 5/8" (w) x 12 5/8" (d) x 1 7/8" (h)

321 mm (w) x 321 mm (d) x 48 mm (h)

Display: 1 1/4" / 32 mm LCD

Connectivity: USB

Options & Accessories:

Power Source: (2) 3V lithium batteries included

Functions: LB / KG / ST Switch, Auto Zero, Auto Off

Model	Description
844KLS	Digital floor scale (single)
64771	Carrying case for 844KL

349KLX

DIGITAL FLOOR SCALE



The 349KLX is a versatile scale that can be used in a host of locations. Users can place the remote display head where it is easiest for them to read the weight. The scale's small footprint makes it a great choice for home use.

Capacity: 400 lb / 180 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size:

12 1/2" (w) x 12" (d) x 1 7/8" (h)

318 mm (w) x 305 mm (d) x 48 mm (h)

Display: 1" / 25 mm LCD

Connectivity: Serial Port

Options & Accessories:

Power Source: 6 AA batteries included, 120V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, Zero, Tare, Auto Zero, Auto Off

Model	Description
349KLXAD	Scale with adapter included
SS-349KLX	2 year extended warranty
ADPT40	Adapter for 349KLX
64771	Carrying case for 349KLX



322KL



★ TOP SELLER

822KL

DIGITAL FLOOR SCALE (PACK OF 2)

This easy to transport scale offers high resolution with a wide platform to make it easier for patients with a wide stance. A versatile scale that is ideal for many uses including traveling healthcare professionals, diet groups, and home use.

Large 1 ½" LCD Display: The large display shows the user an easy-to-read weight measurement.

5-Second Weight Hold: Weight measurement will stay displayed for 5 seconds allowing the user time to stabilize themselves after stepping off the scale.

Low Profile Platform: Low platform height requires only minimal step up, offering patient a comfortable and stable weighing experience.

Capacity	Resolution	Platform Size
400 lb / 180 kg	0.1 lb / 0.05 kg	12 ¾" (w) x 10 ¾" (d) x 1 ½" (h) 349 mm (w) x 273 mm (d) x 38 mm (h)

Display
1 ½" LCD

Options & Accessories:

Model	Description
822KLS	Digital floor scale (single)
64771	Carrying case for 822KL

Power Source: (2) AAA batteries included

Functions: LB / KG Switch, Auto Zero, 5 Second Weight Hold, Auto Off

100LB / 100KG

MECHANICAL FLOOR SCALE -
POUNDS ONLY / KILOGRAMS ONLY
(PACK OF 3)



Simple and lightweight, these mechanical scales are ideal for home use or by visiting nurses.

Capacity:

100LB: 270 lb

100KG: 120 kg

Graduation:

100LB: 1 lb

100KG: 1 kg

Platform Size:

10 ¼" (w) x 10" (d) x 1 7/8" (h)

260 mm (w) x 254 mm (d) x 48 mm (h)

Functions:

100LB: Measures weight in pounds.

100KG: Measures weight in kilograms.

Dial Size: 2" / 51 mm

Options and Accessories:

Model	Description
100LBS	Mechanical floor scale – pounds only (single)
100KGS	Mechanical floor scale – kilograms only (single)
64771	Carrying case for 100LB/100KG

175LB

MECHANICAL FLOOR SCALE -
POUNDS ONLY (PACK OF 2)



When only pounds measurements are required, the 175LB offers a large, easy to read dial on a low profile platform. A great product for home use.

Capacity: 330 lb

Graduation: 1 lb

Platform Size:

11 3/8" (w) x 9 3/4" (d) x 2" (h)

289 mm (w) x 248 mm (d) x 51 mm (h)

Functions: Mechanically measures weight in pounds.

Dial Size: 3" / 89 mm

Options and Accessories:

Model	Description
64771	Carrying case for 175LB

All Health o meter® Professional scales are manufactured in FDA-registered factories.

142KL

MECHANICAL FLOOR SCALE
(PACK OF 2)



For use in the home or clinic, the 142KL is a good option for patients who require a spacious platform and large dial. The scale is also a good choice for weight loss programs.

Capacity: 330 lb / 150 kg

Graduation: 1 lb / 0.5 kg

Platform Size:

11 ½" (w) x 10 ¾" (d) x 2 ¼" (h)

292 mm (w) x 273 mm (d) x 57 mm (h)

Dial Size: 6 ¼" / 159 mm

Options and Accessories:

Functions: Measures weight in pounds and kilograms.

Model	Description
142KLS	Mechanical floor scale (single)

160KL

MECHANICAL FLOOR SCALE
(PACK OF 2)



This classic "Bigfoot" scale has a deep platform with non-slip surface and an 8" easy-to-read raised dial—good for patients with vision challenges or those who prefer more foot room.

Capacity: 400 lb / 180 kg

Graduation: 1 lb / 1 kg

Platform Size:

11" (w) x 12 ½" (d) x 3" (h)

279 mm (w) x 318 mm (d) x 76 mm (h)

Dial Size: 8" / 203 mm, raised 4" / 102 mm

Options & Accessories:

Functions: Measures weight in pounds and kilograms.

Model	Description
160LB	Mechanical floor scale—pounds only (2 pack)
160KG	Mechanical floor scale—kilograms only (2 pack)
160KLS	Mechanical floor scale (single)
160LBS	Mechanical floor scale—pounds only (single)
160KGS	Mechanical floor scale—kilograms only (single)

All Health o meter® Professional scales are manufactured in FDA-registered factories.



Height Measuring



★ TOP SELLER

201HR

WALL-MOUNTED TELESCOPIC METAL HEIGHT ROD

Sleek and robust measuring device requires minimal wall space and has a foldable headpiece. Designed for long life and to withstand dozens of uses every day.

Durable Construction: Aluminum height rod withstands frequent use. Rigorous laboratory tests showed that Health o meter® Professional's 201HR lasts 2 ½ times longer than similar products.

Foot Positioner: Attached footpiece serves as a guide for heel placement to ensure accurate measurements.

Large Measurement Range: The 2" - 90" / 6 cm - 230 cm measuring range provides height readings for a wide range of patients.

Wall Mount Design: Installs easily and quickly with simple tools and provided hardware.

Eye-Level Reading: Superior design ensures that caregivers of any height can easily read even the tallest patient's measurement without a problem.

Measuring Range

2" - 90"
6 cm - 230 cm

Graduation

1/16" / 1 mm

Options & Accessories:

Model	Descrip-
201HR	Wall-mounted telescopic metal height rod
201HR-400	Optional telescopic metal height ro for 400 series

205HR

HIGH-STRENGTH WALL-MOUNTED
HEIGHT ROD



Manufactured with a durable high-strength plastic polymer, this sturdy stadiometer resists bending and ensures that the headpiece slides smoothly. Featuring a double sided measuring tape, it offers a choice of centimeters only or feet and inches.

Measurement Range: 55.5" / 141 cm

Measuring Tape: 0" - 90" / 0 cm - 230 cm

Graduation: 1/16" / 1 mm



PORTROD

WALL-MOUNTED HEIGHT ROD



With its large headpiece and clear indicator the PORTROD stadiometer makes it easy to read a patient's height. Its wall mount design is easy and quick requiring only simple tools and provided hardware.

Measurement Range: 23 5/8" - 83" / 60 cm - 213 cm

Graduation: 1/16" / 1 mm

BABYHR

PEDIATRIC MEASURING ROD



The BABYHR is perfect for practitioners who need a tool to measure infant length. It can be mounted on the wall along an exam table and users can fold the headpiece and footpiece when not in use.

Measurement Range: 3" - 39" / 7 cm - 100 cm

Graduation: 1/16" / 1 mm

WMROD

402KLROD WITH ANCHORS AND SCREWS



Included on the 402KL scale, this metal height rod can be purchased separately to mount on the wall for adult or pediatric use.

Measurement Range: 23 5/8" - 84" / 60 cm - 213 cm

Graduation: 1/8" / 1 mm

* Information about height rods for specific scales can be found on each scale's product page.



Accessories

64769

CARRYING CASE FOR 549KL, 880KL



Features: Nylon bag protects scale during transport. Convenient for mobile use

553CASE

CARRYING CASE FOR 553KL & 550KL



Features: Nylon bag protects scale during transport. Convenient for mobile use

64771

CARRYING CASE FOR 349KLX, 498KL, 752KL, 100LB/KG, 175LB, 800KL, 822KL, 844KL, 8320KL



Features: Reinforced nylon bag protects scale during transport.

2210CART

ROLLING SCALE / MEDICAL DEVICE CART FOR USE WITH 2210KL, 2210KL4, 553KL, 549KL, 550KL



Features: Lockable cabinet and storage shelf.

Cabinet Dimensions:

23" (w) x 16 3/4" (d) x 17" (h)

584 mm (w) x 427 mm (d) x 432 mm (h)

Capacity: 200 lb / 90 kg

Cart Dimensions:

26 1/2" (w) x 18" (d) x 34 1/8" (h)

673 mm (w) x 457 mm (d) x 867 mm (h)

Wheels: 4 Non-conductive wheels (2 with locks)

ADPT30

100-240V POWER ADAPTER FOR USE WITH:
553KL, 1100KL, 2000KL, 2101KL, 2400KL,
2500KL, 2500CKL, 2650KL, 2700KL,
3001KL-AM, 3001KL-AMX

ADPT31

POWER ADAPTER FOR USE WITH: 498KL,
499KL, 500KL, 502KL, 522KL, 524KL,
2842KL

ADPT40

120V POWER ADAPTER FOR USE WITH:
349KLX

ADPT50

100-240V POWER ADAPTER FOR USE WITH:
597KL, 599KL, 600KL, 600KL-BT, ELEVATE-BT,
752KL, 3102KL-AM, ELEVATE-C AND ELEVATE-L

PLUGSET30

INTERNATIONAL WALL RECEPTACLE
ADAPTERS FOR SCALES THAT USE
ADPT30 OR ADPT50 POWER SOURCE

55000

WHEEL BRACKET FOR HEALTH O METER®
PROFESSIONAL 400 SERIES BEAM SCALES

Features: Wheel bracket attaches easily to scale for mobility.

MUACTAPE

MIDDLE UPPER ARM CIRCUMFERENCE TAPE



Features: Measures nutritional status in children
ages 10 - 14 years old.

55070

COUNTERWEIGHTS FOR 450KL
MECHANICAL BEAM SCALE

Features: Increases capacity of 450KL beam scales.

55090

COUNTERWEIGHTS FOR 400KL AND 402KL

Features: Increases capacity of 400KL and 402KL
beam scales.

B2500RAMP

SECONDARY RAMP FOR 2500KL, 2500CKL

B2600RAMP

SECONDARY RAMP FOR 2600KL

PELSTARONE-P

WIRELESS THERMAL PRINTER

Features: Lightweight printer connects to compatible scales via
Pelstar wireless technology. Powered by a rechargeable lithium
battery, the device prints patient weight, height and BMI (if
available) on 58 mm (2 ¼") thermal paper.

BMICHARTKG

BMI WHEEL (METRIC)



Features: Handy tool quickly displays BMI by selecting patient
height and weight. (metric units only)

PRODUCT SPECIFICATIONS

Stand-On Scales

	499KL	599KL	500KL	597KL	600KL/ELEVATE
					
FEATURES					
CAPACITY	500 lb / 220 kg	600 lb / 272 kg	550 lb / 250 kg	600 lb / 272 kg	600 lb / 272 kg
RESOLUTION	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
HEIGHT ROD INCLUDED			•	•	•
DISPLAY SIZE	1" / 25 mm	1" / 25 mm	1" / 25 mm	1" / 25 mm	1" / 25 mm
WHEELS	•	•	•	•	•
CONNECTIVITY/EMR	•	•	•	•	•
POWER ADAPTER INCLUDED		•		•	•
BATTERIES INCLUDED	•		•		•
WARRANTY	2	2	2	2	2
WIRELESS READY			500KL-BT		600KL-BT/ELEVATE-BT
	(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)
PLATFORM SIZE	13 3/4" x 16 1/2" x 2 3/8"	14 1/4" x 14 1/4" x 2 5/8"	13 3/4" x 16 1/2" x 2 3/8"	14 1/4" x 14 1/4" x 2 5/8"	14 1/4" x 14 1/4" x 2 5/8"
PRODUCT FOOTPRINT	14 1/8" x 21 1/8" x 43 1/2"	14 3/4" x 23 1/4" x 42 1/2"	14 1/8" x 21 1/8" x 53 5/8"	14 3/4" x 23 1/4" x 54 1/2"	14 3/4" x 24 1/8" x 54 3/8"
PRODUCT WEIGHT	20 lb	20 lb	21 lb	24 lb	28 lb
FUNCTIONS					
LB/KG CONVERSION	•	•	•	•	•
LB/KG LOCKOUT	•	•	•	•	•
HOLD/RELEASE	•	•	•	•	•
BMI	•	•	•	•	•
ZERO	•	•	•	•	•
TARE	•	•	•	•	•
PRETARE			•		
AUTO ZERO	•	•	•	•	•
AUTO OFF	•	•	•	•	•
REWEIGH			•		
RECALL			•		
EVERLOCK®			•		
ACCESSORIES (OPTIONAL)					
POWER ADAPTER	ADPT31	ADPT50	ADPT 31	ADPT50	ADPT50
HEIGHT ROD	499KLROD	499KLROD			
WRAP-AROUND HANDLEBAR			•		
WIRELESS KIT	•	•	•	•	•
WIRELESS PRINTER	•	•	•	•	•
	PAGE 6	PAGE 6	PAGE 5	PAGE 7	PAGE 9

1100KL

2101KL

3001KL-AM

3001KL-AMX

3102KL-AM



FEATURES

CAPACITY	1000 lb / 454 kg				
RESOLUTION	0.2 lb / 0.1 kg				
DISPLAY SIZE	1 ½" / 38 mm	¾" / 19 mm			
WHEELS	•	•	•	•	•
CONNECTIVITY/EMR	•	•	•	•	•
POWER ADAPTER INCLUDED	•	•	•	•	•
BATTERIES INCLUDED					•
WARRANTY	2	2	2	2	2
ANTIMICROBIAL			•	•	•
ANTIMICROBIAL KEYPAD	•	•	•	•	•

PRODUCT DIMS

	(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)
PLATFORM SIZE	22" x 15 ¾" x 2 ½"	26" x 22" x 2 ¾"	22" x 15 ¾" x 2 ⅝"	22" x 15 ¾" x 2 ⅝"	19 ¼" x 21 ¾" x 1 ¾"
PRODUCT FOOTPRINT	31 ¼" x 31 ¼" x 49 ½"	31" x 27" x 50"	22" x 33 ⅝" x 52"	22" x 34 ½" x 52"	23 ½" x 34 ¼" x 48"
PRODUCT WEIGHT	51 lb	66 lb	42 lb	43 lb	128 lb

FUNCTIONS

LB/KG CONVERSION	•	•	•	•	•
LB/KG LOCKOUT	•	•	•	•	•
HOLD/RELEASE	•	•	•	•	•
BMI	•	•	•	•	
ZERO	•	•	•	•	•
TARE	•	•	•	•	
PRE-TARE	•	•	•	•	•
AUTO ZERO	•	•	•	•	•
AUTO OFF	•	•	•	•	•
REWEIGH	•	•	•	•	•
RECALL	•	•	•	•	•
EVERLOCK®	•	•	•	•	

ACCESSORIES (OPTIONAL)

POWER ADAPTER	ADPT30	ADPT30	ADPT30	ADPT30	ADPT 50
HEIGHT ROD	PROPLUSROD	STROD	3001HR		
WIRELESS KIT	•	•	•	•	•
WIRELESS PRINTER	•	•	•	•	•

	349KLX	498KL	752KL	2842KL
--	--------	-------	-------	--------



FEATURES

- CAPACITY
- RESOLUTION
- DISPLAY SIZE
- CONNECTIVITY/EMR
- POWER ADAPTER INCLUDED
- BATTERIES INCLUDED
- WARRANTY

400 lb / 180 kg	500 lb / 220 kg	600 lb / 272 kg	600 lb / 270 kg
0.2 lb / 0.1 kg			
1" / 25 mm	1" / 25 mm	1" / 25 mm	7/8" / 22 mm
Serial Port	•	•	
		•	•
•	•		
2	2	2	2

PRODUCT DIMS

- PLATFORM SIZE
- PRODUCT FOOTPRINT
- PRODUCT WEIGHT

(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)
12 1/2" x 12" x 1 7/8"	13 3/4" x 16 1/2" x 2 3/8"	14 1/4" x 14 1/4" x 2 5/8"	22 1/4" x 42" x 2 1/2"
12 1/2" x 12" x 1 7/8"	14" x 17 1/4" x 2 3/8"	14 1/2" x 14 1/2" x 2 3/4"	22 1/4" x 42" x 2 1/2"
10 lb	14 lb	14 lb	53 lb

FUNCTIONS

- LB/KG CONVERSION
- LB/KG LOCKOUT
- HOLD/RELEASE
- BMI
- ZERO
- TARE
- AUTO ZERO
- AUTO OFF
- REWEIGH

•	•	•	
•	•	•	•
	•	•	
	•	•	
•	•	•	
•	•	•	•
•	•	•	•
•	•	•	
			•

ACCESSORIES (OPTIONAL)

- POWER ADAPTER
- CARRYING CASE
- HEIGHT ROD
- WIRELESS KIT
- WIRELESS PRINTER

ADPT 40	ADPT 31	ADPT 50	ADPT 31
64771	64771	64771	
*****See pages 51 and 52 for height measuring instruments*****			
	•	•	
	•	•	

	400KL	402KL	402LB	402KG SPECIAL ORDER ONLY	450KL
--	-------	-------	-------	-----------------------------	-------



FEATURES

- CAPACITY
- GRADUATION
- HEIGHT ROD INCLUDED
- ROTATING POISE BARS
- ENGRAVED NUMERALS
- WARRANTY

390 lb / 180 kg	390 lb / 180 kg	400 lb	200 kg	500 lb / 200 kg
¼ lb / 100 g	¼ lb / 100 g	¼ lb	100 g	¼ lb / 100 g
	•	•	•	•
				•
				•
10	10	10	10	10

PRODUCT DIMS

- PLATFORM SIZE
- PRODUCT FOOTPRINT
- PRODUCT WEIGHT

	(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)
PLATFORM SIZE	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"
PRODUCT FOOTPRINT	20 ¼" x 10 ⅜" x 57 ½"	20 ¼" x 10 ⅜" x 57 ½"	20 ¼" x 10 ⅜" x 57 ½"	20 ¼" x 14" x 57 ½"	20 ¼" x 14" x 57 ½"
PRODUCT WEIGHT	24 lb	28 lb	28 lb	28 lb	31 lb

ACCESSORIES (OPTIONAL)

- WHEELS
- HEIGHT ROD
- COUNTERWEIGHTS

WHEELS	55000	55000	55000	55000	55000
HEIGHT ROD	402KLROD or 201HR-400				
COUNTERWEIGHTS	55090	55090			

PAGE 16

PAGE 15

PAGE 15

PAGE 15

PAGE 16

PRODUCT SPECIFICATIONS

Pediatric Scales

	522KL	524KL	553KL	549KL
--	-------	-------	-------	-------



FEATURES

- CAPACITY
- RESOLUTION
- MEASURING TAPE INCLUDED
- DISPLAY SIZE
- CONNECTIVITY/EMR
- BATTERIES INCLUDED
- WARRANTY

CAPACITY	50 lb / 23 kg	50 lb / 23 kg	44 lb / 20 kg	44 lb / 20 kg
RESOLUTION	0-20 lb / 0.2 oz; 20-50 lb / 0.5 oz; 0-9 kg / 5 g; 9-23 kg / 10 g	0-20 lb / 0.2 oz; 20-50 lb / 0.5 oz; 0-9 kg / 5 g; 9-23 kg / 10 g	0.5 oz / 10 g	0.5 oz / 10 g
MEASURING TAPE INCLUDED	•	•	•	•
DISPLAY SIZE	1" / 25 mm	1" / 25 mm	1" / 25 mm	5/8" / 15 mm
CONNECTIVITY/EMR	•	•	•	•
BATTERIES INCLUDED	•	•	•	•
WARRANTY	2	2	2	2

PRODUCT DIMS

- TRAY/SEAT SIZE
- PRODUCT FOOTPRINT
- PRODUCT WEIGHT

	(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)
TRAY/SEAT SIZE	24 1/8" x 14 1/2" x 2 5/8"	24 1/4" x 14 1/4" x 18 1/2"	24 7/8" x 13" x 1 3/8"	19 1/4" x 13" x 2 5/8"
PRODUCT FOOTPRINT	24 1/8" x 21" x 23 1/8"	24 1/2" x 21" x 23 1/8"	24 7/8" x 16 1/2" x 3 7/8"	19 3/8" x 13 1/2" x 4 3/8"
PRODUCT WEIGHT	20 lb	21 lb	10 lb	4 lb

FUNCTIONS

- LB/KG CONVERSION
- LB/KG LOCKOUT
- HOLD/RELEASE
- ZERO
- TARE
- AUTO ZERO
- AUTO OFF

LB/KG CONVERSION	•	•	•	•
LB/KG LOCKOUT	•	•	•	•
HOLD/RELEASE	•	•	•	•
ZERO	•	•	•	•
TARE	•	•	•	•
AUTO ZERO	•	•	•	•
AUTO OFF	•	•	•	•

ACCESSORIES (OPTIONAL)

- POWER ADAPTER
- CARRYING CASE
- HEIGHT ROD
- WIRELESS KIT
- WIRELESS PRINTER

POWER ADAPTER	ADPT31	ADPT31	ADPT30	
CARRYING CASE			553CASE	64769
HEIGHT ROD	BABYHR	BABYHR	BABYHR	BABYHR
WIRELESS KIT	•	•	•	
WIRELESS PRINTER	•	•	•	

PAGE 24

PAGE 24

PAGE 23

PAGE 25

2210KL

2210KL4

8320KL



FEATURES

CAPACITY	45 lb / 20 kg / 20,000 g	45 lb / 20 kg / 20,000 g	33 lb / 15 kg
RESOLUTION	1 g (to 10,000 g); 5 g (10,001 to 20,000 g); 0.1 oz (0-45 lb)	1 g (to 10,000 g); 5 g (10,001 to 20,000 g); 0.1 oz (0-45 lb)	0.5 oz / 10 g
MEASURING TAPE INCLUDED	●	●	
DISPLAY SIZE	7/8" / 22 mm	7/8" / 22 mm	3/4" / 19 mm
BATTERIES INCLUDED	●	●	●
WARRANTY	2	2	2
REMOVABLE TRAY	●	●	

PRODUCT DIMS

TRAY SIZE	(WxDxH) 24" x 15 3/8" x 3"	(WxDxH) 27" x 15 3/8" x 3"	(WxDxH) 23 1/2" x 11" x 7/8"
PRODUCT FOOTPRINT	25 1/2" x 17 3/8" x 5 3/8"	27" x 17 3/8" x 5 3/8"	29 1/8" x 12 3/8" x 2 3/8"
PRODUCT WEIGHT	24 lb	25 lb	5 lb

FUNCTIONS

LB/KG/GM LOCKOUT	●	●	
LB/KG CONVERSION	●	●	●
LB/KG LOCKOUT	●	●	
HOLD/RELEASE			
ZERO	●	●	●
TARE	●	●	●
AUTO ZERO	●	●	●
AUTO OFF	●	●	●
REWEIGH	●	●	
RECALL	●	●	

ACCESSORIES (OPTIONAL)

CARRYING CASE			64771
HEIGHT ROD	BABYHR	BABYHR	BABYHR
OTHER	2210CART	2210CART	

1522KL

1524KL

386S-01

386KGS-01

BABYSLINGSCALE



FEATURES

CAPACITY

130 lb / 65 kg

130 lb / 65 kg

50 lb

25 kg

25 kg

GRADUATION

1 oz / 20 g

1 oz / 20 g

¼ lb

100 kg

100 g

MEASURING TAPE INCLUDED

●

●

●

DIAL SIZE

3" / 76 mm

3" / 76 mm

5 ½" / 140 mm

WARRANTY

2

2

1

1

1

PRODUCT DIMS

TRAY/SEAT SIZE

(WxDxH)

24 ⅞" x 14 ½" x 3 ¼"

(WxDxH)

24 ½" x 14 ¼" x 18 ½"

(WxDxH)

19 ⅜" x 12 ⅜" x 3 ¾"

(WxDxH)

19 ⅜" x 12 ⅜" x 3 ¾"

(WxDxH)

12" x 43" x 1 ⅞"

PRODUCT FOOTPRINT

24 ⅞" x 20" x 22 ½"

24 ½" x 21" x 22 ½"

19 ⅜" x 14 ½" x 5 ⅞"

19 ⅜" x 14 ½" x 5 ⅞"

12" x 43" x 1 ⅞"

PRODUCT WEIGHT

25 lb

32 lb

10 lb

10 lb

2 lb

ACCESSORIES (OPTIONAL)

HEIGHT ROD

BABYHR

BABYHR

BABYHR

BABYHR

BABYHR

PAGE 27

PAGE 27

PAGE 28

PAGE 28

PAGE 28

PRODUCT SPECIFICATIONS

Wheelchair/Chair/Multi-Use

2500KL

2600KL

2610KL

2650KL

2700KL



FEATURES

CAPACITY	1000 lb / 454 kg				
RESOLUTION	0.2 lb / 0.1 kg				
DISPLAY SIZE	1 1/2" / 38 mm	1 1/2" / 38 mm	1 1/2" / 38 mm	3/4" / 19 mm	3/4" / 19 mm
WHEELS	•	•	•	•	•
CONNECTIVITY/EMR	•	•	•	•	•
POWER ADAPTER INCLUDED	•	•	•	•	•
WARRANTY	2	2	2	2	2

PRODUCT DIMS

PLATFORM SIZE	(WxDxH) 29 3/8" x 28 1/4" x 2 3/8"	(WxDxH) 36" x 32 1/4" x 2 1/2"	(WxDxH) 36" x 32 1/4" x 2 1/2"	(WxDxH) 31 1/2" x 31 1/2" x 2"	(WxDxH) 43" x 42" x 2 1/2"
PRODUCT FOOTPRINT	40 1/4" x 34 1/2" x 45 1/2"	48 1/8" x 43 7/8" x 49 3/8"	56 3/4" x 43 7/8" x 49 3/8"	51 3/4" x 33 5/8" x 5 1/4"	69 1/8" x 51 1/2" x 3 3/8"
PRODUCT WEIGHT	81 lb	93 lb	98 lb	82 lb	187 lb

FUNCTIONS

LB/KG CONVERSION	•	•	•	•	•
LB/KG LOCKOUT	•	•	•	•	•
HOLD/RELEASE	•	•	•	•	•
BMI	•	•	•	•	•
ZERO	•	•	•	•	•
TARE	•	•	•	•	•
PRE-TARE	•	•	•	•	•
AUTO ZERO	•	•	•	•	•
AUTO OFF	•	•	•	•	•
REWEIGH	•	•	•	•	•
RECALL	•	•	•	•	•

ACCESSORIES (OPTIONAL)

POWER ADAPTER	ADPT30	ADPT30	ADPT30	ADPT30	ADPT30
WIRELESS PRINTER	•	•	•	•	•
WIRELESS KIT	•	•	•	•	•
RAMP	B2500RAMP	B2600RAMP			

PAGE 32

PAGE 31

PAGE 31

PAGE 33

PAGE 33

	2400KL	2500CKL	2595KL	594KL	445KL
--	--------	---------	--------	-------	-------



FEATURES

CAPACITY	800 lb / 363 kg	1000 lb / 454 kg	600 lb / 270 kg	600 lb / 270 kg	440 lb / 200 kg
RESOLUTION	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	¼ lb / 100 g
DISPLAY SIZE	1 ½" / 38 mm	1 ½" / 38 mm LCD	1 ⅝" / 30 mm	1 ⅜" / 35 mm	
CONNECTIVITY/EMR	•	•	•		
WHEELS	•	•	•	•	•
POWER ADAPTER INCLUDED	•	•	•	•	
WARRANTY	2	2	2	2	2

PRODUCT DIMS

	(WxDxH) 25 ¾" - 50 ⅞" 654 mm - 1273 mm	(WxDxH) 29 ⅜" x 28 ¼" x 2 ⅜"	(WxDxH) 12 ¾" x 10 ½" x 3 ⅞"	(WxDxH) 18" x 18 ½"	(WxDxH) 10 ½" x 5 ¾" x 3"
PLATFORM/ SEAT SIZE					
PRODUCT FOOTPRINT	51 ¼" x 40 ½" x 3 ½"	40 ¼" x 34 ½" x 45 ½"	23 ½" x 37 ½" x 39 ⅞"	23 ¼" x 40 ⅞" x 37 ⅞"	24 ⅜" x 31 ⅝" x 39 ⅞"
PRODUCT WEIGHT	40 lb	121 lb	84 lb	39 lb	60 lb

FUNCTIONS

LB/KG CONVERSION	•	•	•	•	
LB/KG LOCKOUT	•	•	•	•	
HOLD/RELEASE	•	•		•	
BMI	•	•		•	
ZERO	•	•	•	•	
TARE	•	•	•	•	
PRE-TARE	•	•			
AUTO ZERO	•	•		•	
AUTO OFF	•	•	•	•	
REWEIGH	•	•	•		
RECALL	•	•			

ACCESSORIES (OPTIONAL)

POWER ADAPTER	ADPT 30	ADPT 30	ADPT 30	ADPT 30	
WIRELESS PRINTER	•	•	•		•
WIRELESS KIT	•	•	•	•	
RAMP		B2500RAMP			

PRODUCT SPECIFICATIONS

Body Composition Scales

FEATURES

CAPACITY	1000 lb / 453.5 kg
RESOLUTION	0.1 lb / 0.1 kg
BODY FAT ACCURACY	+/-2%
MEASUREMENT TIME	15 sec
POWER ADAPTER	•
6 WINDOW DISPLAY	•
LIMITED TEST RESTRICTIONS	•
OHMS MODULATION LOGIC	•
WARRANTY	1

MEASUREMENTS

WEIGHT	•
% BODY FAT	•
% BODY WATER	•
% MUSCLE MASS	•
BMI (OR BMI %)	•
REE / BMR	•
ACTIVITY-BASED CALORIC BURN RATE	•

IMPACT ON ACCURACY/ PRE-TEST RESTRICTIONS

RECENT FOOD/DRINK	Ok to Test
RECENT EXERCISE	Ok to Test
FULL NIGHT SLEEP	Ok to Test
RECENT ALCOHOL INTAKE	Ok to Test
FULL BLADDER	Ok to Test


REPORTS

PLAIN PAPER	•
ILLUSTRATED	•
OFF-THE-SHELF PRINTER	•

PRODUCT DIMS

	(WxDxH)
PLATFORM SIZE	14" x 16" x 2 7/8"
PRODUCT FOOTPRINT	18" x 23 1/2" x 41 1/4"
PRODUCT WEIGHT	27 lb

PAGE 19



USAGE

ADULT PATIENTS	•	•	•	•
PEDIATRIC PATIENTS			•	
SEGMENTED (INCL. ARMS & LEGS)		•		
UPPER BODY ONLY				•
SHOES OK				•

BCS-G6-adult BCS-G6-limbs BCS-G6-DUO BCS-G61-upper

•	•	•	•
		•	
	•		
			•
			•

PRODUCT SPECIFICATIONS

Professional Home Care

800KL

880KL

844KL

822KL



FEATURES

CAPACITY	397 lb / 180 kg	500 lb / 220 kg	440 lb / 200 kg	400 lb / 180 kg
RESOLUTION	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.1 lb / 0.05 kg	0.1 lb / 0.05 kg
DISPLAY SIZE	1 ½" / 38 mm	1" / 25 mm	1 ¼" / 32 mm	1 ½" / 38 mm
BATTERIES INCLUDED	●	●	●	●
CASE PACK QUANTITY	4	4	2	2
WARRANTY	1	1	1	1

PRODUCT DIMS

	(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)
PLATFORM SIZE	12 ¼" x 12 ¼" x 1 ⅛"	17" x 13" x 1 ¼"	12 ⅝" x 12 ⅝" x 1 ⅞"	14 ¼" x 11 ¾" x 1 ½"
PRODUCT FOOTPRINT	12 ¼" x 12 ¼" x 1 ⅛"	17" x 13" x 1 ¼"	12 ⅝" x 12 ⅝" x 1 ⅞"	14 ¼" x 11 ¾" x 1 ¾"
PRODUCT WEIGHT	4 lb	8 lb	4 lb	4 lb

FUNCTIONS

EMR CAPABLE			●	
AUTO ZERO	●	●	●	●
AUTO OFF	●	●	●	●

ACCESSORIES (OPTIONAL)

CARRYING CASE	64771	64769	64771	64771
---------------	-------	-------	-------	-------

PAGE 40

PAGE 40

PAGE 41

PAGE 43

142KL**160KL****160LB****160KG****FEATURES**

CAPACITY
GRADUATION
DIAL SIZE
CASE PACK QUANTITY
WARRANTY

330 lb / 150 kg	400 lb / 180 kg	400 lb	180 kg
1 lb / 0.5 kg	1 lb / 1 kg	1 lb	1 kg
6 ¼" / 159 mm	8" / 203 mm	8" / 203 mm	8" / 203 mm
2	2	2	2
1	1	1	1

PRODUCT DIMS

PLATFORM SIZE
PRODUCT FOOTPRINT
PRODUCT WEIGHT

(WxDxH)	(WxDxH)	(WxDxH)	(WxDxH)
11 ½" x 10 ¾" x 2 ¼"	11" x 12 ½" x 3"	11" x 12 ½" x 3"	11" x 12 ½" x 3"
11 ¾" x 17 ¾" x 3 ½"	11 ¾" x 18 ¼" x 7 ¼"	11 ¾" x 18 ¼" x 7 ¼"	11 ¾" x 18 ¼" x 7 ¼"
6 lb	9 lb	9 lb	9 lb

PAGE 45**PAGE 45****PAGE 45****PAGE 45**

100LB**100KG****175LB****FEATURES**

CAPACITY
GRADUATION
DIAL SIZE
CASE PACK QUANTITY
WARRANTY

	100LB	100KG	175LB
CAPACITY	270 lb	120 kg	330 lb
GRADUATION	1 lb	0.5 kg	1 lb
DIAL SIZE	2" / 51 mm	2" / 51 mm	3 1/2" / 89 mm
CASE PACK QUANTITY	3	3	2
WARRANTY	1	1	1

PRODUCT DIMS

PLATFORM SIZE
PRODUCT FOOTPRINT
PRODUCT WEIGHT

	100LB	100KG	175LB
PLATFORM SIZE (WxDxH)	10 1/4" x 9 7/8" x 1 7/8"	10 1/4" x 9 7/8" x 1 7/8"	11 3/8" x 9 3/4" x 2"
PRODUCT FOOTPRINT	10 1/4" x 9 7/8" x 2"	10 1/4" x 9 7/8" x 2"	11 3/8" x 9 3/4" x 2 1/4"
PRODUCT WEIGHT	4 lb	4 lb	4 lb

ACCESSORIES (OPTIONAL)

CARRYING CASE

CARRYING CASE	64771	64771	64771
---------------	-------	-------	-------

PAGE 44**PAGE 44****PAGE 44**

PRODUCT SPECIFICATIONS

Specialty Scales

	2000KL	3400KL	3401KL	YG500R
				
FEATURES				
CAPACITY	400 lb / 181 kg	176 oz / 5000 g	176 oz / 5000 g	500 g
RESOLUTION OR GRADUATION	0.2 lb / 0.1 kg	0.05 oz / 1 g	0.05 oz / 1 g	2 g
DISPLAY SIZE	¾" / 19 mm	⅝" / 16 mm	⅝" / 16 mm	
WHEELS	•			
CONNECTIVITY	•			
POWER ADAPTER INCLUDED	•	•	•	
BATTERIES INCLUDED		Built-in Rechargeable	Built-in Rechargeable	
WARRANTY	2	2	2	1
PRODUCT DIMS				
	(WxDxH) Canvas Size 70 ½" x 30"	(WxDxH) Pan size 10" x 9 ⅞" x 2 ½"	(WxDxH) Tray Size 7 ¾" x 9 ¾" x ⅝"	(WxDxH) 6" x 4 ¾"
PLATFORM SIZE	70 ½" x 46 ½" x 60 ½"	11 ½" x 9 ⅞" x 8"	10" x 11 ¼" x 5 ⅝"	6 ½" x 6 ⅞" x 8 ⅝"
PRODUCT FOOTPRINT				
PRODUCT WEIGHT	114 lb	11 lb	10 lb	3 lb
FUNCTIONS				
LB/KG CONVERSION	•	•	•	
LB/KG LOCKOUT	•			
HOLD/RELEASE	•			
BMI	•			
ZERO	•	•	•	
TARE	•	•	•	
AUTO ZERO	•	•	•	
AUTO OFF	•	•	•	
REWEIGH	•			
RECALL	•			
ACCESSORIES (OPTIONAL)				
POWER ADAPTER	ADPT30	3400ADPT	3400ADPT	
OTHER		3400PAN (3-pack pans) 3400TRAY (3-pack trays)	3400PAN (3-pack pans) 3400TRAY (3-pack trays)	

PRODUCT SPECIFICATIONS

Height Measuring Instruments

	205HR	201HR	PORTROD	BABYHR	WMROD
					
FEATURES					
MEASUREMENT RANGE	55.5" / 141 cm	2" - 90" / 6 cm - 230 cm	23 5/8" - 83" / 60 cm - 213 cm	3" - 39" / 7 cm - 100 cm	23 5/8" - 84" / 60 cm - 213 cm
GRADUATION	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/8" / 1 mm
MEASURING TAPE RANGE	0" - 90" / 0 cm - 230 cm				
WARRANTY	2	2	2	2	2
PRODUCT DIMS					
PRODUCT FOOTPRINT	(WxDxH) 2 5/8" x 8" x 60"	(WxDxH) 2 1/4" x 11" x 54"	(WxDxH) 2 3/4" x 8 1/4" x 60"	(WxDxH) 41" x 11 1/2" x 2 1/8"	(WxDxH) 1 1/2" x 13 1/2" x 80 1/4"
	PAGE 48	PAGE 47	PAGE 48	PAGE 48	PAGE 48

Glossary/Symbols

Auto Off

Digital scale will automatically shut down after a designated period of time, reducing electrical usage and sustaining battery life.

Auto Zero

Ability of a digital scale to show a "0" value automatically when powered up.

BMI

A statistical measure that is the ratio between a person's weight and height and interpreted as a BMI score, a.k.a. Quetelet Index. Scales with this feature can calculate BMI score based on digitally measured or manually entered height data.

Connectivity

Ability of the scale to interface with a computer, vital signs monitor or other electronic device to transmit weight, height and/or BMI.

Factory Calibrated

All Health o meter® Professional scales are inspected and factory-calibrated prior to shipment.

Hold / Release

Digital function that allows the user to hold the patient weight until the key is pressed again to release that information.

LB / KG, LB / KG / GM, or G / OZ Conversion

Allows weight to be displayed in different units while patient remains on the scale.

LB / KG Lockout

Allows an administrator to select and use only pounds / inches or kilograms / centimeters, helping to reduce the risk of data being recorded in the incorrect unit of measurement.

LB / KG / ST or LB / KG Switch

Switch located on the scale platform that allows the user to select weight measuring unit.

Live Handrails

Allows patient to hold onto handrails for stability while scale accurately calculates weight.

Low Profile Platform

Low platform height requires a minimal step up, offering patient a comfortable and stable weighing experience.

Motion Detection

Digital feature that enables accurate weight despite movement of the subject on the scale.

Recall

Displays the last weight measurement.

Remote Display

The scale's display head is separate from the weighing platform and can be handheld, mounted on a wall or placed on a table.

Resolution / Graduation

The degree of difference/minimum weight change that a scale maintains when measuring weight, such as fraction of an ounce or pound.

Reweigh

Allows repeated weighing while the patient remains on the scale.

Swivel Display

Display head rotates for the convenience of the caregiver, allowing line-of-sight viewing and information privacy.

Zero

Zeroes the display to set the weight/reference point to zero.

Tare

Subtracts the weight of an object, such as a diaper or wheelchair, from the patient's weight measurement.

Pre-tare

Entering a known weight of an object such as a wheelchair, to be subtracted from the patient's weight.

Everlock®

Entering a known weight of an object such as a wheelchair, to be subtracted from the patient's weight.



EMR Capable - Products bearing this symbol can connect to EMR systems via multiple connectivity solutions. To learn more about connecting Health o meter® Professional scales, see Connectivity Solutions on page 3 or contact Technical Support at 1 (800) 638-3722.



CSA (us) - These products are certified by a CSA accredited third party lab, and have met applicable standards for the United States and Canada.



cMETus - Products have been evaluated for safety by MET, an accredited independent lab, in accordance with the appropriate national standards for the United States and Canada.



Ships By Truck Only - Due to their size, products can only be shipped by a common carrier. Please call 1 (800) 815-6615 for specific shipping information to your facility.



TAA Compliant - The manufacturing of these products meets the rules and requirements of the U.S. Trade Agreements Act.


All Health o meter® Professional products are listed with the US Food & Drug Administration.

All Health o meter® Professional products are manufactured without the use of natural rubber latex.

Global Support

Health o meter® Professional Scales has over 170 skilled professionals based in the United States and Canada who manufacture, sell, service and support our products. Our global headquarters in McCook, IL is home to our highly trained Customer Service and Technical Support teams along with our engineering, logistics and management team. On-site sales and support are provided in partnership with MedTech/MedCare, an industry-recognized leader in medical sales throughout the US and Canada.

Our American engineered products are designed and manufactured to meet the requirements of the International Organization for Standardization (ISO).



Health o meter® Professional Scales takes pride in designing world-class products for healthcare professionals around the globe. Our sales and support network includes select dealers and distributors throughout Latin America, Asia, Africa and the Middle East.



REHENSIVE
onsive
ty
BLE
AVAILABLE
standardization
weigh easier[®]

warranty

SUPPORT

WEIGH

EASIER[®]

reliable TIMELY

quality service

established FLEXIBLE

COMPREHENSIVE

connected reliable

RESPONSIVE precise

accurate easy

support

xible

ONSIVE

lization



Health o meter®

P r o f e s s i o n a l

Contact Information

Main Number: 1 (708) 377-0600

Customer Service: 1 (800) 815-6615

Technical Support: 1 (800) 638-3722

Fax: 1 (708) 377-0601 or 1 (800) 654-7330

Email: homprocs@homscales.com

Corporate Headquarters:

Health o meter® Professional Scales

9500 West 55th St.

McCook, IL 60525-7110 USA