


# Introducing the Patient Transfer Scale

- A transfer board with a built-in medical scale provides a weight measurement when transferring a patient
- Presents an accurate weight when it is critical for safe medication management for time sensitive patients
- Can be utilized as part of an existing patient transfer process



---

Weigh critical or immobile  
patients instantly and accurately

---

*Making Healthcare Professionals' Workflow Weigh**asier***<sup>®</sup>

**Health o meter**<sup>®</sup>  
Professional



UK PATENT  
Patent number 2543835



US PATENT PENDING



EU PATENT PENDING



CANADA PATENT PENDING



HONG KONG PATENT PENDING

# The Patient Transfer Scale: Every Minute Counts



## Key Features / Specifications:

- **Large, clear digital weight display**
- **Portable**
- **Wipe clean surface**
- **Capacity:** 550lb / 250kg
- **Resolution:** 1lb / 0.5kg
- **Power:** Internal rechargeable battery
- **Battery life:** 50 hours
- **Weight:** 25lb / 11.4 kg
- **Dimensions:**  
71" x 27 1/2" x 1" / 1805mm x 700mm x 30 mm
- **Model Number:**  
PTS-1000KL

- When managing a patient with a critical condition like a stroke, **Every Minute Counts.** It is imperative to get an accurate weight quickly when the patient is first admitted to the hospital. The Patient Transfer Scale (PTS) does just that.
- The Patient Transfer Scale was developed and invented by an Emergency Care nurse, with the final design based on feedback from healthcare professionals  
... **For Nurses by Nurses.**
- All departments can benefit from the Patient Transfer Scale, particularly Stroke Units, the Emergency Room, ICU, and Radiology.
- Learn more about how the PTS works Visit <https://www.homscales.com/products/patient-transfer-scale>
- Get to know the inventor of the PTS Gillian Taylor. Read her success stories.  
 **Twitter: @gilliantaym999**  
 **Instagram: @gilliantaylorpts**

**Scott Gottman**  
Senior Product Manager  
(708) 297-6582  
sgottman@homscales.com

**Health o meter® Professional Scales**  
(800) 815-6615  
www.homscales.com