

## THE 2210KL-AM-C

# ANTIMICROBIAL HIGH-RESOLUTION DIGITAL NEONATAL/PEDIATRIC TRAY SCALE

The 2210KL-AM-C neonatal scale was designed for use in the NICU, nursery, delivery room, and pediatric practices, where a high degree of accuracy is critical to the smallest of patients. The scale gives providers and patients safety with an Antimicrobial body, tray, handles, and keypad. The smooth-edged removable tray has a patent-pending mounting system that eliminates breakage, making it easy to clean. The built-in measuring tape allows users to measure both weight and height simultaneously. The scale will quickly lock despite the movement of the infant. The backlit display makes reading a patient's weight easier at almost any angle and light condition.

The high 1 gram resolution scale features a capacity of 20,000 grams/20 kg/45lb. The 2210KL-AM-C's user-friendly digital interface includes breast milk intake, auto-weight hold, hold/release, zero, tare, pre-tare, reweigh, recall, GM/KG/LB conversion and lockout, auto-zero, and auto-off. The

scale's motion detection gives it the ability to weigh accurately despite the movement of the patient. The exclusive Everlock® function allows the measuring units (LB, KG, or G) to be permanently locked on the scale.

Smooth  
Rounded  
Edges



The scale provides easy integration and seamless communication with today's EMR applications. The 2210KL-AM-C includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions. Options such as kilograms only, grams only, and built-in wireless models are available.

**Health o meter®**  
Professional

Health o meter® Professional 2210KL-AM-C

## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 2210KL-AM-C is your choice for accuracy and reliability

**1 GRAM  
RESOLUTION**

**120+ LOCAL  
FIELD SPECIALISTS**

**100%  
CONNECTIVITY**

**#1 TRUSTED  
BRAND**

**100 YEARS  
OF SERVICE**

weigh**easier**®

## TECHNICAL SPECIFICATIONS

### Capacity

45 lb / 20 kg / 20,000 gm

### Resolution

1 g (to 10,000 g);

5 g (10,001 to 20,000 g)

### EMR Connectivity

Via USB, optional Pelstar® wireless technology

### Backlit LCD Display

5 1/4" x 1" (133 mm x 35 mm)

4 levels of brightness

### Batteries

6 D batteries included

### Functions

LB / KG / GM conversion

LB / KG / GM lock out

Zero

Tare

Pre-tare

Reweigh

Recall

Auto zero

Auto off

Motion detection

Everlock®

Hold / release

Variable auto-off Time

Audible/mute sound option

Disable/enable backlight for display

Enabling unit view with unit of measure Locked

Auto-weight hold for 30 to 60 seconds

Breast milk intake

### Tray Size (w x d x h)

24" x 15 1/4" x 2 3/4"

(610 mm x 388 mm x 71 mm)

### Product Footprint

24 3/4" x 17 3/4" x 6"

(630 mm x 450 mm x 153 mm)

### Product Weight

28 lb / 13 kg

### Measuring Tape

0" - 24" / 0 cm - 60 cm

### Portability

Heavy duty handles for ease of moving

**2210CART-AC** - optional acute care cart

**2210CART** - optional rolling cart

### Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

MODEL	DESCRIPTION
2210KL-AM-C	Antimicrobial high resolution digital neonatal / pediatric tray scale, LB & KG
2210KG-AM-C	Antimicrobial high resolution digital neonatal / pediatric tray scale, KG only
2210KL-AM-BT-C	2210KL-AM-C with built-in Pelstar® wireless technology
2210KG-AM-BT-C	2210KG-AM-C with built-in Pelstar® wireless technology
2210CART-AC	Multi-Function Acute Care Cart for 2210KL-AM Series of Scales
2210CART	Rolling Scale/ Medical Device Cart
SS-2210KL	ScaleSurance® 2 year extended warranty for 2210 series*
PELSTARONE-P	Wireless thermal printer



### 2210CART-AC

Multi-Function Acute Care Cart for 2210KL-AM Series of Scales

To learn more, call **1-800-815-6615**  
or visit **homscales.com/2210KL-AM-C**

**Health o meter®**  
Professional

\*Not available in all countries.