

THE 245EHR-3001

DIGITAL HEIGHT ROD FOR THE 3001 SERIES OF SCALES

The Health o meter® Professional digital 245EHR-3001 is an optional state-of-the-art digital height rod for the Antimicrobial 3001 series of scales. The market-leading height rod provides the most accurate height measurement in its class, with a graduation/resolution of 1 mm/0.1" and a full measuring range of 11 ¾" – 79 ½" / 30cm – 202cm and a digital range of 3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm.



Built-In Reset Feature

Proprietary technology eliminates mechanical sensors, reducing failures, increasing accuracy and reliability

The new accessory uses solid-state technology, which provides a facility confidence that a patient's height will be measured accurately, and workflow will not be interrupted. The built-in reset feature increases accuracy and reliability, and the folding headpiece limits the height rod's profile, providing safety to patients and caregivers.



ORDERING INFORMATION

MODEL	DESCRIPTION
245EHR-3001	Optional digital height rod for 3001 Series*
SS-245EHR	ScaleSurance® 2 year extended warranty**

* Compatible with 3001 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 17 or higher.

Not available for AMX series scales. For assistance, contact Customer Service at 800-815-6615.

** Not available in all countries

Health o meter®
Professional

Health o meter® Professional 245EHR-3001

TECHNICAL SPECIFICATIONS

Digital Height Range

3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm

Resolution: 1 mm / 0.1"

Mechanical Height Range

11 ¾" – 79 ½" / 30cm – 202cm

Graduation: 1 mm / 1/16 inch

Product Footprint with 3001KL-AM

(height rod extended)

22" x 33 ½" x 82 ½"

(559 mm x 841 mm x 2096)

Product Footprint with 3001KL-AM

(height rod in home position)

22" x 33 ½" x 50"

(559 mm x 841 mm x 1270 mm)

Product Weight with 3001KL-AM

46 lb / 21 kg

Warranty

2-year limited manufacturer warranty

245EHR-3001 Compatible Scales

3001KL-AM

3001KG-AM

3001KL-AM-BT

3001KG-AM-BT

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/245EHR-3001**