



US Patent #D667323

Includes Connectivity Kit

External adapter allows scale to interface with Welch Allyn Connex® Vital Signs Monitors

Improve workflow with this top-of-the-line physician office scale. The built-in digital height rod eliminates the need to manually enter height for computing BMI. The 350° swivel head allows for line-of-sight viewing and patient privacy.

Advantages

✓ EMR Capable

Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions

✓ 350° Swivel Display

Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy

✓ Built-in Digital Height Rod

The easy-to-use, state of the art digital height rod provides measurements from 23 5/8" - 88" / 60 cm - 222 cm

✓ Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost

✓ Calculates Body Mass Index (BMI)

Scale functions include the ability to calculate the patient's BMI

✓ Patented Design

The 600KL-BT's patented design provides a user-friendly and intuitive interface for measuring patient weight, height and BMI with a convenient rotating display head



US Patent #D667323

Specifications

- Capacity: 660 lb / 300 kg
- Resolution: 0.2 lb / 0.1 kg
- Height Rod Resolution: 0.1" / 0.1 cm
- Platform Size (w x d x h):
14 1/4" x 14 1/4" x 2 5/8"
362 mm x 362 mm x 67 mm
- Product Footprint (w x h x d):
14 3/4" x 24 1/8" x 54 3/8"
375 mm x 613 mm x 1381 mm
- Product Weight: 28 lb / 13 kg
- Digital Height Rod: 23 5/8" - 88" / 60 cm - 222 cm
- Includes connectivity kit to interface with Welch Allyn® vital signs monitors
- Connectivity: USB and built-in Pelstar wireless technology
- 1" / 25 mm LCD display, 35° left / right swivel (350° with height rod fully retracted) for ease of viewing
- Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- 2 wheels for ease of movement
- 100 - 240V adapter (ADPT50) included, 6 C-cell batteries (not included)
- 2 year limited warranty

Model	Description
600KG-BT	Eye level scale with digital height rod & built-in Pelstar wireless technology, KG only
SS-600KL-BT	ScaleSurance 2 year extended warranty for 600KL-BT & 600KG-BT
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America