

THE 245EHR-1100

DIGITAL HEIGHT ROD FOR THE 1100 SERIES OF SCALES

The Health o meter® Professional digital 245EHR-1100 is an optional state-of-the-art digital height rod for the 1100 series of scales. The market-leading height rod provides the most accurate height measurement in its class, with a graduation/resolution of 1mm/0.1" and a full measuring range of 11 ¾" – 79 ½" / 30cm – 202cm and a digital range of 3ft 11.2" – 6ft 7.5" / 120 cm – 202 cm.



Built-In Reset Feature

Proprietary technology eliminates mechanical sensors, reducing failures, increasing accuracy and reliability

The new accessory uses solid-state technology, which provides a facility confidence that a patient's height will be measured accurately, and workflow will not be interrupted. The built-in reset feature increases accuracy and reliability, and the folding headpiece limits the height rod's profile, providing safety to patients and caregivers. The optional digital height rod has quick and easy to assemble brackets.



TECHNICAL SPECIFICATIONS

Digital Height Range

3ft 11.2" – 6ft 7.5" / 120 cm – 202 cm
Resolution: 1 mm / 0.1"

Mechanical Height Range

11 ¾" – 79 ½" / 30cm – 202cm
Graduation: 1mm / 1/16 inch

Product Footprint with 1100KL

(height rod extended)
31 ¼" x 27 ¼" x 82 ½"
(794 mm x 692 mm x 2096 mm)

Product Footprint with 1100KL

(height rod in home position)
31 ¼" x 27 ¼" x 50 ¼"
(794 mm x 692 mm x 1276 mm)

Product Weight with 1100KL

58 lb / 26 kg

Warranty

2-year limited manufacturer warranty

245EHR-1100 Compatible Scales*

1100KL
1100KG
1100KL-BT
1100KG-BT

ORDERING INFORMATION

| MODEL | DESCRIPTION |
|-------------|---|
| 245EHR-1100 | Optional digital height rod for 1100 Series |
| SS-245EHR | ScaleSurance® 2 year extended warranty** |

*Compatible with 1100 series scales with the letter "E" or "U" in the serial number and firmware versions 02.08.00 or higher. For assistance, contact Customer Service at 800-815-6615.

** Not available in all countries

Health o meter®
Professional

Health o meter® Professional 245EHR-1100

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/245EHR-1100**