

# THE 245EHR-2101

## DIGITAL HEIGHT ROD FOR THE 2101 SERIES OF SCALES

The Health o meter® Professional digital 245EHR-2101 is an optional state-of-the-art digital height rod for the 2101 series of scales. The market-leading height rod provides the most accurate height measurement in its class, with a graduation/resolution of 1mm/0.1" and a full measuring range of 11 ¾" – 79 ½" / 30cm – 202cm and a digital range of 3ft 11.2" – 6ft 7.5" / 120 cm – 202 cm.



### Built-In Reset Feature

Proprietary technology eliminates mechanical sensors, reducing failures, increasing accuracy and reliability

The new accessory uses solid-state technology, which provides a facility confidence that a patient's height will be measured accurately, and workflow will not be interrupted. The built-in reset feature increases accuracy and reliability, and the folding headpiece limits the height rod's profile, providing safety to patients and caregivers. The optional digital height rod has quick and easy to assemble brackets.



## TECHNICAL SPECIFICATIONS

### Digital Height Range

3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm  
Resolution: 1 mm / 0.1"

### Mechanical Height Range

11 ¾" – 79 ½" / 30cm – 202cm  
Graduation: 1mm / 1/16 inch

### Product Footprint with 2101KL

(height rod extended)  
31" x 27" x 82 ½"  
(787 mm x 686 mm x 2096 mm)

### Product Footprint with 2101KL

(height rod in home position)  
31" x 27" x 50"  
(787 mm x 686 mm x 1270 mm)

### Product Weight with 2101KL

70 lb / 30 kg

### Warranty

2-year limited manufacturer warranty

### 245EHR-2101 Compatible Scales

2101KL  
2101KG  
2101KL-BT  
2101KG-BT

## ORDERING INFORMATION

MODEL	DESCRIPTION
245EHR-2101	Optional digital height rod for 2101 Series*
SS-245EHR	ScaleSurance® 2 year extended warranty**

\*Compatible with 2101 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 17 or higher. For assistance, contact Customer Service at 800-815-6615.

\*\* Not available in all countries

**Health o meter®**  
Professional

Health o meter® Professional 245EHR-2101

weigh-easier®

To learn more, call **1-800-815-6615**  
or visit **homscales.com/245EHR-2101**