



Includes Connectivity Kit

500KL-BT includes a connectivity kit that allows the scale to interface directly with Welch Allyn[®] vital signs monitors.

The 500KL-BT features built-in wireless technology, no external hardware is required to attach to the scale. Utilizing Pelstar[®] wireless technology, the scale transmits data to Welch Allyn[®] Connex[®] devices or Windows[®] PCs.

Advantages

- ✓ **EMR Connectivity**
Utilizes Pelstar[®] wireless technology to transmit data to Welch Allyn[®] Connex[®] devices or Windows[®] PCs
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-Tare, Hold/Release, Auto Zero, Auto Off, Reweigh, Recall, Parent Child Function) included at no additional cost
- ✓ **Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience
- ✓ **Built-in Height Rod**
The easy-to-use height rod provides height measurements from 30" - 84 1/4" / 76 cm - 214 cm
- ✓ **Top-Selling Digital Scale**
Built-in wireless version of the Health o meter[®] Professional 500KL, the most widely used scale in physician offices nationwide
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI
- ✓ **Everlock[®] (patent pending)**
Exclusive feature allows the measuring unit (LB or KG) to be permanently locked on the scale.



Optional **500HB Handlebar Accessory** can be mounted to your new or existing 500KL, 500KG, 500KL-BT & 500KG-BT.



The 500 Series are the **ONLY** digital Primary Care Scales with a Handlebar Accessory



Specifications

- Capacity: 550 lb / 250 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity: Built-in Pelstar® wireless technology - includes components to wirelessly interface with a Welch Allyn® Connex monitor or Windows® PC
- Display: 1" / 25 mm LCD
- Functions: LB / KG Conversion, LB / KG Lock Out, Body Mass Index (BMI), Zero, Tare, Pre-Tare, Hold / Release, Auto Zero, Auto Off, Reweigh, Recall, Everlock®, Parent / Child Function
- Platform Size (w x d x h): 13 ¾" x 16 ½" x 2 ⅜" / 349 mm x 419 mm x 60 mm
- Product Foot Print (w x d x h): 14 ⅛" x 21 ⅛" x 53 ⅝" / 359 mm x 537 mm x 1362 mm
- Product Weight: 21 lb / 9 kg
- Height Rod: 30" - 84 ¼" / 76 cm - 214 cm
- 2 wheels for ease of movement
- 6 AA batteries included, Optional 100V - 240V adapter (order # ADPT30)
- 2 year limited warranty

Model	Description
500KG-BT	Eye level digital scale with built-in Pelstar® wireless technology, kg only
500HB	Optional handlebar accessory for 500 series scales
SS-500KL	ScaleSurance 2 year extended warranty (500 models)*
ADPT30	Optional 100-240V power adapter