## THE 500KL-BT

# PRIMARY CARE DIGITAL SCALE WITH BUILT-IN PELSTAR® WIRELESS TECHNOLOGY

The number-one-selling physician scale in North America, the all-digital Health o meter® 500KL-BT primary care scale provides the optimum in quality, reliability, and ease of use. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale features a built-in height rod in the scale column, a large non-skid platform, and reinforced base for added stability. The scale gives providers and patients an eye-level weight indicator for patients up to 660 pounds, with a user-friendly digital interface that includes BMI, hold/release, zero, tare, reweigh, recall, KG/LB, clear, and enter functions. The 500KL offers the best in its class extra-wide 8" wheelbase making transport easy. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The first healthcare scale to provide easy integration and seamless communication with today's EMR applications, the 500KL-BT features built-in wireless technology to interface with a computer, EMR software, or other electronic device.

With its low platform height, the 500KL-BT is easy to set up and offers patients a more comfortable and stable weighing experience.

#### **500HB**

"live" wrap-around handlebar

The 500HB is a unique handlebar accessory for the 500KL and is easy to attach to any new scale using the pre-drilled mounting holes and can also be attached to any existing scale using the easy-to-use drill guides. It provides the stability that many patients need without affecting the accuracy of the weight reading.

# Health o meter®

Professional

Health o meter® Professional 500KL-BT

THERE'S STRENGTH IN OUR NUMBERS.

**660 LB** CAPACITY

150,000+ INSTALLED NATIONWIDE

100% CONNECTIVITY

120+ LOCAL
FIELD SPECIALISTS

**8"** WIDE WHEELBASE BEST IN ITS CLASS

#1 TRUSTED
PHYSICIAN SCALE

100 YEARS OF SERVICE

weigheasier®

### **TECHNICAL SPECIFICATIONS**

**ORDERING INFORMATION** 

**DESCRIPTION** 

#### Capacity

660 lb / 300 kg

#### **Resolution**

0.2 lb / 0.1 kg

#### **EMR** connectivity

Built-in Pelstar® wireless technology

#### Display

1" / 25 mm, LCD display

#### **Batteries**

**MODEL** 

500HB

ADPT31

SS-500KL

PLUGSET30

500KG-BT

6 AA batteries included, optional adapter (order# ADPT31)

#### **Functions**

LB/KG conversion LB/KG lock out

BMI

Zero

Tare, Pre-tare

Hold/Release

Auto Zero

Auto Off

Reweigh

Recall

500KG with built-in Pelstar® wireless technology

Optional adapter for 500KL-BT & 500KG-BT

Plug adapters for use outside North America

Optional handlebar accessory for 500 series of scales

ScaleSurance® 2 year extended warranty for 500 series\*

Parent/Child weighing function Everlock® - permanent unit lock

#### Platform Size ( $w \times d \times h$ )

13 ¾" x 16 ½" x 2 ¾"

(349 mm x 419 mm x 60 mm) Dent and rust proof platform

#### **Product Footprint**

14 %" x 21 %" x 53 %"

(359 mm x 537 mm x 1362 mm)

#### **Product Weight**

21 lb / 9 kg

#### **Height Rod**

30"-84 1/4" / 76 cm-214 cm Graduation: 1/8" / 1 mm

#### **Portability**

2 wheels, wide 8" wheelbase for ease of movement

#### **Warranty**

2-year limited manufacturer warranty



Professional

# To learn more, call **1-800-815-6615** or visit **homscales.com/500KL-BT**

