THE 2101 KL

DIGITAL PLATFORM SCALE

The all-digital 2101KL is a comprehensive bariatric scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The extra wide "live" handrails allow the patient to hold on for stability while the scale accurately calculates weight.

Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large non-skid lowprofile platform, and wheels that enable portability.

The scale's high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock[®] function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications, the 2101KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter[®] Professional Connectivity Solutions.



When you add up the numbers on our scale, it's not surprising that the 2101 KL is your choice for a bariatric platform scale.

120+ LOCAL FIELD SPECIALISTS

100% CONNECTIVITY

#1 TRUSTED BRAND

100 YEARS OF SERVICE

weigh**easier**®



Health o meter® Professional 2101 KL

TECHNICAL SPECIFICATIONS

Capacity 1000 lb / 454 kg

Resolution 0.2 lb / 0.1 kg

EMR Connectivity Via USB, optional Pelstar® wireless technology

Display 1 ½" / 38 mm High-Contrast Color TFT-LCD Screen

Power Source 100-240V adapter (ADPT30) included, 6 D batteries not included

Functions

LB/KG conversion LB/KG lock out BMI Zero Tare Pre-tare Hold/release Reweigh Recall Auto zero Auto off Time/date Variable auto off time Audible/mute sound option Parent/child weighing function Everlock[®] - permanent weight lock

Platform Size ($w \ge d \ge h$) 26" x 22" x 2 3/8" (660 mm x 559 mm x 60 mm)

Product Footprint 31" x 27" x 50" (787 mm x 686 mm x 1270 mm)

Product Weight 66 lb / 30 kg

Portability 2 wheels for ease of movement

Warranty 2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
2101 KL	Digital platform scale, LB & KG
2101 KL-BT	2101 KL with Pelstar [®] wireless technology
2101KG	Digital platform scale, kilograms only
2101 KG-BT	2101KG with built-in Pelstar® wireless technology
2101 KLHR	2101KL with mechanical height rod included
2101KGHR	2101KG with mechanical height rod included
SS-2101 KL	ScaleSurance® 2 year extended warranty*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
245EHR-2101	Optional digital height rod for 2101 models**
STROD	Optional mechanical height rod for 2101 models
PELSTARONE-P	Wireless thermal printer

* Not available in all countries

** Compatible with 2101 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 17 or higher. For assistance, contact Customer Service at 800-815-6615.

To learn more, call 1-800-815-6615 or visit homscales.com/2101KL



Digital and Mechanical Height Rod Options Available



