

## **SCALES WITH CONNECTIVITY**

### **Scales Compatible with the Following Connectivity Options:**

- Wired connection to a PC via USB cable
- Wired connection to a Welch Allyn monitor via C-HOMWA-1 adapter
- Wireless connection to PC or Welch Allyn monitor via PELSTARONE wireless kit

| Model Number<br>(Previous versions also compatible.)                                | Scale Description  |
|---|--|
| 522KL / 522KG   | Pediatric Digital Tray Scale   |
| 524KL / 524KG   | Pediatric Digital Seat Scale   |
| 553KL / 553KG   | Pediatric Digital Tray Scale   |
| 597KL / 597KG   | Heavy Duty Eye Level Digital Scale   |
| 599KL / 599KG   | Heavy Duty Waist High Digital Scale  |
| 600KL / 600KG   | Heavy Duty Eye Level Digital Scale with Digital Height Rod   |
| 752KL / 752KG   | Heavy Duty Remote Display Digital Scale  |
| ELEVATE-C / ELEVATE-L /<br>ELEVATE-KG   | Heavy Duty Eye-Level Digital Scale with Digital Height Rod (includes connectivity kit for Welch Allyn devices) |
| 1100KL / 1100KG   | Digital Platform Scale   |
| 2000KL / 2000KG   | Digital In-Bed/Stretcher Scale   |
| 2101KL / 2101KG   | Digital Platform Scale   |
| 2400KL / 2400KG   | Digital Portable Wheelchair Scale  |
| 2500KL / 2500KG   | Digital Wheelchair Ramp Scale  |
| 2595KL  | Digital Chair Scale  |
| 2500CKL / 2500CKG   | Digital Wheelchair Ramp Scale with Fold Away Seat  |
| 2600KL / 2600KG   | Digital Wheelchair Ramp Scale  |
| 2610KL / 2610KG   | Digital Wheelchair Dual Ramp Scale   |
| 2650KL / 2650KG   | Digital Dual Ramp Wheelchair Scale   |
| 2700KL / 2700KG   | Digital Dual Ramp Wheelchair Scale with X-Large Platform   |
| 3001KL-AM / 3001KG-AM   | Antimicrobial Digital Platform Scale   |
| 3001KL-AMX / 3001KG-AMX   | Antimicrobial Digital Platform Scale with Extended Handrails   |
| Compatible discontinued SKUs: 4011, 4021, 4024, 4025, 4026, 4261, 4650, 4700, 753KL |  |

Continued on next page

| <b>Scales with Built-in Pelstar Wireless Technology</b> |  |
|---|--|
| <b>Model Number</b>                                     | <b>Scale Description</b>   |
| 500KL-BT / 500KG-BT                                     | Eye Level digital scale with Built-in Pelstar Wireless Technology  |
| 501KL-BT / 501KG-BT                                     | TAA Compliant Eye Level digital scale with Built-in Pelstar Technology   |
| 600KL-BT / 600KG-BT                                     | Eye Level Scale with Digital Height Rod and Built-in Pelstar Wireless Technology   |
| ELEVATE-BT / ELEVATE-KG-BT                              | Eye Level Scale with Digital Height Rod and Built-in Pelstar Wireless Technology (includes connectivity kit for Welch Allyn devices) |
| 1100KL-BT / 1100KG-BT                                   | Digital Platform Scale with Built-in Pelstar Wireless Technology   |
| 2000KL-BT / 2000KG-BT                                   | Digital In-Bed/Stretcher Scale with Built-in Pelstar Wireless Technology   |
| 2101KL-BT / 2101KG-BT                                   | Digital Platform Scale with Built-in Pelstar Wireless Technology   |
| 2400KL-BT / 2400KG-BT                                   | Digital Portable Wheelchair Scale with Built-in Pelstar Wireless Technology  |
| 2500KL-BT / 2500KG-BT                                   | Digital Wheelchair Ramp Scale with Built-in Pelstar Wireless Technology  |
| 2600KL-BT / 2600KG-BT                                   | Digital Wheelchair Ramp Scale with Built-in Pelstar Wireless Technology  |
| 2610KL-BT / 2610KG-BT                                   | Digital Wheelchair Scale with Dual Ramps and Built-in Pelstar Wireless Technology  |
| 2650KL-BT / 2650KG-BT                                   | Digital Wheelchair Dual Ramp Scale with Built-in Pelstar Wireless Technology   |
| 2700KL-BT / 2700KG-BT                                   | Digital Wheelchair Dual Ramp Scale with Extra Large Platform and Built-in Pelstar Wireless Technology                                |
| 2210KL-AM-BT / 2210KG-AM-BT                             | Antimicrobial High-Resolution Digital Neonatal/Pediatric Tray Scale with Built-in Pelstar Wireless Technology                        |