

## THE PTS-1000KL-C

# PATIENT TRANSFER SCALE

The PTS-1000KL-C is a transfer board with a built-in medical scale that provides a weight measurement when transferring a patient. All departments can benefit from the Patient Transfer Scale, particularly Stroke Units, the Emergency Room, ICU, and Radiology.

Designed by Clinicians for Clinicians the Patient Transfer Scale weighs critical or immobile patients instantly and accurately helping shorten door-to-needle times. It offers a superior alternative to other weighing equipment, such as in-bed scales, that are often not calibrated or inaccurate.

The scale is light and portable and can be utilized as part of an existing patient transfer process when a weight is needed, and will not disrupt the patient flow process.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that Health o meter® Professional Scales is the industry-leading name in medical scales.

**120+** LOCAL  
FIELD SPECIALISTS

**#1** TRUSTED  
BRAND

**100 YEARS**  
OF SERVICE

**Health o meter®**  
Professional

Health o meter® Professional PTS-1000KL-C

weigh**easier**®

## TECHNICAL SPECIFICATIONS

### Capacity

550 lb / 250 kg

### Resolution

1 lb / 0.5 kg

### Display

3" x 1" / 76 mm x 25 mm Backlight LCD

### Power Source

3.7V rechargeable battery

### Functions

LB/KG conversion

Zero

Hold/release

Auto off

### Weighing Mat Size (w x d)

71" x 21"

(1805mm x 533mm)

### Product Footprint (w x d x h)

71" x 27 1/2" x 1"

(1805mm x 700mm x 30mm)

### Product Weight

25 lb / 11.4 kg

### Measuring Tape

Height guide shown on mat 100cm – 180cm  
shown at 5 cm

### Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

MODEL	DESCRIPTION
PTS-1000KL-C	Patient Transfer Scale, LB & KG*
SS-PTS-1000KL	ScaleSurance® 2 year extended warranty for Patient Transfer Scale**



To learn more, call **1-800-815-6615**  
or visit **[homscales.com/PTS-1000KL-C](https://homscales.com/PTS-1000KL-C)**

\* Evaluated to SPE 3000 for sale in Canada

\*\*Not available in all countries.

**Health o meter®**  
Professional