

THE 2210KL-AM

ANTIMICROBIAL HIGH-RESOLUTION DIGITAL NEONATAL/PEDIATRIC TRAY SCALE

The 2210KL-AM neonatal scale was designed for use in the NICU, nursery, delivery room, and pediatric practices, where a high degree of accuracy is critical to the smallest of patients. The scale gives providers and patients safety with an Antimicrobial body, tray, handles, and keypad. The smooth-edged removable tray has a patent-pending mounting system that eliminates breakage, making it easy to clean. The built-in measuring tape allows users to measure both weight and height simultaneously. The scale will quickly lock despite the movement of the infant. The backlit display makes reading a patient's weight easier at almost any angle and light condition.

The high 1 gram resolution scale features a capacity of 20,000 grams/20 kg/45lb.

The 2210KL-AM's user-friendly digital interface includes breast milk intake, auto-weight hold, hold/release, zero, tare, pre-tare, reweigh, recall, GM/KG/LB conversion and lockout, auto-zero, and auto-off. The scale's motion detection gives it the ability to weigh accurately despite the movement of the patient. The exclusive Everlock[®] function allows the measuring units (LB, KG, or G) to be permanently locked on the scale.


2210KL4-AM

Scale with 4-Sided Tray Option

The scale provides easy integration and seamless communication with today's EMR applications. The 2210KL-AM includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter[®] Professional Connectivity Solutions. Options such as kilograms only, grams only, and built-in wireless models are available.

Health o meter[®]
Professional

Health o meter[®] Professional 2210KL-AM



Smooth
Rounded
Edges



**THERE'S STRENGTH IN OUR
NUMBERS.**

When you add up the numbers on our scale, it's not surprising that the 2210KL-AM is your choice for accuracy and reliability


**1 GRAM
RESOLUTION**

**120+ LOCAL
FIELD SPECIALISTS**

**100%
CONNECTIVITY**

**#1 TRUSTED
BRAND**

**100 YEARS
OF SERVICE**



weigh**easier**[®]

TECHNICAL SPECIFICATIONS

Capacity

45 lb / 20 kg / 20,000 gm

Resolution

1 g (to 10,000 g);

5 g (10,001 to 20,000 g)

EMR Connectivity

Via USB, optional Pelstar® wireless technology

Backlit LCD Display

5 ¼" x 1 " (133 mm x 35 mm)

4 levels of brightness

Batteries

6 D batteries included,
adapter ADPT30 included

Functions

LB / KG / GM conversion

LB / KG / GM lock out

Zero

Tare

Pre-tare

Reweigh

Recall

Auto zero

Auto off

Motion detection

Everlock®

Hold / release

Variable auto-off Time

Audible/mute sound option

Disable/enable backlight for display

Enabling unit view with unit of measure locked

Auto weight hold for 30 to 60 seconds

Breast milk intake

Tray Size (w x d x h)

(2210KL-AM) 25 ½" x 15 ¾" x 3 ¾"

(645mm x 403mm x 95mm)

(2210KL4-AM) 27 ½" x 15 ¾" x 3 ¾"

(700 mm x 405 mm x 95 mm)

Product Footprint

(2210KL-AM) 25 ½" x 18 ¼" x 6 ½"

(645mm x 460mm x 160mm)

(2210KL4-AM) 27 ½" x 18 ¼" x 6 ½"

(700 mm x 460mm x 160mm)

Product Weight

(2210KL-AM) 28 lb / 13 kg

(2210KL4-AM) 29 lb / 13.5 kg

Measuring Tape

0" - 24" / 0 cm - 60 cm

Portability

Heavy duty handles for ease of moving
2210CART - optional rolling cart

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
2210KL-AM	Antimicrobial high resolution digital neonatal / pediatric tray scale, LB & KG
2210KG-AM	Antimicrobial high resolution digital neonatal / pediatric tray scale, KG only
2210KL-AM-BT	2210KL-AM with built-in Pelstar® wireless technology
2210KG-AM-BT	2210KG-AM with built-in Pelstar® wireless technology
2210KL4-AM	Antimicrobial high resolution digital neonatal / pediatric 4 sided tray scale, LB & KG
2210KG4-AM	Antimicrobial high resolution digital neonatal / pediatric 4 sided tray scale, KG only
2210KL4-AM-BT	2210KL4-AM with built-in Pelstar® wireless technology
2210KG4-AM-BT	2210KG4-AM with built-in Pelstar® wireless technology
2210CART	Optional rolling cart for 2210 series
SS-2210KL	ScaleSurance® 2 year extended warranty for 2210 series*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer



2210CART
Optional Rolling Cart

To learn more, call **1-800-815-6615**
or visit homscapes.com/2210KL-AM

Health o meter®
Professional

*Not available in all countries.