THE 2900KL-AM

ANTIMICROBIAL WALL-MOUNTED FOLD-UP WHEELCHAIR SCALE

The number-one-selling medical scale company in North America offers an antimicrobial wall-mounted fold-up wheelchair scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The low profile scale extends off the wall 3.5" when

folded up and is excellent for facilities that need to save space. The large easy-to-lift platform accommodates most wheelchairs and provides a comfortable weighing experience for patients. The dual ramps offer easy on, easy off access and the anti-slip caution tape delivers enhanced safety. Other safety features include antimicrobial protection, a raised edge that follows ADA guidelines, a recessed handle, and no pinch points, which provide both the clinician and patient confidence during the weighing process.

Designed to meet the needs of hospitals, clinics, acute, long-term care facilities, doctor's offices, and anywhere space is limited. The scale

gives providers flexibility with a wall-mounted display that can be mounted on either side of the scale. The 2900KL-AM's high-contrast color screen provides a user-friendly interface that includes BMI, hold/ release, reweigh, recall, zero, KG/LB, tare, pretare and more. The exclusive Everlock[®] function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications. The 2900KL-AM includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions



Health o meter® Professional 2900KL-AM

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that the Health o meter[®] Professional Scales is the industryleading name in medical scales.

120+ LOCAL FIELD SPECIALISTS

100% CONNECTIVITY

#1 TRUSTED BRAND

100 YEARS OF SERVICE



weigheasier

TECHNICAL SPECIFICATIONS

Capacity 1000 lb / 454 kg

Resolution 0.2 lb / 0.1 kg

EMR connectivity Via USB, optional Pelstar[®] wireless technology

Display 1 ½" / 38 mm High-Contrast Color TFT-LCD Screen

Power Source 100-240V adapter (ADPT30) included, 6 D batteries not included

Functions

LB/KG conversion LB/KG lock out BMI Zero Tare Pre-tare Hold/release Reweigh Recall Auto zero Auto off Time/date Variable auto off time Audible/mute sound option Parent/child weighing function Everlock[®] - permanent unit lock

Platform Size ($w \times d \times h$)

29 ½" x 32 ½" x 2 ¾" (749 mm x 819 mm x 70 mm) **Ramp Size (w x d)** 6" x 33 ½" (159 mm x 864 mm)

Product Footprint Unfolded $(w \ge d \ge h)$

58" x 37" x 37" (1473 mm x 940 mm x 940)

Product Footprint Folded

(w x d x h) 58" x 3.5" x 37" (1473 mm x 101 mm x 940 mm)

Product Weight 91 lb / 41 kg

Warranty 2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
2900KL-AM	Antimicrobial wall-mounted fold-up wheelchair scale, LB & KG
2900KG-AM	Antimicrobial wall-mounted fold-up wheelchair scale, KG only
2900KL-AM-BT	2900KL-AM with built-in Pelstar® wireless technology
2900KG-AM-BT	2900KG-AM with built-in Pelstar® wireless technology
SS-2900KL	ScaleSurance® 2 year extended warranty for 2900 series*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer



Easy to Lift Platform with a Recessed Handle



To learn more, call **1-800-815-6615** or visit **homscales.com/2900KL-AM**

