

# THE ACUTE CARE LINE

Dedicated to bringing you great value-based scales engineered with a purpose, the Health o meter® Professional hospital line was developed to improve workflow while keeping patients and caregivers safe. The newest state-of-the-art products include a fold-up wheelchair scale, a heavy-duty highly-mobile platform scale, and a high-resolution neonatal tray scale, leading the industry in antimicrobial protection, connectivity, accuracy, reliability, and safety. Other Acute Care scales include a stainless steel wet diaper/lap sponge/organ scale with an easy-to-read backlight and a patient transfer scale that weighs immobile patients instantly, shortening door to needle time.

## Heavy Duty Mobile Antimicrobial Platform Scale (upgraded) **3105KL-AM**

**NEW  
AVAILABLE  
NOW!**

- Capacity of 1000 lb/454 kg
- EMR connectivity, BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare, Everlock®
- Optional electronic and mechanical height rod
- Improved robust structural design and low profile platform
- Options such as KG only and built-in wireless models are available



## High-Resolution Antimicrobial Digital Neonatal/ Pediatric Tray Scale (upgraded)

**2210KL-AM** (2 sided-tray)  
**2210KL4-AM** (4 sided-tray)

**NEW  
AVAILABLE  
NOW!**

- High 1 gram resolution with a capacity of 20,000g / 20kg / 45lb
- EMR connectivity, breast milk intake, auto-weight hold, hold/release, zero, tare, pre-tare, reweigh, recall, GM/KG/LB conversion and lockout, auto-zero, and auto-off
- Round-edged removable tray and backlit display
- Options such as KG only, GM only, LB only and built-in wireless models are available



## Antimicrobial Wall Mounted Fold Down Wheelchair Scale **2900KL-AM**

**NEW  
AVAILABLE  
NOW!**

- Capacity of 1000 lb/454 kg
- Large easy to lift platform conserves space
- EMR connectivity, BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare, Everlock®
- Safety features include anti-slip caution tape, raised edges, recessed handle, and no pinch points
- Options such as KG only and built-in wireless models are available



## Digital Wet Diaper/Lap Sponge/ Organ Scale

**3400KL / 3400KG** (Pan)  
**3401KL** (Tray)

- Stable, stainless steel base and removable stainless steel weighing pan or tray that can be autoclaved for sterilization
- Powered by a long-life rechargeable battery
- Backlit display that's easy to read at any angle and light condition



## Patient Transfer Scale **PTS-1000KL**

- A transfer board with a built-in medical scale that provides a weight measurement when transferring a patient
- Weighs immobile patients instantly and accurately
- It helps shorten door-to-needle times for time-critical patients
- All departments can benefit from the PTS, particularly Stroke Units, the Emergency Room, ICU, and Radiology.



To learn more, call **1-800-815-6615**  
or visit [homscales.com/acutecarescales](https://homscales.com/acutecarescales)

**Health o meter®**  
Professional