

# THE BTKIT-01

## WIRELESS UPGRADE KIT

The Health o meter® Professional BTKIT-01 allows a practitioner to add wireless connectivity to their existing scale. The kit includes a wireless module for installation into the scale's display head and a USB dongle for the scale to communicate with an external device, such as a vital signs monitor or PC.

The BTKIT-01 allows facilities to wirelessly transmit measurements directly into an EMR system to increase workflow efficiency and reduce the risk of error from manual data entry. This direct transmission conveniently eliminates the frustration of tangled wires and the tripping hazard of long cables and cords.

This kit is for facilities that want to connect to Welch Allyn® Connex® devices: Connex® Spot Monitor, Connex® Vital Signs Monitors, Connex® Integrated Wall Systems. This kit can also be used to connect to select EMR systems via a Windows® PC or Ceiba Healthcare's IoMT eConnect Box and Cloud platform. Transmitting data into an EMR system requires that the user's Windows® PC has one of the following systems installed: Allscripts TouchWorks® or Professional™, Midmark® IQmanager® or ChARM Health EHR.

Kit includes:

- USB Wireless Dongle (for use with Welch Allyn monitor)
- USB Cable 1 ft. (for use with the dongle and Welch Allyn monitor)
- Wireless Module (for installation into the scale's display head)

System Requirements When Connecting to a PC:

- CPU of 1.7GHz and above
- Minimum 512MB RAM
- Windows XP/Vista/7
- Available USB 2.0 Port
- Bluetooth® Capable or Bluetooth® Card\*

Communication Range:

- Up to approximately 328ft/100m



## TECHNICAL SPECIFICATIONS

### Warranty

2-year limited manufacturer warranty

### BTKIT-01 Compatible Scales\*\*

500 Series†  
1100 Series  
2101 Series  
2210KL-AM Series  
2210KL4-AM Series  
2400 Series  
2500 Series  
2600 Series  
2610 Series  
2650 Series  
2700 Series  
2900KL-AM Series  
3001KL-AM Series  
3105KL-AM Series

## ORDERING INFORMATION

MODEL	DESCRIPTION
BTKIT-01	Wireless Upgrade Kit

\*Bluetooth® is a registered trademark of the Bluetooth Special Interest Group. While Health o meter® Professional scales use proprietary technologies for reliably communicating with other devices, it can also be used with many Bluetooth® interfaces.

\*\*For assistance with compatible scales, contact Customer Service at 800-815-6615.

†BTKIT-01 is compatible with 500KL scales with a date code with the last two digits of 20 or higher, and 500KG scales with a date code with the last two digits of 21 or higher.

**Health o meter®**  
Professional

Health o meter® Professional BTKIT-01

weigh**easier**®

To learn more, call **1-800-815-6615**  
or visit **homscales.com/BTKIT-01**