

# THE BTKIT-01

## WIRELESS UPGRADE KIT FOR USE WITH WELCH ALLYN® MONITORS

The Health o meter® Professional BTKIT-01 allows a practitioner to add wireless connectivity to their existing scale. The kit includes a wireless module for installation into the scale's display head and a USB dongle for the scale to communicate with a Welch Allyn® Vital Signs Monitor.

The BTKIT-01 allows facilities to wirelessly transmit measurements directly to a Welch Allyn® Monitor to increase workflow efficiency and reduce the risk of error from manual data entry. This direct transmission conveniently eliminates the frustration of tangled wires and the tripping hazard of long cables and cords.

This kit is for facilities that want to connect to Welch Allyn® Connex® devices: Connex® Spot Monitor, Connex® Vital Signs Monitors, Connex® Integrated Wall Systems. From the Welch Allyn device, the data can be displayed, saved, and transferred into EMR systems.\*

Kit includes:

- USB Wireless Dongle (for use with Welch Allyn monitor)
- USB Cable 1 ft. (for use with the dongle and Welch Allyn monitor)
- Wireless Module (for installation into the scale's display head)

Communication Range:

- Up to approximately 328ft/100m

## TECHNICAL SPECIFICATIONS


### Warranty

2-year limited manufacturer warranty



### BTKIT-01 Compatible Scales\*\*

500 Series†  
502 Series  
597 Series  
599 Series  
600 Series  
752 Series  
1100 Series  
1110 Series  
2101 Series  
2210KL-AM Series  
2210KL4-AM Series  
2400 Series  
2500 Series  
2600 Series  
2610 Series  
2650 Series  
2700 Series  
2900KL-AM Series  
3001KL-AM Series  
3001KL-AMX Series  
3105KL-AM Series



## ORDERING INFORMATION

MODEL	DESCRIPTION
BTKIT-01	Wireless Upgrade Kit for Use with Welch Allyn® Monitors

\*Transfer to the EMR requires the EMR vendor to have an existing interface that accepts these parameters from Welch Allyn® devices.

\*\*BTKIT-01 is compatible with versions of these scales that have a date code with the last two digits of 19 or higher. For assistance with compatible scales, contact Customer Service at 800-815-6615.

†BTKIT-01 is compatible with 500KL scales with a date code with the last two digits of 20 or higher, and 500KG scales with a date code with the last two digits of 21 or higher.

**Health o meter®**  
Professional

Health o meter® Professional BTKIT-01

weigh**ea**sier®

To learn more, call **1-800-815-6615**  
or visit **homscales.com/BTKIT-01**