

THE BTKIT-01

WIRELESS UPGRADE KIT FOR USE WITH WELCH ALLYN® MONITORS

The Health o meter® Professional BTKIT-01 allows a practitioner to add wireless connectivity to their existing scale. The kit includes a wireless module for installation into the scale's display head and a USB dongle for the scale to communicate with a Welch Allyn® Vital Signs Monitor.

The BTKIT-01 allows facilities to wirelessly transmit measurements directly to a Welch Allyn® Monitor to increase workflow efficiency and reduce the risk of error from manual data entry. This direct transmission conveniently eliminates the frustration of tangled wires and the tripping hazard of long cables and cords.

This kit is for facilities that want to connect to Welch Allyn® Connex® devices: Connex® Spot Monitor, Connex® Vital Signs Monitors, Connex® Integrated Wall Systems. From the Welch Allyn device, the data can be displayed, saved, and transferred into EMR systems.*

Kit includes:

- USB Wireless Dongle (for use with Welch Allyn monitor)
- USB Cable 1 ft. (for use with the dongle and Welch Allyn monitor)
- Wireless Module (for installation into the scale's display head)

Communication Range:

- Up to approximately 328ft/100m

TECHNICAL SPECIFICATIONS


Warranty

2-year limited manufacturer warranty



BTKIT-01 Compatible Scales**

498 Series
500 Series
502 Series
522 Series
524KL
553 Series
565 Series
597 Series
599 Series
600 Series
752 Series
1100 Series
1110 Series
2000 Series
2101 Series
2210KL-AM Series
2210KL4-AM Series
2400 Series
2500 Series
2600 Series
2610 Series
2650 Series
2700 Series
2900KL-AM Series
3001KL-AM Series
3001KL-AMX Series
3105KL-AM Series



ORDERING INFORMATION

MODEL	DESCRIPTION
BTKIT-01	Wireless Upgrade Kit for Use with Welch Allyn® Monitors

*Transfer to the EMR requires the EMR vendor to have an existing interface that accepts these parameters from Welch Allyn® devices.

**List of compatible scales refers to new, current version scales. For legacy scales already in the field, contact Customer Service at 1-800-815-6615 to confirm compatibility.

Health o meter®
Professional

Health o meter® Professional BTKIT-01

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/BTKIT-01**