


THE 280HR

INFANT MEASURING MAT

The Health o meter® Professional 280HR is a portable Measuring Mat for infants and toddlers. The space-saving mat is soft, lightweight, foldable, and easy to store. The fixed headpiece and sliding foot positioner allow for precise measurement of infants or toddlers while lying down. The mat is easy-to-read with large font printed on both edges and measures in both centimeters and inches. It is also easy to use, with diagram instructions printed on the front of the mat. The 280HR can be conveniently cleaned with most disinfectants and cleaning wipes.



Lightweight & Foldable

Infant Measuring Mat can be rolled up for easy storage & portability



TECHNICAL SPECIFICATIONS

Measuring Range

3.9" – 43.3" / 10cm – 110cm

Graduation

1/4" / 5mm

Product Footprint (l x w)

55" x 12"

1397mm x 305mm

Rolled Up Measurement

4" x 12"

102mm x 305mm

Product Weight

.88 lb / 0.4 kg

Portability

Rolls up for easy transport

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
280HR	Infant Measuring Mat

Health o meter®
Professional

Health o meter® Professional 280HR

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/280HR**