

# THE BTKIT-02

## WIRELESS UPGRADE KIT

The Health o meter® Professional BTKIT-02 allows a practitioner to add wireless connectivity to their existing scale. The kit includes a wireless module for installation into the scale's display head and a USB dongle for the scale to communicate with an external device, such as a vital signs monitor or PC.

The BTKIT-02 allows facilities to wirelessly transmit measurements directly into an EMR system to increase workflow efficiency and reduce the risk of error from manual data entry. This direct transmission conveniently eliminates the frustration of tangled wires and the tripping hazard of long cables and cords.

This kit is for facilities that want to connect to Welch Allyn® Connex® devices: Connex® Spot Monitor, Connex® Vital Signs Monitors, Connex® Integrated Wall Systems. This kit can also be used to connect to select EMR systems via a Windows® PC or Ceiba Healthcare's IoMT eConnect Box and Cloud platform. Transmitting data into an EMR system requires that the user's Windows® PC has one of the following systems installed: Allscripts TouchWorks® or Professional™, Midmark® IQmanager® or ChARM Health EHR.

Kit includes:

- USB Wireless Dongle (for use with Welch Allyn monitor)
- USB Cable 1 ft. (for use with the dongle and Welch Allyn monitor)
- Wireless Module (for installation into the scale's display head)

System Requirements When Connecting to a PC:

- CPU of 1.7GHz and above
- Minimum 512MB RAM
- Windows XP/Vista/7
- Available USB 2.0 Port
- Bluetooth® Capable or Bluetooth® Card\*

Communication Range:

- Up to approximately 328ft/100m

## ORDERING INFORMATION

MODEL	DESCRIPTION
BTKIT-02	Wireless Upgrade Kit

\*Bluetooth® is a registered trademark of the Bluetooth Special Interest Group. While Health o meter® Professional scales use proprietary technologies for reliably communicating with other devices, it can also be used with many Bluetooth® interfaces.

\*\*For assistance with compatible scales, contact Customer Service at 800-815-6615.

†BTKIT-02 is compatible with versions of these scales that have a date code with the last two digits of 19 or higher. For assistance with compatible scales, contact Customer Service at 800-815-6615.

# Health o meter®

P r o f e s s i o n a l

Health o meter® Professional BTKIT-02

## TECHNICAL SPECIFICATIONS

### Warranty

2-year limited manufacturer warranty



### BTKIT-02† Compatible Scales\* \*

597KL/597KG

599KL/599KG

600KL/600KG

752KL/752KG



weigh**easier**®

To learn more, call **1-800-815-6615**  
or visit **homscales.com/BTKIT-02**