

THE BTKIT-100

WIRELESS UPGRADE KIT FOR USE WITH PRE-APPROVED DEVICES AND WINDOWS® / LINUX® APPLICATIONS (USB DONGLE NOT INCLUDED)

The Health o meter® Professional BTKIT-100 allows a practitioner to add wireless connectivity to their new or existing scale. The kit includes a wireless module for installation into the scale's display head for the scale to communicate with Windows® and Linux® devices such as a Windows® PC. The wireless module in the scale allows for direct pairing to pre-approved devices and Windows® / Linux® applications.

The BTKIT-100 allows facilities to wirelessly transmit measurements directly to a myriad of devices to increase workflow efficiency and reduce the risk of error from manual data entry. Once the data reaches the device, it can be conveniently displayed, reviewed, and smoothly transferred into EMR systems.

Kit includes:


- Wireless Module (for installation into the scale's display head)

System Requirements When Connecting to a PC/Device:

- Windows 10/11
- Bluetooth® Capable or Bluetooth® Card*

Communication Range:

- Up to approximately 328ft/100m



TECHNICAL SPECIFICATIONS

Warranty

2-year limited manufacturer warranty

BTKIT-100 Compatible Scales**

498 Series
500 Series
502 Series
522 Series
524KL
553 Series
565 Series
597 Series
599 Series
600 Series
752 Series
1100 Series
1110 Series
2000 Series
2101 Series
2210KL-AM Series
2210KL4-AM Series
2400 Series
2500 Series
2600 Series
2610 Series
2650 Series
2700 Series
2900KL-AM Series
3001KL-AM Series
3001KL-AMX Series
3105KL-AM Series

ORDERING INFORMATION

MODEL	DESCRIPTION
BTKIT-100	Wireless Upgrade Kit for Use with Pre-Approved Devices and Windows® / Linux® Applications, USB Dongle Not Included

*Bluetooth® is a registered trademark of the Bluetooth Special Interest Group. While Health o meter® Professional scales use proprietary technologies for reliably communicating with other devices, it can also be used with many Bluetooth® interfaces.

**BTKIT-100 is not compatible with Welch Allyn® Vital Signs Monitors or BT version scales. List of compatible scales refers to new, current version scales. For assistance with compatible scales, connectivity solutions, or legacy scales already in the field, contact Customer Service at 1-800-815-6615 to confirm compatibility.

Health o meter®
Professional

Health o meter® Professional BTKIT-100

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/BTKIT-100**