



### Advantages

- ✓ **Steel Frame & Powder Coated Fabrication**  
The classic physician scale is built with a robust steel frame and powder-coated column to prevent rust
- ✓ **Exclusive 1 Year Limited Warranty**  
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship
- ✓ **Height Rod Included**  
Includes a convenient height rod so clinicians can measure height and weight all in one place. The easy-to-use height rod provides height measurements from 11  $\frac{7}{8}$ " - 83.5" / 30cm - 212cm



### Specifications

- Capacity: 440 lb / 200 kg
- Graduation: 1/4 lb / 100 g
- Platform Size (w x d x h): 10.75" x 14.75" x 4" /  
273 mm x 375 mm x 101 mm, platform mat: non-slip black
- Fixed Beam, pounds and kilograms can be read on both sides
- Professional Dual Reading Height Rod: 11 7/8" -  
83.5" / 30cm - 212cm
- Optional wheels available (item # 69-00090)
- 1 Year Limited Warranty