

THE BTKIT-03

WIRELESS UPGRADE KIT FOR USE WITH PRE-APPROVED DEVICES AND WINDOWS®/ LINUX® PCS (USB DONGLE INCLUDED)

The Health o meter® Professional BTKIT-03 allows a practitioner to add wireless connectivity to their new or existing scale. The kit includes a wireless module for installation into the scale's display head. The kit also includes a wireless USB dongle that allows for connectivity to pre-approved devices and Windows®/ Linux® PCs. Pre-approved devices include Edan iM3 and MDPro 2500 Vital Signs Monitors.

The BTKIT-03 allows facilities to wirelessly transmit measurements directly to a myriad of devices to increase workflow efficiency and reduce the risk of error from manual data entry. Once the data reaches the device, it can be conveniently displayed, reviewed, and smoothly transferred into EMR systems.

Kit includes:


- USB Wireless Dongle (for use with the PC/device)
- USB Cable 1 ft. (for use with USB dongle)
- Wireless Module (for installation into the scale's display head)

System Requirements When Connecting to a PC/Device:

- Windows 10/11
- Available USB 2.0 Port (for use with included USB dongle)

Communication Range:

- Up to approximately 328ft/100m



TECHNICAL SPECIFICATIONS

Warranty

2-year limited manufacturer warranty

BTKIT-03 Compatible Scales*

498 Series
500 Series
502 Series
522 Series
524KL
553 Series
565 Series
597 Series
599 Series
600 Series
752 Series
1100 Series
1110 Series
2000 Series
2101 Series
2210KL-AM Series
2210KL4-AM Series
2400 Series
2500 Series
2595 Series
2600 Series
2610 Series
2620 Series
2650 Series
2700 Series
2900KL-AM Series
3001KL-AM Series
3001KL-AMX Series
3105KL-AM Series

ORDERING INFORMATION

MODEL	DESCRIPTION
BTKIT-03	Wireless Upgrade Kit for Use with Pre-Approved Devices and Windows®/ Linux® PCs, USB Dongle Included

*BTKIT-03 is not compatible with Welch Allyn® Vital Signs Monitors or BT version scales. List of compatible scales refers to new, current version scales. For assistance with compatible scales, connectivity solutions, or legacy scales already in the field, contact Customer Service at 1-800-815-6615 to confirm compatibility.

Health o meter®
Professional

Health o meter® Professional BTKIT-03

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/BTKIT-03**