

THE BTPRINTERKIT-01

WIRELESS THERMAL PRINTER KIT

The Health o meter® Professional BTPRINTERKIT-01 allows practitioners to add wireless connectivity to their existing scale that communicates directly with a thermal wireless printer. The kit includes a wireless module for installation into the scale's display head that is specifically paired to the included wireless printer.

Kit includes:

- Wireless Module (for installation into the scale's display head)
- Thermal printer
- Battery charger for printer
- Roll of BPA-free thermal printer paper
- AC to USB wall power adapter

Communication Range:

- Up to approximately 328ft/100m



TECHNICAL SPECIFICATIONS

Type of Paper

58 mm (2 1/4") Thermal Paper
(15240mm/50' Roll)

Print Width

48 mm / 1.9"

Printing Speed

90 mm/sec

Power Source

Rechargeable Lithium battery, 7V/1500mAh

Product Footprint

4" x 3" x 1 7/8"
(104 mm x 75 mm x 45 mm)

Product Weight

134 g / 1/3 lb

Warranty

2-year limited manufacturer warranty

BTPRINTERKIT-01 Compatible Scales*

500 Series†
502 Series
1100 Series
2101 Series
2210KL-AM Series
2210KL4-AM Series
2400 Series
2500 Series
2600 Series
2610 Series
2650 Series
2700 Series
2900KL-AM Series
3001KL-AM Series
3001KL-AMX Series
3105KL-AM Series

ORDERING INFORMATION

MODEL	DESCRIPTION
BTPRINTERKIT-01	Wireless Thermal Printer Kit

*Kit is not compatible with BT version scales already in the field. For assistance with compatible scales, contact Customer Service at 800-815-6615.

†BTPRINTERKIT-01 is compatible with 500KL scales with a date code with the last two digits of 20 or higher, and 500KG scales with a date code with the last two digits of 21 or higher.

Health o meter®
Professional

Health o meter® Professional BTPRINTERKIT-01

weigh**ea**sier®

To learn more, call **1-800-815-6615**
or visit **hom**scales.com**/BTPRINTERKIT-01**