

THE BTPRINTERKIT-01

WIRELESS THERMAL PRINTER KIT

The Health o meter® Professional BTPRINTERKIT-01 allows practitioners to add wireless connectivity to their existing scale that communicates directly with a thermal wireless printer. The kit includes a wireless module for installation into the scale's display head that is specifically paired to the included wireless printer.

Kit includes:

- Wireless Module (for installation into the scale's display head)
- Thermal printer
- Battery charger for printer
- Roll of BPA-free thermal printer paper
- AC to USB wall power adapter

Communication Range:

- Up to approximately 328ft/100m



TECHNICAL SPECIFICATIONS

Type of Paper

58 mm (2 1/4") Thermal Paper
(15240mm/50' Roll)

Print Width

48 mm / 1.9"

Printing Speed

90 mm/sec

Power Source

Rechargeable Lithium battery, 7V/1500mAh

Product Footprint

4" x 3" x 1 7/8"
(104 mm x 75 mm x 45 mm)

Product Weight

134 g / 1/3 lb

Warranty

2-year limited manufacturer warranty

BTPRINTERKIT-01 Compatible Scales*

498 Series
500 Series
502 Series
522 Series
524KL
553 Series
565 Series
1100 Series
1110 Series
2000 Series
2101 Series
2210KL-AM Series
2210KL4-AM Series
2400 Series
2500 Series
2600 Series
2610 Series
2650 Series
2700 Series
2900KL-AM Series
3001KL-AM Series
3001KL-AMX Series
3105KL-AM Series

weigh easier®

ORDERING INFORMATION

MODEL	DESCRIPTION
BTPRINTERKIT-01	Wireless Thermal Printer Kit

*BTPRINTERKIT-01 is not compatible with BT version scales already in the field. List of compatible scales refers to new, current version scales. For assistance with compatible scales, connectivity solutions, or legacy scales already in the field, contact Customer Service at 1-800-815-6615 to confirm compatibility.

Health o meter®
Professional

Health o meter® Professional BTPRINTERKIT-01

To learn more, call **1-800-815-6615**
or visit **hom scales.com/BTPRINTERKIT-01**