

# THE BTPRINTERKIT-01

## WIRELESS THERMAL PRINTER KIT

The Health o meter® Professional BTPRINTERKIT-01 allows practitioners to add wireless connectivity to their existing scale that communicates directly with a thermal wireless printer. The kit includes a wireless module for installation into the scale's display head that is specifically paired to the included wireless printer.

Kit includes:

- Wireless Module  
(for installation into the scale's display head)
- Thermal printer
- Battery charger for printer
- Roll of BPA-free thermal printer paper
- AC to USB wall power adapter

Communication Range:

- Up to approximately 328ft/100m



### TECHNICAL SPECIFICATIONS

**Type of Paper**

58 mm (2 ¼") Thermal Paper  
(15240mm/50' Roll)

**Print Width**

48 mm / 1.9"

**Printing Speed**

90 mm/sec

**Power Source**

Rechargeable Lithium battery, 7V/1500mAh

**Product Footprint**

4" x 3" x 1 7/8"  
(104 mm x 75 mm x 45 mm)

**Product Weight**

134 g / 1/3 lb

**Warranty**

2-year limited manufacturer warranty

**BTPRINTERKIT-01 Compatible Scales\***

- 498 Series
- 500 Series
- 502 Series
- 522 Series
- 524KL
- 553 Series
- 565 Series
- 1100 Series
- 1110 Series
- 2000 Series
- 2101 Series
- 2210KL-AM Series
- 2210KL4-AM Series
- 2400 Series
- 2500 Series
- 2600 Series
- 2610 Series
- 2620 Series
- 2650 Series
- 2700 Series
- 2900KL-AM Series
- 3001KL-AM Series
- 3001KL-AMX Series
- 3105KL-AM Series

### ORDERING INFORMATION

MODEL	DESCRIPTION
BTPRINTERKIT-01	Wireless Thermal Printer Kit

\*BTPRINTERKIT-01 is not compatible with BT version scales already in the field. List of compatible scales refers to new, current version scales. For assistance with compatible scales, connectivity solutions, or legacy scales already in the field, contact Customer Service at 1-800-815-6615 to confirm compatibility.