

# THE BTPRINTERKIT-01

## WIRELESS THERMAL PRINTER KIT

The Health o meter® Professional BTPRINTERKIT-01 allows practitioners to add wireless connectivity to their existing scale that communicates directly with a thermal wireless printer. The kit includes a wireless module for installation into the scale's display head that is specifically paired to the included wireless printer.

Kit includes:

- Wireless Module  
(for installation into the scale's display head)
- Thermal printer
- Battery charger for printer
- Roll of BPA-free thermal printer paper
- AC to USB wall power adapter

Communication Range:

- Up to approximately 328ft/100m



## TECHNICAL SPECIFICATIONS

### Type of Paper

58 mm (2 ¼") Thermal Paper  
(15240mm/50' Roll)

### Print Width

48 mm / 1.9"

### Printing Speed

90 mm/sec

### Power Source

Rechargeable Lithium battery, 7V/1500mAh

### Product Footprint

4" x 3" x 1 7/8"  
(104 mm x 75 mm x 45 mm)



### Product Weight

134 g / 1/3 lb

### Warranty

2-year limited manufacturer warranty

### BTPRINTERKIT-01 Compatible Scales\*

498 Series  
500 Series  
502 Series  
522 Series  
524KL  
553 Series  
565 Series  
1100 Series  
1110 Series  
2000 Series  
2101 Series  
2210KL-AM Series  
2210KL4-AM Series  
2400 Series  
2500 Series  
2600 Series  
2610 Series  
2620 Series  
2650 Series  
2700 Series  
2900KL-AM Series  
2920KL-AM Series  
3001KL-AM Series  
3001KL-AMX Series  
3105KL-AM Series

## ORDERING INFORMATION

MODEL	DESCRIPTION
BTPRINTERKIT-01	Wireless Thermal Printer Kit

\*BTPRINTERKIT-01 is not compatible with BT version scales already in the field. List of compatible scales refers to new, current version scales. For assistance with compatible scales, connectivity solutions, or legacy scales already in the field, contact Customer Service at 1-800-815-6615 to confirm compatibility.

**Health o meter®**  
Professional

Health o meter® Professional BTPRINTERKIT-01

weigh**easier**®

To learn more, call **1-800-815-6615**  
or visit **homscales.com/BTPRINTERKIT-01**