


Introducing the Patient Transfer Scale

- A transfer board with a built-in medical scale provides a weight measurement when transferring a patient
- A **kilograms only Patient Transfer Scale** offers a safe accurate solution for facilities standardizing on metric-only measurements
- Can be utilized as part of an existing patient transfer process



Weigh critical or immobile patients instantly and accurately

*Making Healthcare Professionals' Workflow Weigh**asier***

Health o meter
Professional



UK PATENT
Patent number 2543835



US PATENT PENDING



EU PATENT PENDING



CANADA PATENT PENDING



HONG KONG PATENT PENDING

The Patient Transfer Scale: Every Minute Counts



Key Features / Specifications:

- **Large, clear digital weight display**
- **Portable**
- **Wipe clean surface**
- **Capacity:** 250kg
- **Resolution:** 0.5kg
- **Power:** Internal rechargeable battery
- **Battery life:** 50 hours
- **Weight:** 25lb / 11.4 kg
- **Dimensions:**
- 71" x 27 1/2" x 1" / 1805mm x 700mm x 30 mm
- **Model Number:**
PTS-1000KG
PTS-1000KG-C

- When managing a patient with a critical condition like a stroke, **Every Minute Counts**. It is imperative to get an accurate weight quickly when patient weights are crucial and where dosages of critical medications are based on kilogram measurements. The Patient Transfer Scale (PTS) does just that.
- Slide Sheets should be utilized with the PTS when transferring a patient to minimize the risk of injuries and ensure safety for both patients and healthcare providers. To learn about SlideMate™ Slide Tubes visit [https:// bridgehcsa.com/safe-patient-handling/slidemate-single-patient-use-slide-tubes/](https://bridgehcsa.com/safe-patient-handling/slidemate-single-patient-use-slide-tubes/)
- Learn more about how the PTS works Visit <https://www.homscales.com/products/patient-transfer-scale>
- Get to know the inventor of the PTS Gillian Taylor. Read her success stories.
 Twitter: @gilliantaym999
 Instagram: @gilliantaylorpts
- The Patient Transfer Scale was developed and invented by an Emergency Care nurse, with the final design based on feedback from healthcare professionals . . . **For Nurses by Nurses.**
- All departments can benefit from the Patient Transfer Scale, particularly Stroke Units, the Emergency Room, ICU, and Radiology.

Scott Gottman
Director of Marketing and Product Development
(708) 297-6582
sgottman@homscales.com

Health o meter® Professional Scales
(800) 815-6615
www.homscales.com