

## THE PTS-1000KG

# PATIENT TRANSFER SCALE, KG ONLY

The PTS-1000KG is a transfer board with a built-in medical scale that provides a weight measurement when transferring a patient. A kilograms only Patient Transfer Scale offers a safe solution for facilities standardizing on metric only measurements. All departments can benefit from the Patient Transfer Scale, particularly Stroke Units, the Emergency Room, ICU, and Radiology.

Designed by Clinicians for Clinicians the Patient Transfer Scale weighs critical or immobile patients instantly and accurately helping shorten door-to-needle times. A kilograms only Patient Transfer Scale is especially useful for Stroke Units where quick and accurate patient weights are crucial and where dosages of critical medications are based on kilogram measurements. It offers a superior alternative to other weighing equipment, such as in-bed scales, that are often not calibrated or inaccurate.

The scale is light and portable and can be utilized as part of an existing patient transfer process when a weight is needed, and will not disrupt the patient flow process.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that Health o meter® Professional Scales is the industry-leading name in medical scales.

**120+** LOCAL  
FIELD SPECIALISTS

**#1** TRUSTED  
BRAND

**100** YEARS  
OF SERVICE



**Health o meter®**  
Professional

Health o meter® Professional PTS-1000KG

weigh**easier**®

# TECHNICAL SPECIFICATIONS

### Capacity

250 kg

### Resolution

0.5 kg

### Display

3" x 1" / 76 mm x 25 mm Backlight LCD

### Power Source

3.7V rechargeable battery

### Functions

Zero  
Hold/release  
Auto off

### Weighing Mat Size (w x d)

71" x 21"  
(1805mm x 533mm)

### Product Footprint (w x d x h)

71" x 27 1/2" x 1"  
(1805mm x 700mm x 30mm)

### Product Weight

25 lb / 11.4 kg

### Measuring Tape

Height guide shown on mat 100cm – 180cm  
shown at 5 cm

### Warranty

2-year limited manufacturer warranty

# ORDERING INFORMATION

MODEL	DESCRIPTION
PTS-1000KG	Patient Transfer Scale, KG only
PTS-1000KG-C	Patient Transfer Scale, KG only, Canada only
SS-PTS-1000KL*	ScaleSurance® 2 year extended warranty for Patient Transfer Scale



Always use a Slide Sheet  
with the Patient Transfer Scale!



Check out the Single Patient Use  
and Reusable Slide Sheets.

To learn more, call **1-800-815-6615**  
or visit **hom scales.com/PTS-1000KG**



\*Not available in all countries.

**Health o meter®**  
Professional