

THE ACUTE CARE LINE

The Health o meter® Professional Acute Care value-based scales are designed to improve workflow and ensure safety for patients and caregivers. They offer a range of state-of-the-art products, such as a fold-up wheelchair scale, heavy-duty platform scale, and high-resolution neonatal scale, which are known for their antimicrobial protection, accuracy, reliability, and safety features.


Heavy Duty Mobile Antimicrobial Platform Scale (upgraded) **3105KL-AM**

- Capacity of 1000 lb/454 kg
- EMR connectivity, BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare, Everlock®
- Optional electronic and mechanical height rod
- Improved robust structural design and low profile platform
- Options such as KG only and built-in wireless models are available



High-Resolution Antimicrobial Digital Neonatal/ Pediatric Tray Scale (upgraded) **2210KL-AM** (2 sided-tray) **2210KL4-AM** (4 sided-tray) **2210CART-AC** (acute care cart)

- High 1 gram resolution with a capacity of 20,000g / 20kg / 45lb
- EMR connectivity, breast milk intake, auto-weight hold, hold/release, zero, tare, pre-tare, reweigh, recall, GM/KG/LB conversion and lockout, auto-zero, and auto-off
- Round-edged removable tray and backlit display
- Optional electronic height rod
- Options such as KG only, GM only, LB only, and built-in wireless models are available



Antimicrobial Wall Mounted Fold Down Wheelchair Scale **2900KL-AM**

- Capacity of 1000 lb/454 kg
- Large easy to lift platform conserves space
- EMR connectivity, BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare, Everlock®
- Safety features include anti-slip caution tape, raised edges, recessed handle, and no pinch points
- Options such as KG only and built-in wireless models are available



Digital Wet Diaper/Lap Sponge/ Organ Scale **3400KL / 3400KG** (Pan) **3401KL** (Tray)

- Stable, stainless steel base and removable stainless steel weighing pan or tray that can be autoclaved for sterilization
- Powered by a long-life rechargeable battery
- Backlit display that's easy to read at any angle and light condition



Patient Transfer Scale

PTS-1000KL/ PTS-1000KG (New KG only model)

- A transfer board with a built-in medical scale that provides a weight measurement when transferring a patient
- Weighs immobile patients instantly and accurately
- It helps shorten door-to-needle times for time-critical patients
- All departments can benefit from the PTS, particularly Stroke Units, the Emergency Room, ICU, and Radiology
- KG only model is a safe solution for facilities standardizing on metric only measurements for critical medication dosages based on kilograms



To learn more, call **1-800-815-6615**
or visit **homscales.com/acutecarescales**

Health o meter®
Professional