



Advantages



Moveable Arm & Foot Rest

The flip-arm and foot rest allows for easier patient access and patient comfort



Motion Detection

Scale has the ability to weigh accurately despite movement of the subject on the scale



Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, Zero Out/Tare, Auto Off, Reweigh) and more!



| Model | Shipping Dims. (L x W x H) / Cube | Weight (LB) | UPC Code |
|--------|-------------------------------------|-------------|--------------|
| 2595KL | 32.625" x 24.875" x 41.875" / 19.67 | 100.25 lb | 892076002626 |



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Seat Size: 18 1/4"(w) x 15 3/4" (d) x / 464 mm (w) x 400 mm (d)
- Display: 1 1/8" LCD
- Wheels: 2
- Connectivity: via optional Pelstar® wireless technology
- Power Source: 110-240V Adapter (included), 6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out /Tare, Auto Off, Variable auto-off Time, Enabling unit view with unit of measure locked, Hold/Release, Reweigh, Pre-Tare, Recall, Everlock®
- 2 Year Limited Warranty

| Model | Description |
|-----------|---|
| 2595KL | Digital Chair Scale - pounds & kilograms |
| 2595KG | Digital Chair Scale - kilograms only |
| ADPT30 | Optional adapter for 2595 series of scales |
| PLUGSET30 | Plug adapters for use outside North America |

Scale also includes connectivity capabilities. To learn about connectivity options, visit hom scales.com/ConnectivitySolutions