

THE 2600ADAKIT

ADA COMPLIANCE KIT FOR 2600 SERIES OF SCALES

The Health o meter® Professional 2600ADAKIT is an optional ADA Compliance Kit for the 2600 series of scales. Certified for Americans with Disabilities Act (ADA) compliance, the kit includes new ramps, new feet, and platform and ramp edge guards with anti-slip caution tape on the ramps to provide enhanced safety. The 2600ADAKIT is easy to install on any existing 2600 series scale, making it an ideal choice for doctor's offices, hospitals, clinics, and acute or long-term care facilities wanting to comply with the latest MDE standards.



See Ordering Table Below for White Scales



2600ADAKIT for 2600 Series of Scales, Designed to Meet the US Access Boards Standards: Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which is enforced by the Department of Justice.

- 32" wide x 40" deep platform
- 2" edge protection on platform
- 2" edge protection on ramps
- 1:8 running slope

Reviewed for ADA compliance and Platinum certified by Corada



REVIEWED FOR ADA COMPLIANCE
CERTIFICATION NO. 4477888

ORDERING INFORMATION

MODEL	DESCRIPTION
2600ADAKIT	ADA Compliance Kit For Gray 2600 Series Scales
2600ADAKIT-WHITE	ADA Compliance Kit For White 2600 Series Scales

*2600ADAKIT is compatible with 2600KL series or 2610 series scales that have date codes ending in 14 or higher. For assistance with compatible scales, contact Customer Service at 1-800-815-6615.

Health o meter®
Professional

Health o meter® Professional 2600ADAKIT

TECHNICAL SPECIFICATIONS

Edge Guard Size

40" x 1/8" x 2"
(1016 mm x 3 mm x 51 mm)

Ramp Size

32 1/4" x 16"
(819 mm x 406 mm)

Platform Size with 2600KL

40" x 32 1/4" x 1 7/8"
(1016 mm x 819 mm x 48 mm)

Product Footprint with 2600KL

77 1/4" x 43 7/8" x 48 3/4"
(1962 mm x 1114 mm x 1238 mm)

Product Footprint with 2600KL (Folded)

77 1/4" x 9 1/2" x 53"
(1962 mm x 241 mm x 1346 mm)

Product Weight with 2600KL

126 lb / 57 kg

Warranty

2-year limited manufacturer warranty

2600ADAKIT Compatible Scales*

2600KL
2600KG
2600KL-BT
2600KG-BT
2610KL
2610KG
2610KL-BT
2610KG-BT

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/2600ADAKIT**