

weigh**easier**®



Health o meter®

P r o f e s s i o n a l

Product Catalog



why choose us?

weigh**easier**®

What sets us apart is that everything we do - from product design to technical support - is done with our customers in mind. Each decision and each response is formed by first asking **What will be weigh**easier** for our customer?** – whether the customer is a physician, healthcare system, distributor, or sales team.

#1 Medical Scale Brand in North America

Health o meter® has been a trusted brand since 1919. Today, Health o meter® Professional Scales is the U.S. market leader and the fastest growing brand in Medical Scales.* We offer:

- High quality products that serve the entire continuum of medical care
- The industry's largest network of trained and qualified field representatives
- Direct toll-free number for Technical Support with customer advocates based in the USA
- Fast shipments and industry-leading fill rates
- Standardization programs
- No-cost product evaluations

100% Focus on Healthcare

Because all we do is support the healthcare market, we understand the impact of healthcare reform, Value-Based Purchasing, and other industry-specific initiatives. We understand that providers are overwhelmed with patient care and are focused on operating efficiently and with an eye on costs. We concern ourselves with patient safety, supply chain efficiency, and technology standardization because our customers are concerned about those things. We design our products to help patients and the providers and institutions that care for them.

Products Designed in the USA

Our US based engineers and product specialists design high-value, high-quality devices for users around the globe.

Rated #1 by Medical Distributors & Healthcare Professionals!

*As reported by Clarivate's US Medical Supply Distribution data.

table of contents



Safe Patient Handling
& Pressure Injury
Prevention Products

02

ScaleSurance® Extended
Warranty

03

Connectivity Solutions

04-55

Products

- ▲ Stand-On
- ▲ Body Composition Scale
- ▲ Pediatric
- ▲ Wheelchair / Chair / Multi-Use
- ▲ Specialty
- ▲ Professional Home Care
- ▲ Height Measuring
- ▲ Accessories
- ▲ Automated Calibrator

56

Glossary

Symbols

57

Global Support

ISO 9001:2015
ISO 13485:2016
CERTIFIED



Extended Warranty Program

Health o meter® Professional Scales offers an extended warranty that protects a facility's scale investment. This program extends the standard Limited Warranty that is available with select Health o meter® Professional products.

Benefits of an Extended Warranty:

- Peace of mind against the financial hardship of unexpected equipment replacement costs
- Quick product replacement – including free outbound shipping
- Hassle-free coverage provided by the manufacturer, not a third party warranty company
- Immediate, live interaction with product specialists for troubleshooting, diagnostics, or EMR connectivity assistance

A fraction of the cost of buying a new scale!

Additional Details:

- Can be purchased with a new scale or for a facility's existing scale before the current warranty has expired.
- Available on Health o meter® Professional digital scales and digital height rods that have a 2 year limited warranty and G6 Body Comp Scales

**For more information visit
www.homscales.com/scalesurance**

Not available in all countries

**To purchase an extended warranty or for more information,
contact your medical supply distributor today!**

Connectivity Solutions

Health o meter® Professional Scales provide seamless connectivity solutions for transmitting weight, height, and Body Mass Index (BMI) measurements to vital signs devices, PCs and EMR systems. This automatic data transfer minimizes transcription errors and is the fastest quickest, easiest way to record patient information. Device connectivity enhances healthcare facility workflow efficiency, patient safety, and EMR compliance.

Our open-architecture, non-propriety connectivity solutions are designed for universal integration with any device or EMR company that would like to establish direct interfaces between our scales and their systems. We will work with any EMR company at the request of our customers to facilitate direct interfaces with their systems. Our compatibility extends to various systems, including:



The following solutions are available for connecting Health o meter® Professional digital scales with vital signs devices, PCs, and EMR systems:

"BT" Model Scales with Built-in Pelstar® Wireless Technology

These scales include wireless connectivity and a wireless USB Dongle for facilities that want to transmit data to Welch Allyn® Connex® devices: Connex® Spot Monitor, Connex® Vital Signs Monitors, and Connex® Integrated Wall Systems.

BTKIT-01

Wireless upgrade kit is for Health o meter® Professional compatible digital scales and allows practitioners to add wireless connectivity to their new or existing scale. This wireless upgrade kit is for facilities that want to connect to Welch Allyn® Connex® devices: Connex® Spot Monitor, Connex® Vital Signs Monitors, and Connex® Integrated Wall Systems.

BTKIT-03

Wireless upgrade kit is for Health o meter® Professional compatible digital scales and allows practitioners to add wireless connectivity to their new or existing scale. This wireless upgrade kit is for facilities that want to connect to pre-approved devices and Windows®/ Linux® PCs. Pre-approved devices include Edan iM3, MDPro 2500 Vital Signs Monitors, and the Schiller DS20 Diagnostic Station.

BTKIT-100

Wireless upgrade kit is for Health o meter® Professional compatible digital scales and allows practitioners to add wireless connectivity to their new or existing scale. This wireless upgrade kit is for facilities that want to connect to pre-approved devices and Windows®/ Linux® applications. Pre-approved devices include Ceiba Healthcare's IoMT eConnect Box which can interface to several EMR Systems including Epic® and Cerner®.

For more information visit
www.homscales.com/ConnectivitySolutions



Stand-On



**Optional 500HB
Handlebar
Accessory can
be mounted
to your new or
existing 500KL,
500KG, 500KL-BT,
and 500KG-BT**

Power Source: 6 AA batteries included,
100-240V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out,
BMI, Zero, Tare, Pre-tare, Hold / Release, Auto Zero,
Auto Off, Reweigh, Recall, Parent / Child Weighing
Function, Variable Auto Off Time, Audible / Mute
Sound Option, Enabling Unit View with Unit of
Measure Locked, Everlock®

Wheels: 2 wheels,
8" wide
wheelbase

Connectivity: optional
Pelstar® wireless
technology

★ TOP SELLER

500KL



Videos & More!

EYE-LEVEL DIGITAL SCALE

The 500KL is the #1 selling physician scale in North America. With a built-in height rod, the scale enhances workflow by measuring weight and height and calculating Body Mass Index without moving the patient. Optional handlebars provide stability to keep both patients and caregivers safe.

Built-in Height Rod: The easy-to-use height rod provides height measurements from 30" - 84 1/4" / 76 cm - 214 cm.

Calculates BMI: Scale functions include the ability to calculate the patient's BMI.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Dent/Rust Proof Low Profile Platform: Non-skid platform requires only minimal step up, offering the patient a comfortable and stable weighing experience.



Capacity	Resolution	Platform Size (w x d x h)
660 lb / 300 kg	0.2 lb / 0.1 kg	13 3/4" x 16 1/2" x 2 3/8" 349 mm x 419 mm x 60 mm

Height Rod	Height Rod Graduation	Display
30" - 84 1/4" 76 cm - 214 cm	1/8" / 1 mm	1" / 25 mm LCD

Options & Accessories:

Model	Description
500KG	Scale in kilograms only
500KLHB	500KL with handlebars included
500KLAD	500KL with adapter included
500KGAD	500KG with adapter included
500KL-BT	500KL with built-in Pelstar® wireless technology
500KG-BT	500KG with built-in Pelstar® wireless technology
500HB	Optional handlebar accessory for 500 series scales
SS-500KL	2 year extended warranty
ADPT31	Adapter for 500 series scales
BTPRINTERKIT-01	Wireless thermal printer kit



Optional 500HB Handlebar Accessory can be mounted to your 502KL, 502KG, 502KL-BT, and 502KG-BT



Power Source: 6 AA batteries included, 100-240V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Auto Zero, Auto Off, Reweigh, Recall, Parent / Child Weighing Function, Variable Auto Off Time, Audible / Mute Sound Option, Enabling Unit View with Unit of Measure Locked, Everlock®

Wheels: 2 wheels, 8" wide wheelbase

Connectivity: optional Pelstar® wireless technology

502KL

PRIMARY CARE DIGITAL SCALE
WITH DIGITAL HEIGHT ROD



Videos & More!

The newest member of the #1 selling physician scale in North America. With a market-leading built-in digital height rod, the scale enhances workflow by measuring weight and height and calculating Body Mass Index without moving the patient. Optional handlebars provide stability to keep both patients and caregivers safe.

Built-in Digital Height Rod: The easy-to-use height rod provides a full measuring range of 11 ¾" – 79 ½" / 30cm – 202 cm and a digital range of 3ft 11.2" – 6ft 7.5" / 120 cm – 202 cm.

Calculates BMI: Scale functions include the ability to calculate the patient's BMI.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Dent/Rust Proof Low Profile Platform: Non-skid platform requires only minimal step up, offering the patient a comfortable and stable weighing experience.



Capacity	Resolution	Platform Size (w x d x h)
660 lb / 300 kg	0.2 lb / 0.1 kg	13 ¾" x 16 ½" x 2 ¾" 349 mm x 419 mm x 60 mm

Digital Height Rod	Mechanical Height Rod	Display
3ft 11.2" – 6ft 7.5" 120 cm – 202 cm Resolution: 1 mm / 0.1"	11 ¾" – 79 ½" 30 cm – 202 cm Graduation: 1/16" / 1 mm	1" / 25 mm LCD

Options & Accessories:

Model	Description
502KG	Scale in kilograms only
502KL-BT	502KL with built-in Pelstar® wireless technology
502KG-BT	502KG with built-in Pelstar® wireless technology
500HB	Optional handlebar accessory for 502 series scales
SS-502KL	2 year extended warranty
ADPT31	Adapter for 502 series scales
BTPRINTERKIT-01	Wireless thermal printer kit

499KL

WAIST-HIGH DIGITAL SCALE



**Optional
499HB Handlebar
Accessory!**

The 499KL has a waist-high easy to read display and the scale is an attractive, robust option for use in medical offices. With an optional height rod and optional handlebars, the 499KL is an economical choice for almost any environment.

Capacity: 460 lb / 209 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

13 3/4" x 16 1/2" x 2 3/8"

349 mm x 419 mm x 60 mm

Display: 1" / 25 mm LCD

Wheels: 2 wheels, 8" wide wheelbase

Power Source: 6 AA batteries

included, 100-240V adapter

(not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function, Everlock®

Options & Accessories:

Model	Description
499KG	Scale in kilograms only
499KLHR	499KL with height rod included
499KGHR	499KG with height rod included
499KLAD	499KL with adapter included
499KGAD	499KG with adapter included
499KLHRAD	499KL with height rod and adapter included
499KGHRAD	499KG with height rod and adapter included
499HB	Optional handlebar accessory for 499 series scales
SS-499KL	2 year extended warranty
499KLROD	Height rod for 499 series scales
ADPT31	Adapter for 499 series scales



599KL

HEAVY DUTY WAIST-HIGH DIGITAL SCALE



Great for doctor offices and clinics with minimal floor space, the 599KL has a high capacity with waist-high display that can swivel 350°. For providers wishing to measure height, an optional height rod is available.

Capacity: 660 lb / 300 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

14 1/4" x 14 1/4" x 2 5/8"

362 mm x 362 mm x 67 mm

Display: 1" / 25 mm LCD, 350° Swivel

Connectivity: USB, optional Pelstar wireless technology

Wheels: 2

Power Source:

100-240V adapter included,

6 C batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Options & Accessories:

Model	Description
599KG	Scale in kilograms only
599KLHR	599KL with height rod included
599KGHR	599KG with height rod included
SS-599KL	2 year extended warranty
499KLROD	Height rod for 599KL
ADPT50	Adapter included, item # for replacement or spare



500KG



ING YOUR WEIGHT
OR YOU?

ward to keep
weight you've lost?

Waist matters

col
know
prev
pres



600KL

HEAVY DUTY EYE-LEVEL DIGITAL SCALE
WITH DIGITAL HEIGHT ROD

Improve workflow with this top-of-the-line physician office scale. The built-in digital height rod eliminates the need to manually enter height for computing BMI. The 350° swivel head allows for line-of-sight viewing and patient privacy.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Automatically Calculates Weight, Height and BMI: Improve efficiency and accuracy of vital sign measurements.

Built-in Digital Height Rod: The easy-to-use, state of the art digital height rod provides measurements from 23 5/8" - 88" / 60 cm - 222 cm.

Patented Design: The 600KL's patented design provides a user-friendly and intuitive interface for measuring patient weight, height and BMI with a convenient rotating display head.



Power Source: 100-240V adapter included, 6 C batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out , BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Wheels: 2 **Connectivity:** USB, optional Pelstar® wireless technology

Capacity

660 lb / 300 kg

Resolution

0.2 lb / 0.1 kg

Platform Size (w x d x h)

14 1/4" x 14 1/4" x 2 5/8"
362 mm x 362 mm x 67 mm

Height Rod

23 5/8" - 88"
60 cm - 222 cm

Height Rod Resolution

0.1" / 0.1 cm

Display

1" / 25 mm LCD
350° Swivel

Options & Accessories:

Model	Description
600KG	Scale in kilograms only
600KL-BT	600KL with built-in Pelstar® wireless technology
600KG-BT	600KG with built-in Pelstar® wireless technology
ADPT50	Adapter included, item # for replacement or spare
SS-600KL	2 year extended warranty

498KL

REMOTE DISPLAY DIGITAL SCALE



The 498KL is a versatile choice for almost any clinical environment. Its portability and remote display head allow the scale to be moved from one location to another.

Capacity: 660 lb / 300 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

13 3/4" x 16 1/2" x 2 3/8"

349 mm x 419 mm x 60 mm

Display: 1" / 25 mm LCD

Connectivity: via optional Pelstar® wireless technology

Power Source: 6 AA batteries included, 100-240V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Reweigh, Auto Zero, Auto Off, Parent / Child Weighing Function, Variable Auto Off Time, Audible / Mute Sound Option, Enabling Unit View with Unit of Measure Locked, Everlock®

Options & Accessories:

Model	Description
498KG	Scale in kilograms only
498KLAD	498KL with adapter included
SS-498KL	2 year extended warranty
ADPT31	Adapter for 498KL
64771	Carrying case for 498KL
BTprinterKIT-01	Wireless thermal printer kit



Optional Height Rods, see page 48

752KL

HEAVY DUTY REMOTE DISPLAY DIGITAL SCALE



With a high 660 lb capacity, the 752KL is a great supplement to accommodate ambulatory bariatric patients. The remote display and portable base make it convenient to move the scale anywhere it is needed.

Capacity: 660 lb / 300 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

14 1/4" x 14 1/4" x 2 5/8"

362 mm x 362 mm x 67 mm

Display: 1" / 25 mm LCD

Connectivity: USB, optional Pelstar® wireless technology

Power Source: 100-240V adapter included, 6 C-cell batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off, Parent / Child Weighing Function

Options & Accessories:

Model	Description
752KG	Scale in kilograms only
SS-752KL	2 year extended warranty
ADPT50	Adapter included, item # for replacement or spare
64771	Carrying case for 752KL



Optional Height Rods, see page 48

★ TOP SELLER

349KLX



DIGITAL FLOOR SCALE

The 349KLX is a top selling scale that can be used in a host of locations, from MD offices to clinics to home care. Users can place the remote display head where it is easiest for them to read the weight. The scale's small footprint also makes it a great choice for weight loss clubs and home use.

Low Profile Platform: Low platform height offers patient a comfortable and stable weighing experience.

Remote Display: The remote display provides versatility with the ability to be mounted on a wall or placed on a table.

Power Source: 6 AA batteries included, 120V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, Zero, Tare, Hold / Release, Reweigh, Auto Zero, Auto Off, Variable auto off time, Parent / Child Weighing Function, Enabling Unit View with Unit of Measure Locked, Enable / Disable Backlight Display with 4 Levels of Brightness, Everlock®

Capacity

400 lb / 180 kg

Resolution

0.2 lb / 0.1 kg

Platform Size (w x d x h)

12 ½" x 12" x 1 ⅞"
318 mm x 305 mm x 48 mm

Display

1" / 25 mm LCD

Options & Accessories:

Model	Description
349KGX	Scale in kilograms only
349KLXAD	349KLX with adapter included
SS-349KLX	2 year extended warranty
ADPT40	Adapter for 349KLX
64771	Carrying case for 349KLX

Optional Height Rods, see page 46

3001KL-AM



★ TOP SELLER

3001KL-AM & 3001KL-AMX

Optional Digital and
Mechanical Height
Rods Available



ANTIMICROBIAL DIGITAL PLATFORM SCALE

Incredibly robust, these platform scales feature feet and wheels specially designed to withstand frequent use and regular movement. Antimicrobial coating on the handrails, transport handle, and keypad help reduce the spread of harmful microbes. The 3001KL-AMX offers extended handrails that provide extra support and help patients maintain balance when stepping onto the platform.

Infection Prevention: Silver ionic Antimicrobial Powder Coating suppresses the growth of bacteria and helps protect against cross-contamination.

Fall Prevention: "Live" handrails are advanced technology that allow the patient to hold on for stability while the scale accurately calculates weight.


Low Profile Platform: Requires only minimal step up, offering the patient a comfortable and stable weighing experience.

Durability: Convenient transport handle allows the user to easily move the scale to where it is needed in the facility. Rugged wheels allow for smooth and easy portability.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.



3001KL-AMX
ANTIMICROBIAL DIGITAL
PLATFORM SCALE WITH
EXTENDED HANDRAILS



Power Source: 100-240V
adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion,
LB / KG Lock Out, BMI, Zero, Tare, Pre-tare,
Hold / Release, Recall, Reweigh, Auto Zero,
Auto Off, Time / Date, Variable Auto Off
Time, Audible / Mute Sound Option, Parent /
Child Weighing Function, Everlock®

Wheels: 2 **Connectivity:** USB, optional
Pelstar® wireless
technology

Capacity	Resolution	Platform Size (w x d x h)
1000 lb / 454 kg	0.2 lb / 0.1 kg	22" x 15 3/4" x 2 5/8" 559 mm x 400 mm x 67 mm

Display

Digit Size: 3/8" / 10 mm,
Screen Size: 2 3/4" x 1 1/2" High-Contrast Color TFT-LCD Screen

Options & Accessories:

Model	Description
3001KG-AM	Assembled scale in kilograms only
3001KL-AM-BT	3001KL-AM with built-in Pelstar® wireless technology
3001KG-AM-BT	3001KG-AM with built-in Pelstar® wireless technology
3001KL-AMHR	3001KL-AM with mechanical height rod
3001KG-AMHR	3001KG-AM with mechanical height rod
3001KL-AMX	Assembled scale with extended handrails
3001KG-AMX	3001KL-AMX in kilograms only
245EHR-3001	Digital height rod for 3001KL-AM & 3001KG-AM
3001HR	Mechanical height rod for 3001KL-AM & 3001KG-AM
SS-3001KL-AM	2 year extended warranty
SS-3001KL-AMX	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTPRINTERKIT-01	Wireless thermal printer kit

★ TOP SELLER

1100KL

Optional Digital and Mechanical Height Rods Available



DIGITAL PLATFORM SCALE WITH EXTRA WIDE HANDRAILS

This Certified ADA Compliant platform scale supports up to 1000 pounds, and offers “live” handrails on both sides of the slip resistant low profile platform to provide a stable and comfortable weighing experience for patients of any size. The keypad features antimicrobial protection. Ideal for any type of facility including primary care, rehab or hospital.

Fall Prevention: “Live” handrails are advanced technology that allow the patient to hold on for stability while the scale accurately calculates weight.

Low Profile Platform: Requires only minimal step up, offering the patient a comfortable and stable weighing experience.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Complies with MDE Standards: Slip resistant platform, handrail supports on both sides of platform, and 34” min. – 38” max. handrail height from platform



Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

Platform Size (w x d x h)

22" x 15 3/4" x 2 1/2"
559 mm x 400 mm x 64 mm

Display

Digit Size: 3/8" / 10 mm,
Screen Size: 2 3/4" x 1 1/2" High-Contrast Color TFT-LCD Screen

Options & Accessories:

Model	Description
1100KG	Scale in kilograms only
1100KL-BT	1100KL with built-in Pelstar® wireless technology
1100KG-BT	1100KG with built-in Pelstar® wireless technology
1100KLHR	1100KL with mechanical height rod included
1100KGHR	1100KG with mechanical height rod included
1100KL-EHR	1100KL with digital height rod included
1100KL-BT-EHR	1100KL-EHR with built-in Pelstar® wireless technology
245EHR-1100	Digital height rod for 1100 series scales
PROPLUSROD	Mechanical height rod for 1100 series scales
SS-1100KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTPRINTKIT-01	Wireless thermal printer kit

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Parent / Child Weighing Function, Everlock®

Wheels: 2 **Connectivity:** USB, optional Pelstar® wireless technology

2101KL

DIGITAL PLATFORM SCALE



This Certified ADA Compliant bariatric scale features “live” handrails on both sides of the platform that allow the patient to hold on for stability while taking an accurate weight. The keypad features antimicrobial protection. Suitable for use in MD offices, long term care, hospitals or anywhere there is adequate space for the device.

Capacity: 1,000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

26" x 22" x 2 ¾"

660 mm x 559 mm x 60 mm

Display: Digit Size: ¾" / 10 mm,

Screen Size: 2 ¾" x 1 ½" High-Contrast

Color TFT-LCD Screen

Connectivity: USB, optional Pelstar®

wireless technology

Power Source: 100-240V adapter included,
6 D batteries (not included)

Functions: LB / KG Conversion,
LB / KG Lock Out, BMI, Zero, Tare, Pre-tare,
Hold / Release, Recall, Reweigh, Auto Zero,
Auto Off, Time / Date, Variable Auto Off Time,
Audible / Mute Sound Option, Parent / Child
Weighing Function, Everlock®

Options & Accessories:

Model	Description
2101KG	Scale in kilograms only
2101KL-BT	2101KL with built-in Pelstar® wireless technology
2101KG-BT	2101KG with built-in Pelstar® wireless technology
2101KLHR	2101KL with mechanical height rod included
2101KGHR	2101KG with mechanical height rod included
245EHR-2101	Digital height rod for 2101 series scales
STROD	Mechanical height rod for 2101 series scales
SS-2101KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTTPRINTERKIT-01	Wireless thermal printer kit



1110KL

HIGH CAPACITY REMOTE DISPLAY DIGITAL SCALE WITH LARGE BASE



A versatile high-capacity remote digital scale designed to be customizable to accomodate any patient or facility. The keypad features antimicrobial protection.

Capacity: 1000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

22" x 15 ¾" (d) x 2 ⅝"

559 mm x 597 mm x 121 mm

Display: Digit Size: ¾" / 10 mm,

Screen Size: 2 ¾" x 1 ½" High-Contrast

Color TFT-LCD Screen

Connectivity: Via USB, optional

Pelstar® wireless technology

Power Source: 100-240V adapter included,
6 D batteries (not included)

Functions: LB / KG Conversion,
LB / KG Lock Out, BMI, Zero, Tare, Pre-tare,
Hold / Release, Recall, Reweigh, Auto Zero,
Auto Off, Time / Date, Variable Auto Off Time,
Audible / Mute Sound Option, Parent / Child
Weighing Function, Everlock®

Options and Accessories:

Model	Description
1110KG	Scale in kilograms only
1110KL-BT	1110KL with built-in Pelstar® wireless technology
1110KG-BT	1110KG with built-in Pelstar® wireless technology
1110HB	Handlebars for the 1110 series scales
245EHR-1110	Digital height rod for 1110 series scales (must have an 1110HB)
201HR-1110	Mechanical height rod for 1110 series scales (must have an 1110HB)
SS-1110KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTTPRINTERKIT-01	Wireless thermal printer kit





3105KL-AM

★ TOP SELLER

3105KL-AM



Optional Digital and Mechanical Height Rods Available



HEAVY-DUTY ANTIMICROBIAL
PLATFORM SCALE

The 3105KL-AM has a robust, innovative design with enhanced safety and mobility features, "live" handrails allow the patient to hold onto handrails for stability while the scale accurately calculates weight. Low profile platform requires only minimal step up, offering the patient a comfortable and stable weighing experience.

Safety Features: Low-profile platform, "live" handlebars, antimicrobial protection on the head, keypad, transport handle, handrails, body, and kickplate, and an innovative design that eliminates places for dirt and contaminants.

Portability: 2 smooth-rolling wheels and easy grip handle make it easy to move from room to room.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.



Capacity	Resolution	Platform Size (w x d x h)
1000 lb / 454 kg	0.2 lb / 0.1 kg	19 3/4" x 21 1/4" x 2 1/8" 502 mm x 540 mm x 54 mm

Display

Digit Size: 3/8" / 10 mm,
Screen Size: 2 3/4" x 1 1/2" High-Contrast Color TFT-LCD Screen

Options & Accessories:

Model	Description
3105KG-AM	Scale in kilograms only
3105KL-AM-BT	3105KL-AM with built-in Pelstar® wireless technology
3105KG-AM-BT	3105KG-AM with built-in Pelstar® wireless technology
245EHR-3105	Digital height rod for 3105KL-AM series scales
201HR-3105	Mechanical height rod for 3105KL-AM series scales
SS-3105KL-AM	2 year extended warranty
ADPT30	Adapter for 3105KL-AM series scales
BTPRINTERKIT-01	Wireless thermal printer kit

PowerSource: 6 D batteries included,
100-240V adapter not included

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Parent / Child Weighing Function, Everlock®

Connectivity: USB, optional Pelstar® wireless technology



★ TOP SELLER

402KL

MECHANICAL BEAM SCALE

The classic physician scale is built with a robust steel frame and powder coated column to prevent rust. Includes a convenient height rod so clinicians can measure height and weight all in one place.

Height Rod Included: The easy-to-use height rod provides height measurements from 11 $\frac{7}{8}$ " - 84" / 30 cm - 213 cm. Headpiece can be folded down when not in use.

Exclusive 10 Year Limited Warranty: Health o meter® Professional's industry-leading warranty protects against defects of materials or workmanship.

Capacity	Graduation	Platform Size (w x d x h)
390 lb / 180 kg (Counterweights raise to 490 lb / 210 kg)	$\frac{1}{4}$ lb / 100 g	10 $\frac{1}{2}$ " x 14" x 3 $\frac{1}{4}$ " 267 mm x 356 mm x 83 mm
Height Rod	Height Rod Graduation	
11 $\frac{7}{8}$ " - 84" 30 cm - 213 cm	$\frac{1}{8}$ " x 1 mm	

Options & Accessories:

Model	Description
402KLCW	402KL with counterweights included
402KLWH	402KL with wheels included
402KLWHCW	402KL with wheels and counterweights included
402KG50	402KL with 50 gram graduation (special order)
402LB	Scale in pounds only (400 lb capacity)
402LBWH	Scale in LB only with wheels included
402KG	Scale in kilograms only (special order)
55090	Counterweights for 402KL
55000	Wheel brackets for 402KL, 402KG & 402LB

Functions: Measures weight in pounds and kilograms. Measures height in inches and centimeters.

400KL

MECHANICAL BEAM SCALE



A traditional mechanical scale with an industry-leading 10 year warranty. Features a unique coating inside the column to prevent rust when used in gyms, locker rooms and humid environments.

Capacity: 390 lb / 180 kg
(Counterweights raise capacity to 490 lb / 210 kg)

Graduation: ¼ lb / 100 g

Platform Size (w x d x h):

10 ½" x 14" x 3 ¼"

267 mm x 356 mm x 83 mm

Options & Accessories:

Model	Description
400KLCW	Scale with counterweights included
400KLWH	Scale with wheels included
400KLWHCW	Scale with wheels and counterweights included
402KLROD	Optional metal height rod for 400KL
55090	Counterweights for 400KL
55000	Wheel brackets for 400KL

Functions: Measures weight in pounds and kilograms.

450KL

MECHANICAL BEAM SCALE WITH ROTATING POISE BARS



The 450KL is the top of the line mechanical scale and a market leader in doctor offices for decades. Special features include long-lasting engraved numerals on rotating (LB or KG) poise bars and included counterweights that increase capacity to 500 lb / 200 kg.

Capacity: 500 lb / 200 kg (with included counterweights)

Graduation: ¼ lb / 100 g

Height Rod:

11 ⅞" - 84" / 30 cm - 213 cm

Height Rod Graduation: 1/8" x 1 mm

Platform Size (w x d x h):

10 ½" x 14" x 3 ¼"

267 mm x 356 mm x 83 mm

Options & Accessories:

Model	Description
450KLWH	Scale with wheels included
55000	Wheel brackets for 450KL

Functions: Measures weight in pounds and kilograms. Measures height in inches and centimeters.



BCS-G6

BODY COMPOSITION SCALE

Affordable, accurate and attractive. The perfect body composition scale for practices desiring consistent results and accuracy regardless of when it's used. With no pre-test restrictions and proprietary Ohms Modulation Logic technology, G6 scales eliminate the Body Fat Bounce Effect and can be used with patients at any time of day.

15 Second Measurement Time: On-screen display delivers rapid access to weight, body fat %, body water %, lean mass, BMR & BMI. A full illustrated printout can be generated in under a minute.

Large Platform: The spacious platform allows patients to find a comfortable and secure stance.

Few Pretest Restrictions: Analysis can be performed at any time of day and under almost any circumstances, resulting in the greatest flexibility of any body composition analyzer.

High Capacity: The scale features an industry-leading 1000 lb capacity that provides the flexibility and accuracy to weigh any patient.



Capacity	Resolution	Platform Size (w x d x h)
1000 lb / 453.5 kg	0.1 lb / 0.1 kg	14" x 16" x 2 7/8" 356 mm x 406 mm x 48 mm

Display

6-window display (1/2" digits) including weight, % body fat, % body water, % muscle mass, BMI, REE

Optional printer, sold separately, not available through Health o meter® Professional Scales. Visit homscases.com/bcs-g6 for BCS-G6 compatible printer list.

Power Source: 110-220V power adapter (included)

Functions: 15 Second Measurement Time, Body Fat, Body Water, Fat Free Mass, BMI, Basal Metabolic Rate, Visceral Fat, Measures Weight in Pounds and Kilograms

BCS-G6

Users can choose from 4 graphic patient report templates including Adult / Bariatric, Adult Spanish, Pediatric and Fitness versions. These personalized printouts include detailed information about a patient's body composition and educational content regarding the health risks associated with being overweight. Providing patients with a personalized risk assessment is an effective way to engage them in their care and motivate them to take control of their health.



Waiting room brochures (BCS-BROCHURE) are effective marketing tools for practices that charge patients for the G6 analysis. The flyers come in a handy stand and provide a way to inform patients that a Body Composition Analysis is available and can be performed during their visit. The brochure highlights the measurements provided by the analysis and what they can learn from the results.

Model	Description
BCS-G6-ADULT	The whole-body adult only scale is ideal for practices that exclusively see patients 18 and over
BCS-G6-DUO	The whole-body adult/pediatric version is a good choice for providers who see adults and children
BCS-G6-LIMBS	The segmented version provides adults with body composition information about their torso and each of their limbs
BCS-G61-UPPER	The upper-body-only scale provides a partial body composition analysis for adults who do not want to remove shoes and socks
BCS-G6-ADULT-SP	The whole body, adult-only scale is ideal for practices that exclusively see adults and who want a Spanish language report
SS-BCS-G6-ADULT	2 year extended warranty for all versions of the G6 scale
IPO-ADU	Patient Report Template - Adult (5-page report, 25 reports/package)
IPO-PED	Patient Report Template - Pediatric (7-page report, 25 reports/package)
IPO-ATH	Patient Report Template - Athletic (3-page report, 50 reports/package)
IPO-ADU-SP	Spanish Language Patient Report Template - Adult (5-page report, 25 reports/package)
BCS-BROCHURE	Waiting room patient brochure (50/package)
BCS-POWER	Adapter included, item # for replacement or spare



Pediatric

★ TOP SELLER

553KL

DIGITAL PEDIATRIC TRAY SCALE WITH EXTRA-WIDE TRAY

The 553KL is Health o meter® Professional's most popular infant scale and perfect for newborns, infants and toddlers. The scale operates on batteries or AC power so you can use it anywhere. Practices benefit from one-button operation that makes it easy to use while ensuring that the baby remains safely on the scale.

Measuring Tape Included: Easy-to-read measuring tape allows user to measure both weight and height on the same device.

Motion Detection: Scale weighs accurately despite movement of the child on the scale, and locks onto the measurement.



553KL
with 553HR

553KL
with 553EHR



2210CART



Capacity	Resolution	Tray Size (w x d x h)
44 lb / 20 kg	0.5 oz / 10 g	24 7/8" x 13" x 1 3/8" 632 mm x 330 mm x 35 mm

Measuring Tape	Display
0" - 24" 0 cm - 58 cm	1" / 25 mm LCD

Options & Accessories:

Model	Description
553KLCT	Scale with rolling cart included
553KG	Scale in kilograms only
553KL-HR	553KL with attachable mechanical height rod
553KG-HR	553KG with attachable mechanical height rod
553KL-EHR	553KL with attachable digital height rod
553KG-EHR	553KG with attachable digital height rod
SS-553KL	2 year extended warranty
553HR	Attachable mechanical height rod for 553 series scales
553EHR	Attachable digital height rod for 553 series scales
2210CART	Rolling cart for 553KL & 553KG
553CASE	Carrying case for 553KL & 553KG
ADPT31	Adapter for 553KL & 553KG
BTPRINTERKIT-01	Wireless thermal printer kit

Power Source: 6 AA batteries included, 100-240V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lockout, Zero, Tare, Pre-Tare, Hold / Release, Reweigh, Recall, Auto Zero, Auto Off, Motion Detection, Variable Auto-Off Time, Enabling Unit View with Unit of Measure Locked, Auto-Weight Hold for 30 to 60 Seconds, Everlock®

Connectivity: optional Pelstar® wireless technology

522KL

DIGITAL PEDIATRIC TRAY SCALE



Capacity: 50 lb / 23 kg

Resolution: 0 - 20 lb / 0.2 oz; 20 - 50 lb / 0.5 oz; 0 - 9 kg / 5 g; 9 - 23 kg / 10 g

Tray Size (w x d x h):

24 1/8" x 14 1/2" x 2 5/8"

613 mm x 368 mm x 67 mm

Measuring Tape: 0" - 23" / 0 cm - 58 cm

Display: 1" / 25 mm LCD

Connectivity: optional Pelstar®

wireless technology

Power Source: 6 AA batteries included, 100-240V adapter (not included)

Options & Accessories:

Model	Description
522KG	Scale in kilograms only
522KL-HR	522KL with attachable mechanical height rod
522KG-HR	522KG with attachable mechanical height rod
522KL-EHR	522KL with attachable digital height rod
522KG-EHR	522KG with attachable digital height rod
SS-522KL	2 year extended warranty
522HR	Attachable mechanical height rod for 522 series scales
522EHR	Attachable digital height rod for 522 series scales
ADPT31	Adapter for 522 series scales
BTPRINTERKIT-01	Wireless thermal printer kit



Functions: LB / KG Conversion, LB / KG Lock Out, Hold / Release, Zero, Tare, Pre-Tare, Reweigh, Recall, Auto Zero, Auto Off, Audible Mute / Sound Option, Variable Auto-Off Time, Enabling Unit View with Unit of Measure Locked, Auto Weight Hold for 30-60 Seconds, Motion Detection, Everlock®

524KL

DIGITAL PEDIATRIC SEAT SCALE



Capacity: 50 lb / 23 kg

Resolution: 0 - 20 lb / 0.2 oz; 20 - 50 lb / 0.5 oz; 0 - 9 kg / 5 g; 9 - 23 kg / 10 g

Seat Size (w x d x h):

24 1/4" x 14 1/4" x 18 1/2"

616 mm x 362 mm x 470 mm

Display: 1" / 25 mm LCD

Connectivity: optional Pelstar® wireless technology

Options & Accessories:

Model	Description
524KG	Scale in kilograms only
SS-524KL	2 year extended warranty
ADPT31	Adapter for 524KL & 524KG



Power Source: 6 AA batteries included, 100-240V adapter (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, Hold / Release, Zero, Tare, Pre-Tare, Reweigh, Recall, Auto Zero, Auto Off, Audible Mute / Sound Option, Variable Auto-Off Time, Enabling Unit View with Unit of Measure Locked, Auto Weight Hold for 30-60 Seconds, Motion Detection, Everlock®

549KL

DIGITAL PEDIATRIC TRAY SCALE



Lightweight and battery-powered with basic functionality, this scale is a good choice for practices with occasional infant patients or for visiting nurses.

Capacity: 44 lb / 20 kg

Resolution: 0.5 oz / 10 g

Tray Size (w x d x h):

19 1/4" x 13" x 2 5/8"

489 mm x 330 mm x 67 mm

Display: 5/8" / 15 mm LCD

Power Source: 9V battery included

Options & Accessories:

Functions: LB / KG Conversion,

LB / KG Lock Out, Tare, Auto Zero,

Auto Off, Auto Hold, Variable Auto Off

Time, Enabling Unit View with Unit of

Measure Locked, Auto-Weight Hold for

30 to 60 Seconds, Motion Detection,

Everlock®

Model	Description
549KG	Scale in kilograms only
2210CART	Rolling Cart for 549KL & 549KG
64769	Carrying case for 549KL
SS-549KL	2 year extended warranty

550KL

INFANT TO TODDLER, 2 IN 1 DIGITAL SCALE



For practices looking to weigh both infants and toddlers, this dual-purpose scale can be used with the infant weighing tray or toddler step-on platform. An easy-to-use scale that is lightweight and portable with simple functionality and a digital display.

Capacity: 44 lb / 20 kg

Resolution: 0.5 oz / 10 g

Tray Size (w x d x h):

23 1/2" x 8" x 2"

597 mm x 203 mm x 51 mm

Options & Accessories:

Platform Size (w x d x h):

12 1/2" x 8 5/8" x 1 1/4"

318 mm x 219 mm x 32 mm

Display: 1/2" / 13 mm digit height LCD

Power Source: 9V battery included

Functions: Auto Off, Hold, Tare, Unit

Model	Description
SS-550KL	2 year extended warranty
553CASE	Carrying case for 550KL
2210CART	Rolling cart for 550KL
eBABYHR	Digital height rod for 550KL (purchase with 550EHRBRACKET to mount on scale)
550EHRBRACKET	550KL mounting bracket for eBABYHR (purchase with eBABYHR to mount on scale)

8320KL

DIGITAL PORTABLE PEDIATRIC TRAY SCALE



8320KL
Closed

Excellent choice for visiting nurses who need a lightweight pediatric scale. The folding scale is easy to use and compact to carry in its optional bag.

Capacity: 33 lb / 15 kg

Resolution: 0.5 oz / 10 g

Tray Size (w x d x h):

23 1/2" x 11" x 7/8"

597 mm x 279 mm x 22 mm

Display: 3/4" / 19 mm

Options & Accessories:

Power Source: 2 AAA batteries included

Functions: LB / KG Conversion, Zero,

Tare, Auto Zero, Auto Off

Model	Description
SS-8320KL	2 year extended warranty
64771	Carrying case for 8320KL

2210KL-AM

ANTIMICROBIAL HIGH RESOLUTION
DIGITAL NEONATAL / PEDIATRIC TRAY
SCALE



2210KL-AM with
2210EHR



2210CART-AC

2210KL4-AM

ANTIMICROBIAL HIGH RESOLUTION
DIGITAL NEONATAL/PEDIATRIC SCALE
WITH 4-SIDED TRAY



2210KL4-AM

The 2210KL-AM neonatal scale was designed for use in the NICU, delivery room, nursery and pediatric practices where a high degree of accuracy is critical to ensure that babies are gaining weight. A convenient mobile acute care cart designed exclusively for the 2210KL-AM series of scales is available to easily move the scale from location to location.

Capacity: 45 lb / 20 kg / 20,000 gm

Resolution: 1 g (to 10,000 g);
5 g (10,001 to 20,000 g)

Tray Size (w x d x h):

24" x 15 1/4" x 2 3/4"

610 mm x 388 mm x 71 mm

Display: 5 1/4" x 1 1/8" (133 mm x 35 mm),
4 levels of brightness

Connectivity: via USB, optional Pelstar®
wireless technology

Measuring Tape: 0" - 24" / 0 cm - 60 cm

Power Source: 6 D batteries included,
100-240V adapter included

Functions: LB / KG / GM Conversion,
LB / KG / GM Lock Out, Zero, Tare,
Pre-tare, Reweigh, Recall, Auto Zero,
Auto Off, Motion Detection, Everlock®,
Hold / Release, Variable Auto Off Time,
Audible / Mute Sound Option, Disable
/ Enable Backlight for Display, Auto
Weight Hold for 30 to 60 seconds, Breast
Milk Intake

Options & Accessories:

Model	Description
2210KG-AM	Scale in kilograms only
2210KL-AM-BT	2210KL-AM with built-in Pelstar® wireless technology
2210KG-AM-BT	2210KG-AM with built-in Pelstar® wireless technology
2210KL-AM-C	Scale, Canada only
2210KG-AM-C	Scale in kilograms only, Canada only
2210KL-AM-BT-C	2210KL-AM with built-in Pelstar® wireless technology, Canada only
2210KG-AM-BT-C	2210KG-AM with built-in Pelstar® wireless technology, Canada only
2210EHR	Attachable digital height rod for 2210KL-AM (2-sided tray scale)
SS-2210KL	2 year extended warranty
2210CART-AC	Acute care cart for 2210KL-AM series of scales
ADPT30	Adapter included, item # for replacement or spare
BTprinterKIT-01	Wireless thermal printer kit

The 2210KL4-AM neonatal scale comes with a 4-sided weighing tray, which provides an additional level of safety and protection for the baby during the weighing process.

Capacity: 45 lb / 20 kg / 20,000 gm

Resolution: 1 g (to 10,000 g);
5 g (10,001 to 20,000 g)

Tray Size (w x d x h):

26 1/2" x 15 1/4" x 2 3/4"

674 mm x 388 mm x 71 mm

Display: 5 1/4" x 1 1/8" (133 mm x 35 mm),
4 levels of brightness

Connectivity: via USB, optional Pelstar®
wireless technology

Measuring Tape: 0" - 24" / 0 cm - 60 cm

Options & Accessories:

Power Source: 6 D batteries included,
100-240V adapter included

Functions: LB / KG / GM Conversion,
LB / KG / GM Lock Out, Zero, Tare,
Pre-tare, Reweigh, Recall, Auto Zero,
Auto Off, Motion Detection, Everlock®,
Hold / Release, Variable Auto Off Time,
Audible / Mute Sound Option, Disable /
Enable Backlight for Display, Enabling
Unit View with Unit of Measure Locked,
Auto Weight Hold for 30 to 60 Seconds,
Breast Milk Intake

Model	Description
2210KG4-AM	Scale in kilograms only
2210KL4-AM-BT	2210KL4-AM with built-in Pelstar® wireless technology
2210KG4-AM-BT	2210KG4-AM with built-in Pelstar® wireless technology
2210KL4-AM-C	Scale, Canada only
2210KG4-AM-C	Scale in kilograms only, Canada only
2210KL4-AM-BT-C	2210KL4-AM with built-in Pelstar® wireless technology, Canada only
2210KG4-AM-BT-C	2210KG4-AM with built-in Pelstar® wireless technology, Canada only
SS-2210KL	2 year extended warranty
2210CART-AC	Acute care cart for 2210KL-AM series of scales
ADPT30	Adapter included, item # for replacement or spare
BTprinterKIT-01	Wireless thermal printer kit

565KL

HIGH RESOLUTION DIGITAL
BABY SCALE



565KL
with 565EHR

The 565KL was designed for use in critical settings such as the NICU, nursery, delivery room, and pediatric practices. Its exceptional accuracy is especially crucial when caring for the tiniest patients.

- Capacity:** 44 lb / 20 kg
- Resolution:** 0 – 44 lb x 0.1 oz / 0 – 6,000 g x 1g / 6,001 g – 20,000 g x 2 g
- Tray Size (w x d x h):**
24 7/8" x 9 7/8" x 2 3/4"
632 mm x 250 mm x 70 mm
- Display:** Digit Size: 3/4" / 19 mm,
Screen Size: 3" x 1" LCD display
- Connectivity:** optional Pelstar® wireless technology
- Options & Accessories:**
- Measuring Tape:** 0" - 24 5/8" / 0 cm - 62.5 cm
- Power Source:**
100-240V Adapter included, 6 AA batteries not included
- Functions:** LB / KG Conversion,
LB / KG Lock Out, Zero, Tare,
Pre-tare, Reweigh, Recall, Auto Zero, Auto Off, Motion Detection, Everlock®, Hold / Release, Variable Auto Off Time, Audible / Mute Sound Option, Auto Weight Hold for up to 5 Minutes

Model	Description
565KG	Scale in kilograms only
565KL-C	Scale, Canada only
565KG-C	Scale in kilograms only, Canada only
565EHR	Attachable digital height rod for 565 series scales
2210CART	Rolling cart for 565KL & 565KG
SS-565KL	2 year extended warranty
69-00001	Adapter included, item # for replacement or spare
BTTPRINTERKIT-01	Wireless thermal printer kit

1522KL

MECHANICAL PEDIATRIC BEAM SCALE



This high capacity pediatric tray scale is perfect for small infants or young children who can sit independently. The mechanical poise bar does not require electricity, allowing the scale to be used anywhere.

- Capacity:** 135 lb / 62 kg
- Graduation:** 1 oz / 20 g
- Tray Size (w x d x h):**
24" x 16 1/2" x 2 7/8"
610 mm x 420 mm x 70 mm
- Measuring Tape:** 0" - 23 3/4" / 0 cm - 60 cm
- Functions:** Weighs in pounds and kilograms, measures length in inches and centimeters

1524KL

MECHANICAL PEDIATRIC SEAT SCALE



The built-in chair on the 1524KL is convenient for use with infants who are unable to sit; it can also help small children keep still for weight measurements. The mechanical poise bar does not require electricity, allowing the scale to be used anywhere.

- Capacity:** 130 lb / 65 kg
- Graduation:** 1 oz / 20 g
- Seat Size (w x d x h):**
24 1/4" x 14 1/4" x 18 1/2"
622 mm x 362 mm x 470 mm
- Functions:** Measures weight in pounds and kilograms

386S-01 / 386KGS-01

MECHANICAL TRAY SCALE
POUNDS ONLY / KILOGRAMS ONLY



These economical mechanical tray scales are simple to use and include a built-in measuring tape. The sturdy body and convenient handle make it a great choice for visiting nurses.

Capacity:

386S-01: 50 lb

386KGS-01: 25 kg

Graduation:

386S-01: ¼ lb

386KGS-01: 100 g

Tray Size (w x d x h): 19 ¾" x 12 ¾" x 3 ¾"

492 mm x 314 mm x 95 mm

Dial Size: 3" / 76 mm

Measuring Tape: 0" - 18" / 0 cm - 46 cm

Functions:

386S-01: Measures weight in pounds.

386KGS-01: Measures weight in kilograms.

BABYSLINGSCALE

HANGING PEDIATRIC SCALE



This scale is an economical, portable option for use with young children that are able to hold their heads up. The large dial makes it easy to read.

Capacity: 25 kg

Graduation: 100 g

Dial Size: 5 ½" / 140 mm

Functions: Measures weight in kilograms

Options and Accessories:

Model	Description
BABYSLING	Additional sling / trouser

3400KL / 3401KL

DIGITAL WET DIAPER / LAP SPONGE /
ORGAN SCALE



The solid, industrial design of the 3400KL includes a stable, stainless steel base and removable, stainless steel weighing pan (3401KL comes with a flat tray) that can be autoclaved for sterilization. Powered by a long-life rechargeable battery and backlit for any lighting level, this portable scale can be used in any environment, from OR to ED to NICU and beyond.

Capacity: 176 oz / 5,000 g

Resolution: 0.05 oz / 1 g

Base Size (w x d x h):

11 1/2" x 9 7/8" x 5 1/4"

292 mm x 250 mm x 133 mm

Pan Size (3400KL/ 3400KG) (w x d x h):

10" x 9 7/8" x 2 1/2"

254 mm x 250 mm x 62 mm

Tray Size (3401KL) (w x d x h):

7 3/4" x 9 3/4" x 5/8"

197 mm x 248 mm x 16 mm

Display: 5/8" / 16 mm LCD

Power Source: 6V rechargeable

battery included, 100 - 240V AC

battery charger included

Functions: Zero, Tare, Auto Zero,

Auto Off, g / oz Conversion,

Low Battery Indicator, Backlit Display

Options & Accessories:

Model	Description
3400KG	Scale with pan, kilograms only
3401KL	Scale with tray
3400KG-C	Scale with pan, kilograms only, Canada only
SS-3400KL	2 year extended warranty
SS-3400KG	2 year extended warranty
SS-3400KG-C	2 year extended warranty
SS-3401KL	2 year extended warranty
3400PAN	3 pack of stainless steel pans
3400TRAY	3 pack of stainless steel trays
3400ADPT	Battery charger for 3400KL, 3400KG, & 3401KL
20-00018	Battery charger for 3400KG-C



YG500R

MECHANICAL DIAPER /
SMALL PLATFORM SCALE



Perfect for facilities that occasionally need to weigh small items, this cost-effective scale is a great choice.

Capacity: 500 g

Graduation: 2 g

Platform Size (w x d x h):

6 1/2" x 4 3/4"

152 mm x 121 mm

Dial Size: 6" / 152 mm

Functions: Measures weight in grams



**Wheelchair /
Chair / Multi-Use**

★ TOP SELLER

2600KL / 2610KL



Optional 2600ADAKIT
ADA Compliant Upgrade
Kit can be installed on
your existing 2600KL
and 2610KL series of
scales to bring your
scale into compliance
with the new MDE
standards

Power Source: 100-240V adapter included,
6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock
Out, BMI, Zero, Tare, Pre-tare, Hold / Release,
Recall, Reweigh, Auto Zero, Auto Off, Time /
Date, Variable Auto Off Time, Everlock®,
Audible / Mute Sound Option

DIGITAL WHEELCHAIR RAMP SCALE /
DIGITAL WHEELCHAIR DUAL RAMP SCALE

Our most popular wheelchair scale is also the most versatile, and is a good choice for any facility. The platform folds for convenient storage and the casters (the only scale in the market that has them!) make transport a breeze. For easy on/off access the 2610KL offers dual ramps.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Casters and Folding Ramps for Easy Portability and Storage:

4 wheels, including 2 casters, afford ease of movement and the folding ramps make it simple to store.

Wide Platform for Wheelchair or Bariatric Use: Wide platform accommodates wheelchairs of all sizes and provides a comfortable weighing experience for bariatric patients.

Calculates BMI: Scale functions include the ability to calculate patient's BMI.



Capacity	Resolution	Platform Size (w x d x h)
1000 lb / 454 kg	0.2 lb / 0.1 kg	36" x 32 1/4" x 2 1/2" 914 mm x 819 mm x 64 mm

Display

Digit Size: 3/8" / 10 mm, Screen Size: 2 3/4" x 1 1/2" High-Contrast Color TFT-LCD Screen


Wheels: 4,
including 2 casters

Connectivity:
USB, optional Pelstar® wireless technology

Options & Accessories:

Model	Description
2600KG	Scale in kilograms only
2610KG	Scale in kilograms only
2600KL-BT	2600KL with built-in Pelstar® wireless technology
2610KL-BT	2610KL with built-in Pelstar® wireless technology
2600KG-BT	2600KG with built-in Pelstar® wireless technology
2610KG-BT	2610KG with built-in Pelstar® wireless technology
2620KL	ADA compliant scale
2620KG	ADA compliant scale in kilograms only
SS-2600KL	2 year extended warranty
SS-2610KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
2600ADAKIT	ADA Compliance Kit for gray 2600KL & 2610KL series scales
2600ADAKIT-WHITE	ADA Compliance Kit for white 2600KL & 2610KL series scales
B2600RAMP	Secondary ramp in white for 2600KL & 2600KG series scales
BTprinterKIT-01	Wireless thermal printer kit

2620KL



[Learn More](#)



Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option

DIGITAL WHEELCHAIR SCALE WITH DUAL RAMPS AND X-LARGE PLATFORM, **CERTIFIED FOR ADA COMPLIANCE**

Designed to meet the MDE Standards for Accessible Medical Diagnostic Equipment and is 3rd party independently certified in compliance with accessibility codes and standards. The scale includes critical features that meet ADA safety requirements and enhanced safety features to ensure a comfortable weighing experience for patients who use wheelchairs.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Complies with MDE Standards: Dual ramps adhere to a 1:8 slope ratio. The ramps and extra-large 32" x 40" platform feature a 2" edge protection.

Safety Features: Ramps and extra-large platform feature edge protection to prevent roll-offs, and anti-slip caution tape on the ramps enhances safety providing both the clinician and patient confidence during the weighing experience.



Capacity	Resolution	Platform Size
1000 lb / 454 kg	0.2 lb / 0.1 kg	40" x 32 1/4" x 1 7/8" 1016 mm x 819 mm x 48 mm

Display

Digit Size: 3/8" / 10 mm, Screen Size: 2 3/4" x 1 1/2" High-Contrast Color TFT-LCD Screen

Wheels: 4,
including 2 casters

Connectivity:
USB, optional Pelstar® wireless technology

Options & Accessories:

Model	Description
2620KG	Scale in kilograms only
SS-2620KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTPRINTERKIT-01	Wireless thermal printer kit

2500KL

DIGITAL WHEELCHAIR RAMP SCALE



This single-ramp wheelchair scale is perfect for use with standard size wheelchairs. To weigh patients in a seated position, a chair accessory is also available, item # CHAIRACC.

Capacity: 1,000 lb / 454 kg
Resolution: 0.2 lb / 0.1 kg
Platform Size:
29 3/8" (w) x 28 1/4" (d) x 2 3/8" (h)
746 mm (w) x 718 mm (d) x 60 mm (h)
Display: Digit Size: 3/8" / 10 mm
Screen Size: 2 3/4" x 1 1/2" High-Contrast
Color TFT-LCD Screen
Connectivity: USB, optional Pelstar® wireless technology

Wheels: 2
Power Source: 100-240V adapter included, 6 D batteries (not included)
Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Everlock®

Options and Accessories:

Model	Description
2500KG	Scale in kilograms only
2500KL-BT	2500KL with built-in Pelstar® wireless technology
2500KG-BT	2500KG with built-in Pelstar® wireless technology
SS-2500KL	2 year extended warranty
CHAIRACC	Optional fold away seat
ADPT30	Adapter included, item # for replacement or spare
B2500RAMP	Secondary ramp for 2500KL and 2500KG
BTTPRINTERKIT-01	Wireless thermal printer kit



2500CKL

DIGITAL WHEELCHAIR RAMP SCALE WITH FOLD AWAY SEAT



Facilities looking for a scale that can be used with wheelchairs and takes advantage of a built-in chair, should consider the 2500CKL, which features both a ramp and a folding patient seat.

Platform Capacity: 1000 lb / 454 kg
Seat Capacity: 419 lb / 190 kg
Resolution: 0.2 lb / 0.1 kg
Platform Size (w x d x h):
29 3/8" x 28 1/4" x 2 3/8"
746 mm x 718 mm x 60 mm
Seat Size (w x d x h): 18" x 15" x 18 1/2"
457 mm x 381 mm x 470 mm
Display: Digit Size: 3/8" / 10 mm,
Screen Size: 2 3/4" x 1 1/2" High-Contrast
Color TFT-LCD Screen

Connectivity: USB, optional Pelstar® wireless technology
Wheels: 2
Power Source: 100-240V adapter included, 6 D batteries (not included)
Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Audible / Mute Sound Option, Everlock®

Options & Accessories:

Model	Description
2500CKG	Scale in kilograms only
SS-2500CKL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
B2500RAMP	Secondary ramp for 2500CKL
BTTPRINTERKIT-01	Wireless thermal printer kit



Wheelchair Scales/ Multi-use scales

2650KL

DIGITAL WHEELCHAIR DUAL RAMP SCALE



Facilities frequently weighing patients in wheelchairs will appreciate the 2650KL, which offers a remote display head and a medium-sized platform and folds for convenient storage. The dual ramps make it easy to wheel patients on and off the platform. Like all wheelchair scales, the 2650KL can be used for walk-on patients, as well.

Capacity: 1000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

31 ½" x 31 ½" x 2"

800 mm x 800 mm x 51 mm

Display: Digit Size: ¾" / 10 mm,

Screen Size: 2 ¾" x 1 ½" High-Contrast

Color TFT-LCD Screen

Connectivity: USB, optional Pelstar® wireless technology

Options & Accessories:

Wheels: 2

Power Source: 100-240V adapter included, 6 D batteries (not included)


Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option

Model	Description
2650KG	Scale in kilograms only
2650KL-BT	2650KL with built-in Pelstar® wireless technology
2650KG-BT	2650KG with built-in Pelstar® wireless technology
SS-2650KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTPRINTERKIT-01	Wireless thermal printer kit



2700KL

DIGITAL WHEELCHAIR DUAL RAMP SCALE WITH EXTRA LARGE PLATFORM



Featuring the largest platform in the wheelchair scale collection, the 2700KL can easily accommodate even the widest wheelchair. It is a good choice in facilities where the scale will remain in one location.

Capacity: 1000 lb / 454 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

43" x 42" x 2 ½"

1092 mm x 1067 mm x 64 mm

Display: Digit Size: ¾" / 10 mm,

Screen Size: 2 ¾" x 1 ½" High-Contrast

Color TFT-LCD Screen

Connectivity: USB, optional Pelstar® wireless technology

Options & Accessories:

Wheels: 2

Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option

Model	Description
2700KG	Scale in kilograms only
2700KL-BT	2700KL with built-in Pelstar® wireless technology
2700KG-BT	2700KG with built-in Pelstar® wireless technology
SS-2700KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTPRINTERKIT-01	Wireless thermal printer kit



2900KL-AM

ANTIMICROBIAL WALL-MOUNTED FOLD-UP WHEELCHAIR SCALE

The 2900KL-AM is the perfect scale for facilities that need to save space. The low profile scale extends off the wall 3.5" when folded up and has a large easy-to-lift platform that accommodates most wheelchairs.

EMR Capable: Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions.

Easy to Lift Large Platform: large easy-to-lift platform accommodates most wheelchairs and provides a comfortable weighing experience.

Safety Features: Anti-slip caution tape, antimicrobial protection, raised platform edge, recessed handle and no pinch points provide both the clinician and patient confidence during the weighing experience.



Power Source: 100-240V adapter included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option

Connectivity: USB, optional Pelstar® wireless technology

Capacity	Resolution	Platform Size (w x d x h)
1000 lb / 454 kg	0.2 lb / 0.1 kg	29 ½" x 32 ½" x 2 ¾" 749 mm x 819 mm x 70 mm
Display		

Digit Size: ¾" / 10 mm, Screen Size: 2 ¾" x 1 ½" High-Contrast Color TFT-LCD Screen

Options & Accessories:

Model	Description
2900KG-AM	Scale in kilograms only
2900KL-AM-BT	2900KL- AM with built-in Pelstar® wireless technology
2900KG-AM-BT	2900KG-AM with built-in Pelstar® wireless technology
SS-2900KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTPRINTKIT-01	Wireless thermal printer kit

2400KL

DIGITAL PORTABLE WHEELCHAIR SCALE



The 2400KL is the perfect scale for locations where occasional wheelchair-bound patients need to be weighed. The scale folds easily and requires minimal room for storage. The scale includes a handle for portability.

Capacity: 800 lb / 363 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size: Expands and adjusts from 25 3/4" - 50 1/8" / 654 mm - 1273 mm

Display: Digit Size: 3/8" / 10 mm, Screen Size: 2 3/4" x 1 1/2" High-Contrast Color TFT-LCD Screen

Connectivity: USB, optional Pelstar® wireless technology

Wheels: 2 wheels, 2 casters

Power Source: 100-240V adapter included, 6 D batteries (not included)


Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Recall, Reweigh, Auto Zero, Auto Off, Time / Date, Variable Auto Off Time, Everlock®, Audible / Mute Sound Option

Model	Description
2400KG	Scale in kilograms only
2400KL-BT	2400KL with built-in Pelstar® wireless technology
2400KG-BT	2400KG with built-in Pelstar® wireless technology
SS-2400KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTprinterKIT-01	Wireless thermal printer kit



2842KL

LARGE PLATFORM DIGITAL SCALE



With its large platform and stainless steel base, the 2842KL can serve in multiple capacities. It can be used as a basic stand-on scale with an extra-large platform, and it can also be used in shipping/receiving or in laundry areas.

Capacity: 600 lb / 270 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

22 1/4" x 42" x 2 1/2"

565 mm x 1067 mm x 64 mm

Display: 7/8" / 22 mm LCD

Options & Accessories:

Power Source: 100-240V adapter

included, 6 C-cell batteries (not included)

Functions: LB / KG Conversion, Zero, Tare, Motion Detection, Reweigh

Model	Description
SS-2842KL	2 year extended warranty
ADPT31	Adapter included, item # for replacement or spare



2595KL

DIGITAL CHAIR SCALE



This best-selling chair scale has a comfortable seat so patients can sit safely while being weighed. The moveable armrests and footrests allow for easy access on and off the scale. A great scale for use in Long-Term Care centers, rehab facilities, and senior living homes.

Capacity: 600 lb / 272 kg

Resolution: 0.2 lb / 0.1 kg

Seat Size (w x d x h):

18 ¼" (w) x 15 ⅝" (d) x 16 ½"

464 mm x 397 mm x 419 mm

Display: 1 ⅞" / 29 mm LCD

Connectivity: optional Pelstar® wireless technology

Options and Accessories:

Wheels: 2

Power Source: 100-240V adapter included, 6 C-cell batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, Zero, Pre-tare, Hold / Release, Recall, Reweigh, Tare, Auto Off, Variable Auto Off Time, Enabling Unit View with Unit of Measure Locked, Everlock®

Model	Description
2595KG	Scale in kilograms only
SS-2595KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare



594KL

DIGITAL CHAIR SCALE



The 594KL has moveable arm and foot rests to allow for easy patient access on and off this economical chair scale. Two wheels make portability easy for bringing the scale to the patient's location.

Capacity: 600 lb / 270 kg

Resolution: 0.2 lb / 0.1 kg

Seat Size (w x d):

18" x 18 ½"

457 mm x 470 mm

Display: 1" / 25 mm LCD

Wheels: 4 casters

Power Source: 100-240V adapter included, 6 AA batteries (not included)

Functions: LB / KG Conversion, LB / KG Lock Out, BSA, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off

Options and Accessories:

Model	Description
594KL-C	Scale, Canada only
SS-594KL	2 year extended warranty
69-00001	Adapter included, item # for replacement or spare

445KL

MECHANICAL CHAIR SCALE



Nursing homes, rehab hospitals and facilities that treat elderly or unstable patients may find this basic chair scale an appropriate choice for weighing their patients.

Capacity: 440 lb / 200 kg

Graduation: ¼ lb / 100 g

Seat Size (w x d x h):

18 ½" x 17 ¼" x 15 ¾"

470 mm x 438 mm x 400 mm

Wheels: 2

Functions: Measures weight in pounds and kilograms





Specialty

★ TOP SELLER

PTS-1000KG

Patient Transfer Scale, Kilograms Only



Power Source: 3.7V rechargeable battery

Functions: Zero, Hold / Release, Auto Off

Always use a Slide Sheet when using the Patient Transfer Scale



Check out the Single Patient Use and Reusable SlideMate™



Capacity	Resolution	Weighing Mat Size (w x d):
250 kg	0.5 kg	71" x 21" 1805 mm x 533 mm
Measuring Tape		Display

100 cm -180 cm height guide on mat

3" x 1" Backlight LCD Screen

Options & Accessories:

Model	Description
PTS-1000KL	Scale, lbs and kilograms
PTS-1000KL-C	Scale, lbs and kilograms, Canada only
PTS-1000KG-C	KG only scale, Canada only
SS-PTS-1000KL	2 year extended warranty
69-00012	Battery charger included, item # for replacement or spare

3400KL / 3401KL

DIGITAL WET DIAPER / LAP SPONGE / ORGAN SCALE



The solid, industrial design of the 3400KL includes a stable, stainless steel base and removable, stainless-steel weighing pan (3401KL comes with a flat tray) that can be autoclaved for sterilization. Powered by a long-life rechargeable battery and featuring a backlit display for all lighting levels, this portable scale can be used in environments, from OR to ED to NICU and beyond.

Capacity: 176 oz / 5,000 g

Resolution: 0.05 oz / 1 g

Base Size (w x d x h): 11 1/2" x 9 7/8" x 5 1/4"
292 mm x 250 mm x 133 mm

Display: 5/8" / 16 mm LCD

Pan Size (w x d x h): 10" x 9 7/8" x 2 1/2"
254 mm x 250 mm x 64 mm

Options and Accessories:

Tray Size: 7 3/4" (w) x 9 3/4" (d) x 5/8" (h)

197 mm (w) x 248 mm (d) x 16 mm (h)

Power Source: 6V rechargeable battery included, 100 - 240V AC battery charger

Functions: Zero, Tare, Auto Zero, Auto Off, g / oz Conversion, Low Battery Indicator, Backlit Display

Model	Description
3400KG	Scale with pan, kilograms only
3401KL	Scale with tray
3400KG-C	Scale with pan, kilograms only, Canada only
SS-3400KL	2 year extended warranty
SS-3400KG	2 year extended warranty
SS-3400KG-C	2 year extended warranty
SS-3401KL	2 year extended warranty
3400PAN	3 pack of stainless steel pans
3400TRAY	3 pack of stainless steel trays
3400ADPT	Battery charger for 3400KL, 3400KG, & 3401KL
20-00018	Battery charger for 3400KG-C



2000KL

DIGITAL IN-BED / STRETCHER SCALE



This scale can weigh immobile patients up to 400 lb / 181 kg. It has a durable steel structure, fabric sling, and wheels for portability.

Capacity: 400 lb / 181 kg

Resolution: 0.2 lb / 0.1 kg

Stretcher Size (w x d):

70 1/2" x 30"

1791 mm x 762 mm

Display: Digit Size: 3/8" / 10 mm,

Screen Size: 2 3/4" x 1 1/2" High-

Contrast Color TFT-LCD Screen

Connectivity: USB, optional Pelstar®

wireless technology

Options and Accessories:

Wheels: 4 (1 with safety lock)

Power Source: 100-240V adapter

included, 6 D batteries (not included)

Functions: LB / KG Conversion, LB / KG

Lock Out, BMI, Zero, Tare, Pre-tare, Hold /

Release, Recall, Reweigh, Auto Zero, Auto

Off, Time / Date, Variable Auto Off Time,

Everlock®, Audible / Mute Sound Option

Model	Description
2000KG	Scale in kilograms only
2000KL-BT	2000KL with built-in Pelstar® wireless technology
2000KG-BT	2000KG with built-in Pelstar® wireless technology
SS-2000KL	2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
BTTPRINTERKIT-01	Wireless thermal printer kit



YG500R

MECHANICAL DIAPER / SMALL PLATFORM SCALE



Perfect for facilities that occasionally need to weigh small items, this cost-effective scale is a great choice.

Capacity: 500 g

Graduation: 2 g

Platform Size (w x d):

6 1/2" x 4 3/4"

152 mm x 121 mm

Dial Size: 6" / 152 mm

Functions: Mechanically measures weight in grams.



Professional Home Care

810KL

DIGITAL GLASS SCALE WITH ANTI-SLIP
TREAD & BACKLIGHT



The 810KL is a high capacity medical grade home care scale with exclusive anti-slip tread and an easy-to-read backlit display. Designed for use in the bathroom / home, with accuracy and features of a professional medical scale.

Capacity: 440lb / 200 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

13" x 13" x 7/8"

311 mm x 311 mm x 44 mm

Display: 1 7/8" LCD

Options & Accessories:

Power Source: (2) 3V lithium batteries included

Functions: LB / KG Switch, Auto Zero, 5 Second Weight Hold, Auto Off

Model	Description
64771	Carrying case for 810KL

800KL

DIGITAL FLOOR SCALE (PACK OF 4)



Our popular digital bath-style scale is great for use in the home with its simple design and its ability to match any décor. A great option for visiting nurses and health assessment professionals.

Capacity: 397 lb / 180 kg

Resolution: 0.1 lb / 0.05 kg

Platform Size (w x d x h):

12 1/4" x 12 1/4" x 1 3/4"

311 mm x 311 mm x 44 mm

Display: 1 1/2" / 38 mm LCD

Options & Accessories:

Power Source: (2) AA batteries included

Functions: LB / KG Switch, Auto Zero, 5 Second Weight Hold, Auto Off

Model	Description
800KLS	Digital floor scale (single)
64771	Carrying case for 800KL

880KL

HEAVY-DUTY DIGITAL FLOOR SCALE
(PACK OF 4)



With its non-slip surface, the 880KL is the best choice for large or unstable patients. The robust design offers the largest platform available and an extra-large digital display. The display's size makes it easier to read than any of our other bath scales and is great for use in a clinic setting or for home health.

Capacity: 500 lb / 220 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

17" x 13" x 1 1/4"

432 mm x 330 mm x 33 mm

Display: 2 1/2" / 64 mm LCD

Options & Accessories:

Power Source: (2) 3V lithium batteries included

Functions: LB / KG / ST Switch, Auto Off, 10 Second Weight Hold, Zero

Model	Description
880KLS	Heavy-duty digital floor scale (single)
64769	Carrying case for 880KL



All Health o meter® Professional scales are manufactured in FDA-registered factories.

844KL

DIGITAL FLOOR SCALE (PACK OF 2)



The 844KL offers a high capacity paired with the highest resolution in the bath scale line. A good option for patients who need greater precision in their weight measurement.

Capacity: 440 lb / 200 kg

Resolution: 0.1 lb / 0.05 kg

Platform Size:

12 5/8" (w) x 12 5/8" (d) x 1 7/8"

321 mm x 321 mm x 48 mm

Display: 1 1/4" / 32 mm LCD

Options & Accessories:

Power Source: (2) AAA

batteries included

Functions: LB / KG / ST Switch,

Auto Zero, Auto Off, 5 Second

Weight Hold

Model	Description
844KLS	Digital floor scale (single)
64771	Carrying case for 844KL

349KLX

DIGITAL FLOOR SCALE



The 349KLX is a versatile scale that can be used in a host of locations. Users can place the remote display head where it is easiest for them to read the weight. The scale's small footprint makes it a great choice for home use.

Capacity: 400 lb / 180 kg

Resolution: 0.2 lb / 0.1 kg

Platform Size (w x d x h):

12 1/2" (w) x 12" (d) x 1 7/8"

318 mm x 305 mm x 48 mm

Display: 1" / 25 mm LCD

Power Source: 6 AA batteries included, 120V adapter (not included)

Options & Accessories:

Functions: LB / KG Conversion,

LB / KG Lock Out, Zero, Tare, Hold/

Release, Reweigh, Auto Zero, Auto Off,

Variable Auto Off Time, Parent / Child

Weighing Function, Enabling Unit View

with Unit of Measure Locked, Enable /

Disable Backlight Display with 4 Levels

of Brightness, Everlock®

Model	Description
349KLXAD	Scale with adapter included
349KGX	Scale in kilograms only
SS-349KLX	2 year extended warranty
ADPT40	Adapter for 349KLX
64771	Carrying case for 349KLX

322KL



★ TOP SELLER

822KL



DIGITAL FLOOR SCALE (PACK OF 2)

This easy to transport scale offers high resolution with a wide platform to make it easier for patients with a wide stance. A versatile scale that is ideal for many uses including traveling healthcare professionals, diet groups, and home use.

Large 1 ½" LCD Display: The large display shows the user an easy-to-read weight measurement.

5-Second Weight Hold: Weight measurement will stay displayed for 5 seconds allowing the user time to stabilize themselves after stepping off the scale.

Low Profile Platform: Low platform height requires only minimal step up, offering patient a comfortable and stable weighing experience.

Capacity

400 lb / 180 kg

Resolution

0.1 lb / 0.05 kg

Platform Size (w x d x h):

12 ¾" x 10 ¾" x 1 ½"
349 mm x 273 mm x 38 mm

Display

1 ½" LCD

Options & Accessories:

Model	Description
822KLS 64771	Digital floor scale (single) Carrying case for 822KL

Power Source: (2) AAA batteries included

Functions: LB / KG Switch, Auto Zero,
5 Second Weight Hold, Auto Off

100LB / 100KG

MECHANICAL FLOOR SCALE -
POUNDS ONLY / KILOGRAMS ONLY
(PACK OF 3)



Simple and lightweight, these mechanical scales are ideal for home use or by visiting nurses.

Capacity:

100LB: 270 lb

100KG: 125 kg

Graduation:

100LB: 1 lb

100KG: .5 kg

Platform Size (w x d x h):

10 ¼" x 10" x 1 ⅞"

260 mm x 254 mm x 48 mm

Functions:

100LB: Measures weight in pounds.

100KG: Measures weight in kilograms.

Dial Size: 2" / 51 mm

Options and Accessories:

Model	Description
100LBS	Mechanical floor scale – pounds only (single)
100KGS	Mechanical floor scale – kilograms only (single)
64771	Carrying case for 100LB/100KG

175LB

MECHANICAL FLOOR SCALE -
POUNDS ONLY (PACK OF 2)



When only pounds measurements are required, the 175LB offers a large, easy to read dial on a low profile platform. A great product for home use.

Capacity: 330 lb

Graduation: 1 lb

Platform Size (w x d x h):

11 ⅜" x 9 ¾" x 2"

289 mm x 248 mm x 51 mm

Functions: Mechanically measures weight in pounds.

Dial Size: 3" / 89 mm

Options and Accessories:

Model	Description
64771	Carrying case for 175LB

142KL

MECHANICAL FLOOR SCALE
(PACK OF 2)



For use in the home or clinic, the 142KL is a good option for patients who require a spacious platform and large dial. The scale is also a good choice for weight loss programs.

Capacity: 330 lb / 150 kg

Graduation: 1 lb / 0.5 kg

Platform Size:

11 ½" x 10 ¾" x 2 ¼"

292 mm x 273 mm x 57 mm

Dial Size: 6 ¼" / 159 mm

Options and Accessories:

Functions: Measures weight in pounds and kilograms.

Model	Description
142KLS	Mechanical floor scale (single)

160KL

MECHANICAL FLOOR SCALE
(PACK OF 2)



This classic "Bigfoot" scale has a deep platform with non-slip surface and an 8" easy-to-read raised dial—good for patients with vision challenges or those who prefer more foot room.

Capacity: 400 lb / 180 kg

Graduation: 1 lb / 1 kg

Platform Size (w x d x h):

11" x 12 ½" x 3"

279 mm x 318 mm x 76 mm

Dial Size: 8" / 203 mm, raised 4" / 102

mm

Functions: Measures weight in pounds and kilograms.

Model	Description
160LB	Mechanical floor scale—pounds only (2 pack)
160KG	Mechanical floor scale—kilograms only (2 pack)
160KLS	Mechanical floor scale (single)
160LBS	Mechanical floor scale—pounds only (single)
160KGS	Mechanical floor scale—kilograms only (single)

All Health o meter® Professional scales are manufactured in FDA-registered factories.



Height Measuring

201HR

WALL-MOUNTED TELESCOPIC
METAL HEIGHT ROD



★ TOP SELLER

Sleek and robust measuring device requires minimal wall space and has a foldable headpiece. Designed for long life and to withstand dozens of uses every day.

Measurement Range: 2" - 90" / 6 cm - 230 cm

Graduation: 1/16" / 1 mm

205HR

HIGH-STRENGTH WALL-MOUNTED
HEIGHT ROD



Manufactured with a durable high-strength plastic polymer, this sturdy stadiometer resists bending and ensures that the headpiece slides smoothly. Featuring a double sided measuring tape, it offers a choice of centimeters only or feet and inches.

Measurement Range: 55.5" / 141 cm

Measuring Tape: 0" - 90" / 0 cm - 230 cm

Graduation: 1/16" / 1 mm



PORTROD

WALL-MOUNTED HEIGHT ROD



With its large headpiece and clear indicator the PORTROD stadiometer makes it easy to read a patient's height. Its wall mount design is easy and quick requiring only simple tools and provided hardware.

Measurement Range: 23 5/8" - 83" / 60 cm - 213 cm

Graduation: 1/16" / 1 mm

WMROD

402KLROD WITH ANCHORS AND SCREWS



Included on the 402KL scale, this metal height rod can be purchased separately to mount on the wall for adult or pediatric use.

Measurement Range: 11 7/8" - 84" / 30 cm - 213 cm

Graduation: 1/8" / 1 mm

221HR

TAA COMPLIANT PORTABLE STADIOMETER WITH CARRYING CASE



Designed to meet the needs of doctor's offices, schools, home health and traveling clinicians, features an easy to assemble/ disassemble structure and a large sturdy anti-slip base. The long headpiece adds comfort for patients of all sizes and the large view finder provides an easy-to-read accurate measurement.

Measurement Range: 7" - 82.6" / 17.7 cm - 210 cm

Graduation: 1/16" / 1 mm

Carrying Case Size (w x d x h): 25" x 5" x 15" / 635 mm x 127 mm x 381 mm



280HR

INFANT MEASURING MAT



The space-saving mat is soft, lightweight, foldable, and easy to store. The fixed headpiece and sliding foot positioner allow for precise measurement of infants or toddlers while lying down.

Measurement Range: 3.9" - 43.3" / 10 cm - 110 cm

Graduation: 1/4" / 5 mm

EBABYHR

DIGITAL PEDIATRIC MEASURING ROD



The eBABYHR makes measuring infant length fast and easy with clear digital measurements, large head and foot paddles, and convenient Hold function.

Measurement Range: 13 3/4" - 31 1/2" / 35 cm - 80 cm

Resolution: 1/16" / 0.1 cm

BABYHR

MECHANICAL PEDIATRIC MEASURING ROD



The BABYHR is perfect for practitioners who need a tool to measure infant length. It can be mounted on the wall along an exam table and users can fold the headpiece and footpiece when not in use.

Measurement Range: 3" - 39" / 7 cm - 100 cm

Graduation: 1/16" / 1 mm

* Information about height rods for specific scales can be found on each scale's product page.



Accessories

64769

CARRYING CASE FOR 549KL, 880KL



Features: Nylon bag protects scale during transport. Convenient for mobile use

553CASE

CARRYING CASE FOR 553KL & 550KL



Features: Nylon bag protects scale during transport. Convenient for mobile use

64771

CARRYING CASE FOR 349KLX, 498KL, 752KL, 100LB/KG, 175LB, 800KL, 822KL, 844KL, 8320KL



Features: Reinforced nylon bag protects scale during transport.

2210CART-AC

ACUTE CARE CART FOR 2210KL-AM / 2210KL4-AM SERIES OF SCALES



Features: 2 easy grip handles and customized plates for secure scale placement of the 2210KL-AM series of scales

Storage Shelf (w x d x h):

21 1/4" x 13" x 3"

540 mm x 330 mm x 76 mm

Capacity: 77 lb / 35 kg

Cart Dimensions (w x d x h):

30" x 22" x 38"

762 mm x 59 mm x 965 mm

Wheels: 4 large anti-static casters with locks

2210CART

ROLLING SCALE / MEDICAL DEVICE CART FOR USE WITH 549KL, 550KL, 553KL, 565KL, 2210KL-AM, 2210KL4-AM



Features: Lockable cabinet and storage shelf.

Cabinet Dimensions (w x d x h):

23" x 16 3/4" x 17"

584 mm x 427 mm x 432 mm

Capacity: 200 lb / 90 kg

Cart Dimensions (w x d x h):

26 1/2" x 18" x 34 1/8"

673 mm x 457 mm x 867 mm

Wheels: 4 casters (2 with locks)



2600ADAKIT

ADA COMPLIANCE KIT FOR GRAY 2600KL AND 2610KL SERIES SCALES, CERTIFIED FOR ADA COMPLIANCE



2600ADAKIT-WHITE
Available for
White Scales

Features: Upgrade kit includes new ramps, feet, edge guards, and platform extension that meet MDE requirements

Edge Guard Dimensions:

40" x 1/8" x 2"

819 mm x 406 mm x 51 mm

Ramp Dimensions:

32 1/4" x 16"

819 mm x 406 mm

Platform Dimensions (with 2600KL):

40" x 32 1/4" x 1 7/8"

1016 mm x 819 mm x 48 mm

ADPT30

100-240V POWER ADAPTER FOR USE WITH:
1100KL, 1110KL, 2000KL, 2101KL, 2210KL-AM,
2400KL, 2500KL, 2500CKL, 2595KL, 2600KL,
2610KL, 2650KL, 2700KL, 2900KL-AM,
3001KL-AM, 3001KL-AMX, 3105KL-AM

ADPT31

100-240V POWER ADAPTER FOR USE WITH:
498KL, 499KL, 500KL, 502KL, 522KL, 524KL,
553KL, 2842KL

ADPT40

120V POWER ADAPTER FOR USE WITH:
349KLX

ADPT50

100-240V POWER ADAPTER FOR USE WITH:
599KL, 600KL, 600KL-BT, ELEVATE-BT, 752KL,
AND ELEVATE-C

69-00001

100-240V POWER ADAPTER FOR USE WITH:
594KL, 565KL

69-00012

100-240V BATTERY CHARGER FOR USE WITH:
PTS-1000KL

55000

WHEEL BRACKET FOR HEALTH O METER®
PROFESSIONAL 400 SERIES BEAM SCALES

Features: Wheel bracket attaches easily to scale for mobility.

55070

COUNTERWEIGHTS FOR 450KL
MECHANICAL BEAM SCALE

Features: Increases capacity of 450KL beam scales.

55090

COUNTERWEIGHTS FOR 400KL AND 402KL

Features: Increases capacity of 400KL and 402KL beam scales.

B2500RAMP

SECONDARY RAMP FOR 2500KL, 2500CKL

B2600RAMP-WHITE

SECONDARY RAMP FOR WHITE 2600KL, 2610KL

BTPRINTERKIT-01

WIRELESS THERMAL PRINTER KIT

Features: Lightweight printer connects to compatible scales via Pelstar wireless technology. Powered by a rechargeable lithium battery, the device prints patient weight, height and BMI (if available) on 58 mm (2 ¼") thermal paper.

BMICHARTKG

BMI WHEEL (METRIC)



Features: Handy tool quickly displays BMI by selecting patient height and weight. (metric units only)

MUACTAPE

MIDDLE UPPER ARM CIRCUMFERENCE TAPE



Features: Measures nutritional status in children ages 10-14 years old.



Automated Calibrator

CALIBOT

Automated Calibrator

The Automated Calibrator is a robust innovative calibration system that allows biomedical engineers and other medical technicians to calibrate medical patient scales quickly and safely without having to transport and stack physical calibration weights on the scale. Eliminating the cumbersome physical calibration weights from the calibration process reduces the time needed for routine scale maintenance and eliminates the safety risks associated with handling heavy dead weights.



CALIBOT-AC



Capacity

Resolution

550 lb / 250 kg

0.01 lb / 0.005 kg

Platform Size CALIBOT-AC (w x d):

32" x 27"
813 mm x 686 mm

Product Weight CALIBOT-AC

97 lb / 44 kg

Platform Size CALIBOT-PC (w x d):

17" x 17"
432 mm x 432 mm

Product Weight CALIBOT-PC

65 lb / 30 kg

Options & Accessories:

Model	Description
CALIBOT-AC	Acute Care Calibrator
CALIBOT-PC	Primary Care Calibrator



CALIBOT-PC

CALIBOT not available in all countries.

Glossary/Symbols

Audible / Mute Sound Option

The display can be set to stay silent or beep while operating the scale.

Auto Off

Digital scale will automatically shut down after a designated period of time, reducing electrical usage and sustaining battery life.

Auto Weight Hold

Scale continues displaying the weight measurement for a set period of time after the patient is off the scale.

Auto Zero

Ability of a digital scale to show a "0" value automatically when powered up.

Breast Milk Intake

Scale calculates the patient's breast milk intake after a feeding.

BMI

A statistical measure that is the ratio between a person's weight and height and interpreted as a BMI score, a.k.a. Quetelet Index. Scales with this feature can calculate BMI score based on digitally measured or manually entered height data.

Connectivity

Ability of the scale to interface with a computer, vital signs monitor or other electronic device to transmit weight, height and/or BMI.

Enabling Unit View with Unit of Measure Locked

Allows the user to briefly display the weight in the alternate weight measuring unit.

Everlock®

An exclusive feature that allows a scale to be permanently locked in kilogram or pound units.

Factory Calibrated

All Health o meter® Professional scales are inspected and factory-calibrated prior to shipment.

Hold / Release

Digital function that allows the user to hold the patient weight until the key is pressed again to release that information.

LB / KG, LB / KG / GM, or G / OZ Conversion

Allows weight to be displayed in different units while patient remains on the scale.

LB / KG Lockout

Allows an administrator to select and use only pounds / inches or kilograms / centimeters, helping to reduce the risk of data being recorded in the incorrect unit of measurement.

LB / KG /ST or LB / KG Switch

Switch located on the scale platform that allows the user to select weight measuring unit.

Live Handrails

Allows patient to hold onto handrails for stability while scale accurately calculates weight.

Low Profile Platform

Low platform height requires a minimal step up, offering patient a comfortable and stable weighing experience.

Motion Detection

Digital feature that enables accurate weight despite movement of the subject on the scale.

Parent / Child Weighing Function

The scale's Tare function provides the convenience of weighing a small child or infant while being held by a parent.

Pre-tare

Entering a known weight of an object such as a wheelchair, to be subtracted from the patient's weight.

Recall

Displays the last weight measurement.

Remote Display

The scale's display head is separate from the weighing platform and can be handheld, mounted on a wall or placed on a table.

Resolution / Graduation

The degree of difference/minimum weight change that a scale maintains when measuring weight, such as fraction of an ounce or pound.

Reweigh

Allows repeated weighing while the patient remains on the scale.

Swivel Display

Display head rotates for the convenience of the caregiver, allowing line-of-sight viewing and information privacy.

Tare

Subtracts the weight of an object, such as a diaper or wheelchair, from the patient's weight measurement.

Variable Auto-off Time

Ability to choose the period of inactivity time before the scale powers off.

Zero

Zeroes the display to set the weight/reference point to zero.



cMETus - Products have been evaluated for safety by MET, an accredited independent lab, in accordance with the appropriate national standards for the United States and Canada.



Ships By Truck Only - Due to their size, products can only be shipped by a common carrier. Please call 1 (800) 815-6615 for specific shipping information to your facility.



TAA Compliant - The manufacturing of these products meets the rules and requirements of the U.S. Trade Agreements Act.


All Health o meter® Professional products are listed with the US Food & Drug Administration.

All Health o meter® Professional products are manufactured without the use of natural rubber latex.


Global Support

We have over 170 skilled professionals based in the United States and Canada who manufacture, sell, service and support our products. Our global headquarters in McCook, IL is home to our highly trained Customer Service and Technical Support teams along with our engineering, logistics and management team. On-site sales and support are provided in partnership with MedTech/MedCare, an industry-recognized leader in medical sales throughout the US and Canada.

Our American engineered products are designed and manufactured to meet the requirements of the International Organization for Standardization (ISO).



We take pride in designing world-class products for healthcare professionals around the globe. Our sales and support network includes select dealers and distributors throughout Latin America, Asia, Africa and the Middle East.





Health o meter®

P r o f e s s i o n a l

Contact Information

Main Number: 1 (708) 377-0600

Customer Service: 1 (800) 815-6615

Technical Support: 1 (800) 638-3722

Fax: 1 (708) 377-0601 or 1 (800) 654-7330

Email: homprocs@homscales.com

Corporate Headquarters:

Pelstar LLC

9500 West 55th St.

McCook, IL 60525-7110 USA

www.homscales.com

Health o meter® Professional products are manufactured, designed and owned by Pelstar LLC. Health o meter® is a registered trademark of Sunbeam, Inc. used under license. ©2025 Pelstar® LLC – PRODUCT CAT-ONLINE 20250728