

Baxter

CONNECTING
HEALTH O METER®
PROFESSIONAL SCALES
& VITAL SIGNS DEVICES



HEALTH O METER®
PROFESSIONAL SCALE



CONNECTIVITY
KIT







VITAL SIGNS
DEVICE



ELECTRONIC
MEDICAL RECORD

The following table outlines the specific **Health o meter®** Professional Scales that can be used with/are compatible with both the **Connex** Vital Signs Monitor 6000 Series (including the **Connex** Integrated Wall System) and the **Connex** Spot Monitor. **Health o meter®** Professional Scales are sold through major medical supply distributors in the U.S. and Canada.

HEALTH O METER® PROFESSIONAL SCALES

Model Number	Scale Description	Connection Type	Scale Family	
522KL	Pediatric Digital Tray Scale	Bluetooth®	Pediatric	
524KL	Pediatric Digital Seat Scale	Bluetooth	Pediatric	
553KL	Pediatric Digital Tray Scale	Bluetooth	Pediatric	
565KL	High-Resolution Digital Baby Scale	Bluetooth	Pediatric	
2210KL-AM	Antimicrobial High-Resolution Digital Neonatal/Pediatric Tray Scale	Bluetooth	Pediatric	
2210KL4-AM	Antimicrobial High-Resolution Digital Neonatal/Pediatric Four-Sided Tray Scale	Bluetooth	Pediatric	
498KL	Remote Display Digital Scale	Bluetooth	Stand-on	
500KL	Eye Level Digital Scale	Bluetooth	Stand-on	
502KL	Digital Physician Scale with Integrated Digital Height Rod	Bluetooth	Stand-on	
599KL	Heavy Duty Waist Level Digital Scale	Bluetooth	Stand-on	
600KL	Heavy Duty Eye Level Digital Scale with Digital Height Rod	Bluetooth	Stand-on	
752KL	Heavy Duty Remote Display Digital Scale	Bluetooth	Stand-on	
1100KL	Digital Platform Scale	Bluetooth	Stand-on/Bariatric	
1110KL	High-Capacity Large Base Remote Display Scale	Bluetooth	Stand-on/Bariatric	
2101KL	Digital Platform Scale	Bluetooth	Stand-on/Bariatric	
3001KL-AM	Digital Platform Antimicrobial Scale	Bluetooth	Stand-on/Bariatric	
3105KL-AM	Antimicrobial Heavy Duty Platform Scale	Bluetooth	Stand-on/Bariatric	
2500KL	Digital Wheelchair Ramp Scale	Bluetooth	Wheelchair/Bariatric	
2600KL	Digital Wheelchair Ramp Scale	Bluetooth	Wheelchair/Bariatric	
2610KL	Digital Wheelchair Dual Ramp Scale	Bluetooth	Wheelchair/Bariatric	
2620KL	Digital Wheelchair Scale With Dual Ramps And Extra-Large Platform, Certified for ADA Compliance	Bluetooth	Wheelchair/Bariatric	
2650KL	Digital Dual Ramp Wheelchair Scale	Bluetooth	Wheelchair/Bariatric	
2700KL	Digital Dual Ramp Wheelchair Scale with Extra Large Platform	Bluetooth	Wheelchair/Bariatric	
2900KL-AM	Antimicrobial Wall-Mounted Fold-Up Wheelchair Scale	Bluetooth	Wheelchair/Bariatric	
2400KL	Digital Portable Wheelchair Scale	Bluetooth	Wheelchair	
2595KL	Digital Chair Scale	Bluetooth	Chair	
2000KL	Digital In-Bed/Stretcher Scale	Bluetooth	Specialty	

*KG only models also available

Users have the option of connecting their scales to **Welch Allyn devices wirelessly with **Health o meter®** Professional Scales **Bluetooth** Wireless Solution. Options and connection type are version-specific for each scale model.

SCALE CONNECTIVITY KITS FOR THE FOLLOWING VITAL SIGNS DEVICES:

HEALTH O METER® WIRELESS CONNECTIVITY KIT	
<p>CONNEX Spot Monitor</p>  A white medical monitor with a digital display showing vital signs, mounted on a stand with a coiled cable.	<p>#BTKIT-01 Wireless Connectivity Kit Connex wireless module for installation into scale, wireless USB dongle (connects to device) and user instructions</p>
<p>CONNEX Vital Signs Monitor</p>  A white medical monitor with a digital display showing vital signs, mounted on a stand with a coiled cable.	<p>#BTKIT-01 Wireless Connectivity Kit Connex wireless module for installation into scale, wireless USB dongle (connects to device) and user instructions</p>
<p>CONNEX Integrated Wall System</p>  A white medical monitor with a digital display showing vital signs, mounted on a wall with various cables and connectors.	<p>#BTKIT-01 Wireless Connectivity Kit Connex wireless module for installation into scale, wireless USB dongle (connects to device) and user instructions</p>

ACQUISITION AND SET UP OF SCALE CONNECTIVITY:

- **Health o meter®** Professional scales and Connectivity Kits are sold through medical supply distributors in the U.S. and Canada.
- **Connex** Vital Signs Monitors require software version 1.7 or later to access the scale connectivity feature.
- **Connex** Spot Monitor devices require software version 1.2 or later to access the scale connectivity feature.
- **Connex** Vital Signs Monitors and **Connex** Spot Monitors are not compatible with scales that communicate via RS-232.
- Connectivity kits described above are compatible with new, current version scales. For legacy scales already in the field, contact **Health o meter®** Professional Customer Service at 1-800-815-6615 to confirm compatibility. Connectivity is also available on “BT” model scales.

WELCH ALLYN DEVICES FUNCTIONING WITH HEALTH O METER® PROFESSIONAL SCALES:

- **Health o meter®** Professional Digital Scales: Measure weight, allow for manual entry of height (600KL and 502KL will capture height digitally) and calculate BMI (where height is digitally captured or can be manually entered on the scale — exceptions are pediatric scales).
- **Welch Allyn** Connectivity Kits: Have the ability to transfer values for weight, height and BMI from a **Health o meter®** Professional Digital Scale to a **Welch Allyn** device.
- **Connex** Vital Signs Monitor devices with software version 1.7 or later and **Connex** Integrated Wall Systems, as well as **Connex** Spot Monitor devices with software version 1.2 or later: Accept and display weight, height and BMI when transmitted from the scale. In the Spot and Office profiles, BMI is calculated based on weight and height values that you enter manually on the device or that an attached weight scale transfers to the monitor.

EHR CONNECTIVITY:

EHR connectivity requires existing **Welch Allyn** device interface built by the specific EHR vendor for customer's specific EHR product and release. This interface must include each specific category (height, weight and BMI) for that information to transfer. Contact your connected device manager to research customer-specific scenarios.

FOR ADDITIONAL QUESTIONS REGARDING HEALTH O METER® PROFESSIONAL SCALES AND ASSOCIATED CONNECTIVITY KITS IN THE U.S. AND CANADA, PLEASE CONTACT:

Pelstar LLC / Health o meter® Professional Scales
9500 W. 55th St. McCook, IL 60525
Toll Free: +1 (800) 815-6615
Telephone: +1 (708) 377-0600
Fax: +1 (708) 377-0601
www.homscales.com

Rx only. For safe and proper use of the products mentioned herein, please refer to the Operators Manual or Instructions for Use.

Baxter.com

Baxter, Connex, Hillrom and Welch Allyn are trademarks of Baxter International, Inc. or its subsidiaries.

Bluetooth is a trademark of Bluetooth SIG, Inc.

Health o meter is a trademark of Sunbeam Products, Inc., and is used under license by Pelstar, LLC.

US-FLC45-230009 (v2.0) 06/2025