


Safe Accurate Accessible ADA Compliant Scales

Health o meter® Professional Scales ADA-compliant scales are thoughtfully engineered to support facilities in meeting accessibility standards. These scales comply with the latest Medical Diagnostic Equipment (MDE) standards for accessible diagnostic equipment and are the only ADA Compliant scales in the market that are independently certified by a third party.

Wheelchair Scale MDE Requirements:

- ① 32" wide by 40" deep platform
- ② 2" edge protection on platform
- ③ 2" edge protection on ramps
- ④ 1:8 running ramp slope



2101KL
Digital Platform Scale



REVIEWED FOR ADA COMPLIANCE
CERTIFICATION N° 5849459

1100KL
Digital Platform Scale with Extra-wide Handrails



REVIEWED FOR ADA COMPLIANCE
CERTIFICATION N° 5776764

Stand-on Scale MDE Requirements:

- ① Slip resistant platform
- ② Handrail supports on both sides of platform
- ③ 34" – 38" max. handrail height from platform



Learn More About our ADA Compliant Scales and the New MDE Standards