

Scales Compatible with Wireless Connectivity Kits

These scales have connectivity capabilities via various wireless kits. New scales can be purchased along with a kit, or kits can be purchased separately as after-market add-ons.* See kit descriptions below and compatibility information.

Select scales can be purchased with factory installed wireless modules. See available models on the following page.

Model Number**	Scale Description
498 Series	Remote Display Digital Scale
500 Series	Eye Level Digital Scale
502 Series	Primary Care Scale with Integrated Digital Height Rod
522 Series	Pediatric Digital Tray Scale
553 Series	Pediatric Digital Tray Scale
524KL	Pediatric Digital Seat Scale
565 Series	High-Resolution Digital Baby Scale
597 Series	Heavy Duty Eye Level Digital Scale
599 Series	Heavy Duty Waist High Digital Scale
600 Series	Heavy Duty Eye Level Digital Scale with Digital Height Rod
752 Series	Heavy Duty Remote Display Digital Scale
1100 Series	Digital Platform Scale
1110 Series	High-Capacity Remote Display Digital Scale with Large Base
2000 Series	Digital In-Bed/Stretchers Scale
2101 Series	Digital Platform Scale
2210KL-AM Series	Antimicrobial High-Resolution Digital Neonatal/Pediatric Tray Scale
2210KL4-AM Series	Antimicrobial High-Resolution Digital Neonatal/Pediatric 4-sided Tray
2400 Series	Digital Portable Wheelchair Scale
2500 Series	Digital Wheelchair Ramp Scale
2595 Series	Digital Chair Scale
2600 Series	Digital Wheelchair Ramp Scale
2610 Series	Digital Wheelchair Dual Ramp Scale
2620 Series	Digital Wheelchair Scale With Dual Ramps And Extra-Large Platform, ADA Compliant
2650 Series	Digital Dual Ramp Wheelchair Scale
2700 Series	Digital Dual Ramp Wheelchair Scale with X-Large Platform
2900KL-AM Series	Antimicrobial Wall-Mounted Fold-Up Wheelchair Scale
2920KL-AM Series	Antimicrobial Wall-Mounted Fold-Up Wheelchair Scale, ADA Compliant
3001KL-AM Series	Antimicrobial Digital Platform Scale
3001KL-AMX Series	Antimicrobial Digital Platform Scale with Extended Handrails
3105KL-AM Series	Heavy-Duty Antimicrobial Platform Scale

Connectivity Kit Descriptions

- BTKIT-01 Wireless Upgrade Kit for Use with Welch Allyn® Monitors
 - Compatible Welch Allyn devices are Connex Spot Vital Signs Monitors, Connex 6000 Vital Signs 6000 Series, and Connex Integrated Wall Systems
- BTKIT-03 Wireless Upgrade Kit for Use with Pre-Approved Devices and Windows® PCs (USB dongle included)
 - Purchase BTKIT-03 for connection via USB dongle inserted into the PC/Device
- BTKIT-100 Wireless Upgrade Kit for Use with Pre-Approved Devices and Windows® Applications (USB dongle not included)
 - Purchase BTKIT-100 for connection via the device's internal Bluetooth
- BTPRINTERKIT-01 Wireless Thermal Printer Kit
 - Kit includes components to allow scale to communicate with thermal wireless printer. (Note: 597, 599, 600, and 752 series scales are not compatible with this kit)

*Kits must be ordered as a separate line item with a new scale, or the kit can be bought as an after-market add-on.

**List of compatible scales refers to new, current version scales. For legacy scales already in the field, see page 3 for compatibility information or contact Customer Service at 1-800-815-6615.

Scales with Built-in Pelstar Wireless Technology

These scales have a factory-installed internal wireless module; no external hardware needs to be attached to the scale. The additional included wireless USB Dongle is used to transmit data to a Welch Allyn device (Welch Allyn Connex Spot Monitors, Connex Vital Signs 6000 Series, and Connex Integrated Wall Systems).

Model Number	Scale Description
500KL-BT / 500KG-BT	Eye Level digital scale with Built-in Pelstar Wireless Technology
502KL-BT / 502KG-BT	Primary Care Scale with Integrated Digital Height Rod with Built-in Pelstar Wireless Technology
600KL-BT / 600KG-BT	Eye Level Scale with Digital Height Rod and Built-in Pelstar Wireless Technology
1100KL-BT / 1100KG-BT	Digital Platform Scale with Built-in Pelstar Wireless Technology
1110KL-BT / 1110KG-BT	High-Capacity Remote Display Digital Scale with Large Base with Built-in Pelstar Wireless Technology
2000KL-BT / 2000KG-BT	Digital In-Bed/Stretchers Scale with Built-in Pelstar Wireless Technology
2101KL-BT / 2101KG-BT	Digital Platform Scale with Built-in Pelstar Wireless Technology
3105KL-AM-BT / 3105KG-AM-BT	Heavy-Duty Antimicrobial Platform Scale with Pelstar Wireless Technology
2400KL-BT / 2400KG-BT	Digital Portable Wheelchair Scale with Built-in Pelstar Wireless Technology
2500KL-BT / 2500KG-BT	Digital Wheelchair Ramp Scale with Built-in Pelstar Wireless Technology
2600KL-BT / 2600KG-BT	Digital Wheelchair Ramp Scale with Built-in Pelstar Wireless Technology
2610KL-BT / 2610KG-BT	Digital Wheelchair Scale with Dual Ramps and Built-in Pelstar Wireless Technology
2650KL-BT / 2650KG-BT	Digital Wheelchair Dual Ramp Scale with Built-in Pelstar Wireless Technology
2700KL-BT / 2700KG-BT	Digital Wheelchair Dual Ramp Scale with Extra Large Platform and Built-in Pelstar Wireless Technology
2900KL-AM-BT / 2900KG-AM-BT	Antimicrobial Wall-Mounted Fold-Up Wheelchair Scale with Pelstar Wireless Technology
3001KL-AM-BT / 3001KG-AM-BT	Antimicrobial Digital Platform Scale with Pelstar Wireless Technology
3001KL-AMX-BT / 3001KG-AMX-BT	Antimicrobial Digital Platform Scale with Extended Handrails with Pelstar Wireless Technology
2210KL-AM-BT / 2210KG-AM-BT	Antimicrobial High-Resolution Digital Neonatal/Pediatric Tray Scale with Built-in Pelstar Wireless Technology
2210KL4-AM-BT / 2210KG4-AM-BT	Antimicrobial High-Resolution Digital Neonatal/Pediatric 4-sided Tray with Built-in Pelstar Wireless Technology

Legacy Scales Compatible with Wireless Connectivity Kits

Select scales already in use in the field may be compatible. Find the date code on your scale's product label to determine connectivity kit compatibility. Date code column indicates the versions of the scale with kit compatibility.

Model Number**	Scale Description	Date Code
498 Series	Remote Display Digital Scale	4123, or with the last two digits of 24 or higher
500 Series	Eye Level Digital Scale	500KL: Last two digits of 20 or higher; 500KG: last two digits of 21 or higher
502 Series	Primary Care Scale with Integrated Digital Height Rod	Last two digits of 19 or higher
522 Series	Pediatric Digital Tray Scale	3423, 3723, or with the last two digits of 24 or higher
553 Series	Pediatric Digital Tray Scale	2224, 2724, 3224, 3324, 3824, 4224, 4724, 5124, or with the last two digits of 25 or higher
524KL	Pediatric Digital Seat Scale	3723, or with the last two digits of 24 or higher
565 Series	High-Resolution Digital Baby Scale	Last two digits of 19 or higher
597 Series	Heavy Duty Eye Level Digital Scale	Last two digits of 19 or higher. Not compatible with BTPRINTINTERKIT-01
599 Series	Heavy Duty Waist High Digital Scale	Last two digits of 19 or higher. Not compatible with BTPRINTINTERKIT-01
600 Series	Heavy Duty Eye Level Digital Scale with Digital Height Rod	Last two digits of 19 or higher. Not compatible with BTPRINTINTERKIT-01
752 Series	Heavy Duty Remote Display Digital Scale	Last two digits of 19 or higher. Not compatible with BTPRINTINTERKIT-01
1100 Series	Digital Platform Scale	Last two digits of 19 or higher
1110 Series	High-Capacity Remote Display Digital Scale with Large Base	Last two digits of 19 or higher
2000 Series	Digital In-Bed/Stretcher Scale	Last two digits of 19 or higher
2101 Series	Digital Platform Scale	Last two digits of 19 or higher
2210KL-AM Series	Antimicrobial High-Resolution Digital Neonatal/Pediatric Tray Scale	Last two digits of 19 or higher
2210KL4-AM Series	Antimicrobial High-Resolution Digital Neonatal/Pediatric 4-sided Tray	Last two digits of 19 or higher
2400 Series	Digital Portable Wheelchair Scale	Last two digits of 19 or higher
2500 Series	Digital Wheelchair Ramp Scale	Last two digits of 19 or higher
2595 Series	Digital Chair Scale	Last two digits of 25 or higher
2600 Series	Digital Wheelchair Ramp Scale	Last two digits of 19 or higher
2610 Series	Digital Wheelchair Dual Ramp Scale	Last two digits of 19 or higher
2620 Series	Digital Wheelchair Scale With Dual Ramps And Extra-Large Platform, ADA Compliant	Last two digits of 19 or higher
2650 Series	Digital Dual Ramp Wheelchair Scale	Last two digits of 19 or higher
2700 Series	Digital Dual Ramp Wheelchair Scale with X-Large Platform	Last two digits of 19 or higher
2900KL-AM Series	Antimicrobial Wall-Mounted Fold-Up Wheelchair Scale	Last two digits of 19 or higher
2920KL-AM Series	Antimicrobial Wall-Mounted Fold-Up Wheelchair Scale, ADA Compliant	Last two digits of 19 or higher
3001KL-AM Series	Antimicrobial Digital Platform Scale	Last two digits of 19 or higher
3001KL-AMX Series	Antimicrobial Digital Platform Scale with Extended Handrails	Last two digits of 19 or higher
3105KL-AM Series	Heavy-Duty Antimicrobial Platform Scale	Last two digits of 19 or higher

Date code can be found on the product label on the back of the display head.

