

THE 2920KL-AM

ANTIMICROBIAL WALL-MOUNTED FOLD-UP WHEELCHAIR SCALE, CERTIFIED FOR ADA COMPLIANCE

As North America's leading medical scale manufacturer, we proudly offer a Certified for Americans with Disabilities Act (ADA) compliance, antimicrobial, wall-mounted wheelchair scale designed for optimal performance in demanding healthcare environments. Engineered to support up to 1,000 pounds (454 kilograms), this scale combines exceptional durability, reliability, and user-friendly operation.

Space-Saving Design

- Wall-mounted with foldable platform
- Extends only 4.7 inches from the wall when folded
- Display can be mounted on either side for flexibility
- Ideal for hospitals, clinics, acute, long-term care, and doctors' offices

Caregiver & Patient Friendly Platform

- Spacious and easy to lift
- Fits most wheelchairs for safe, comfortable weighing

Enhanced Safety Features

- Antimicrobial surface protection
- Raised platform edge
- Recessed handle
- No-pinch-point design



The scale provides easy integration and seamless communication with today's EMR applications. The 2920KL-AM includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

The 2920KL-AM is designed to meet the US Access Board's Standards. Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which are enforced by the Department of Justice.

- 32" wide x 40" deep platform
- 2" edge protection on platform
- 2" edge protection on ramps
- 1:8 running slope



Health o meter®
Professional

Health o meter® Professional 2920KL-AM

Reviewed for
**ADA compliance and
Platinum certified by Corada**

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that the 2920KL-AM provides a safe, accessible weighing experience.



100% CERTIFIED FOR ADA COMPLIANCE

2" EDGE PROTECTION

32" X 40" LOW PROFILE PLATFORM

1:8 RUNNING SLOPE

#1 TRUSTED BRAND

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm
Screen Size: 2 3/4" x 1 1/2" High-Contrast Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Pre-tare
Hold/release
Reweigh
Recall
Auto zero
Auto off
Time/date
Variable auto off time
Audible/mute sound option
Parent/child weighing function
Everlock® - permanent unit lock

Platform Size (w x d x h)

40" x 32" x 1 3/5"
(1016 mm x 819 mm x 38 mm)

Ramp Size (w x d)

12 1/4" x 32"
(318 mm x 813 mm)

Product Footprint Unfolded (w x d x h)

77 1/2" x 37 7/8" x 37"
(1968 mm x 962 mm x 940 mm)

Product Footprint Folded (w x d x h)

77 1/2" x 4 3/4" x 37"
(1968 mm x 114 mm x 940 mm)

Product Weight

118 lb / 53.5 kg

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
2920KL-AM	Antimicrobial wall-mounted fold-up wheelchair scale, ADA Compliant, LB & KG
2920KG-AM	Antimicrobial wall-mounted fold-up wheelchair scale, ADA Compliant, KG only
SS-2920KL-AM	ScaleSurance® 2 year extended warranty for 2920 series*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscapes.com/ConnectivitySolutions



Easy to Lift Platform with a Recessed Handle

To learn more, call **1-800-815-6615**
or visit homscapes.com/2920KL-AM



*Not available in all countries.

Health o meter®
Professional