

Health o meter®



Professional

Safe Accurate Accessible ADA Compliant Scales

Health o meter® Professional Scales ADA compliant scales are thoughtfully engineered to support facilities in meeting accessibility standards. These scales comply with the latest Medical Diagnostic Equipment (MDE) standards for accessible diagnostic equipment and are reviewed for ADA Compliance and Platinum Certified by Corada.

Wheelchair Scale MDE Requirements:

- 1 32" wide by 40" deep platform
- 2 2" edge protection on platform
- 3 2" edge protection on ramps
- 4 1:8 running ramp slope



Stand-on Scale MDE Requirements:

- 1 Slip resistant platform
- 2 Handrail supports on both sides of platform (*Handrail along length of platform for wheelchair scales*)
- 3 34" – 38" max. handrail height from platform



Learn More About our Certified ADA Compliant Scales and the New MDE Standards

