

THE 2622ADAKIT

ADA WHEELCHAIR & STAND-ON COMPLIANCE KIT FOR WHITE 2600 AND 2610 SERIES OF SCALES, CERTIFIED FOR ADA COMPLIANCE

The Health o meter® Professional 2622ADAKIT is an optional ADA Compliance Kit for the white 2600 series of scales. The kit converts a Health o meter® Professional 2600 series wheelchair scale into a dual purpose ADA compliant stand-on and wheelchair scale. Certified for Americans with Disabilities Act (ADA) compliance, the 2622ADAKIT includes new ramps, new feet, new antimicrobial handrail and platform and ramp edge guards with anti-slip caution tape on the ramps for enhanced safety and accessibility. The upgrade also replaces the existing handrail and fixed display head with a wall-mounted display head. It's easy to install on any existing white 2600KL and 2610KL series of scales, making it an ideal choice for doctor's offices, hospitals, clinics, and acute or long-term care facilities looking to meet the latest U.S. Access Board Standards for both stand-on and wheelchair scales.

Kit Includes:

- 2 ramps with edge guards
- platform edge guard
- handrail
- wall mount bracket for display head
- load cell connection cable
- required hardware



The 2622ADAKIT is designed to meet the US Access Board's Standards. Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which are enforced by the Department of Justice.

- 32" wide x 40" deep platform
- 2" edge protection on platform
- 2" edge protection on ramps
- 1:8 running slope
- Slip-resistant platform
- Handrail support along the length of the platform
- 34" – 38" max. handrail height



ORDERING INFORMATION

MODEL	DESCRIPTION
2622ADAKIT	ADA Wheelchair & Stand-on Compliance Kit for White 2600 and 2610 Series

*For assistance with compatible scales, contact Customer Service at 1-800-815-6615.

Health o meter®

Professional

Health o meter® Professional 2622ADAKIT

Reviewed for
**ADA compliance and
Platinum certified by Corada**



TECHNICAL SPECIFICATIONS

Edge Guard Size

40" x 1/8" x 2"
(1016 mm x 3 mm x 51 mm)

Ramp Size

32 1/4" x 16"
(819 mm x 406 mm)

Platform Size with 2600KL

40" x 32 1/4" x 1 7/8"
(1016 mm x 819 mm x 48 mm)

Product Footprint with 2600KL

77 1/2" x 45 1/8" x 39 1/2"
(1816 mm x 1147 mm x 1003mm)

Product Footprint with 2600KL (Folded)

77 1/2" x 9 1/2" x 44"
(1816 mm x 241 mm x 1118 mm)

Product Weight with 2600KL

126 lb / 57 kg

Handrail Height

37.5" / 952 mm

Warranty

2-year limited manufacturer warranty

2622ADAKIT Compatible Scales *

2600KL
2600KG
2600KL-BT
2600KG-BT
2610KL
2610KG
2610KL-BT
2610KG-BT

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/2622ADAKIT**