

Product Specification Guide

Stand-On Scales

	499KL	500KL	502KL	599KL	600KL
					
FEATURES					
Capacity	460 lb / 209 kg	660 lb / 300 kg	660 lb / 300 kg	660 lb / 300 kg	660 lb / 300 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Height Rod Included		•	•		•
Display Size	1" / 25 mm	1" / 25 mm	1" / 25 mm	1" / 25 mm	1" / 25 mm
Wheels	•	•	•	•	•
EMR Capable		Optional	Optional	•	•
Power Adapter Included				•	•
Batteries Included	•	•	•		
Warranty (years)	2	2	2	2	2
Wireless Connectivity		Optional	Optional	Optional	Optional
PRODUCT DIMS (w x d x h)					
Platform Size	13 ¾" x 16 ½" x 2 ⅝"	13 ¾" x 16 ½" x 2 ⅝"	13 ¾" x 16 ½" x 2 ⅝"	14 ¼" x 14 ¼" x 2 ⅝"	14 ¼" x 14 ¼" x 2 ⅝"
Product Footprint	14 ⅞" x 21 ⅞" x 43 ½"	14 ⅞" x 21 ⅞" x 53 ⅝"	14 ⅞" x 21 ⅞" x 53 ⅝"	14 ¾" x 23 ¼" x 42 ½"	14 ¾" x 24 ⅞" x 54 ⅜"
Product Weight	21 lb / 9 kg	21 lb / 9 kg	22.4 lb / 10.2 kg	20 lb / 9 kg	28 lb / 13 kg
FUNCTIONS					
LB/KG Conversion	•	•	•	•	•
LB/KG Lockout	•	•	•	•	•
Hold/ Release	•	•	•	•	•
BMI	•	•	•	•	•
Zero	•	•	•	•	•
Tare	•	•	•	•	•
Pretare		•	•		
Auto Zero	•	•	•	•	•
Auto Off	•	•	•	•	•
Reweigh		•	•		
Recall		•	•		
Everlock®	•	•	•		
Parent/ Child Weighing	•	•	•	•	•
Variable Auto Off Time	•	•	•		
Audible/ Mute Sound	•	•	•		
Enabling Unit View		•	•		
ACCESSORIES (Optional)					
Power Adapter	ADPT31	ADPT31	ADPT31	ADPT50	ADPT50
Height Rod	499KLROD			499KLROD	
Handrail	499HB	500HB	500HB		
Wireless Kit		•	•	•	•
Wireless Printer		•	•		

	1100KL	1110KL	2101KL	3001KL-AM	3001KL-AMX	3105KL-AM
						
FEATURES						
Capacity	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Display Size	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm
Wheels	•		•	•	•	•
EMR Capable	•	•	•	•	•	•
Power Adapter Included	•	•	•	•	•	
Batteries Included						•
Warranty (years)	2	2	2	2	2	2
Antimicrobial				•	•	•
Antimicrobial Keypad	•	•	•	•	•	•
Wireless Connectivity	Optional	Optional	Optional	Optional	Optional	Optional
ADA Compliant	•		•			
PRODUCT DIMS (w x d x h)						
Platform Size	22" x 15 ¾" x 2 ½"	22" x 15 ¾" x 2 ⅝"	26" x 22" x 2 ⅝"	22" x 15 ¾" x 2 ⅝"	22" x 15 ¾" x 2 ⅝"	19 ¾" x 21 ¼" x 2 ⅝"
Product Footprint	31 ¾" x 33 ½" x 49 ¾"	22" x 23 ½" x 4 ¾"	31" x 27" x 50"	22" x 34 ½" x 52"	22" x 34 ½" x 52"	23 ¼" x 28" x 52 ¼"
Product Weight	52 lb / 23 kg	27 lb / 12 kg	66 lb / 30 kg	42 lb / 19 kg	43 lb / 19.5 kg	86 lb / 39 kg
FUNCTIONS						
LB/KG Conversion	•	•	•	•	•	•
LB/KG Lockout	•	•	•	•	•	•
Hold/ Release	•	•	•	•	•	•
BMI	•	•	•	•	•	•
Zero	•	•	•	•	•	•
Tare	•	•	•	•	•	•
Pretare	•	•	•	•	•	•
Auto Zero	•	•	•	•	•	•
Auto Off	•	•	•	•	•	•
Reweigh	•	•	•	•	•	•
Recall	•	•	•	•	•	•
Everlock®	•	•	•	•	•	•
Parent/ Child Weighing	•	•	•	•	•	•
Time/ Date	•	•	•	•	•	•
Variable Auto Off Time	•	•	•	•	•	•
Audible/ Mute Sound	•	•	•	•	•	•
Enabling Unit View	•	•	•	•	•	•
ACCESSORIES (Optional)						
Power Adapter	ADPT30	ADPT30	ADPT30	ADPT30	ADPT30	ADPT30
Height Rod	PROPLUSROD or 245EHR-3105	201HR-1110 or 245EHR-1110	STROD or 245EHR-2101	3001HR or 245EHR-3001		201HR-3105 or 245EHR-3105
Wireless Kit	•	•	•	•	•	•
Wireless Printer	•	•	•	•	•	•

*2622KL also meets stand-on scale ADA requirements — see Wheelchair & Chair Scales for full specifications.

	349KLX	498KL	752KL	2842KL
--	--------	-------	-------	--------



FEATURES

Capacity	400 lb / 180 kg	660 lb / 300 kg	660 lb / 300 kg	600 lb / 270 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Display Size	1" / 25 mm	1" / 25 mm	1" / 25 mm	¾" / 22 mm
EMR Capable	Optional		•	
Power Adapter Included			•	•
Batteries Included	•	•		
Warranty (years)	2	2	2	2
Wireless Connectivity	Optional		Optional	

PRODUCT DIMS (w x d x h)

Platform Size	12 ½" x 12" x 1 ⅞"	13 ¾" x 16 ½" x 2 ⅜"	14 ¼" x 14 ¼" x 2 ⅝"	22 ¼" x 42" x 2 ½"
Product Footprint	12 ½" x 12" x 1 ⅞"	13 ¾" x 16 ½" x 2 ⅜"	14 ½" x 14 ½" x 2 ¾"	22 ¼" x 42" x 2 ½"
Product Weight	10 lb / 4.5 kg	14 lb / 6.4 kg	14 lb / 6.4 kg	53 lb / 24 kg

FUNCTIONS

LB/KG Conversion	•	•	•	•
LB/KG Lockout	•	•	•	•
Hold/ Release	•	•	•	•
BMI		•	•	
Zero	•	•	•	•
Tare	•	•	•	•
Pretare		•		•
Auto Zero	•	•	•	
Auto Off	•	•	•	•
Reweigh	•	•		•
Recall				•
Everlock®	•	•		•
Parent/ Child Weighing	•	•	•	
Variable Auto Off Time	•	•		•
Audible/ Mute Sound		•		
Enabling Unit View	•	•		•

ACCESSORIES (Optional)

Power Adapter	ADPT40	ADPT31	ADPT50	ADPT31
Carrying Case	64771	64771	64771	
Wireless Kit		•	•	
Wireless Printer		•		

Mechanical Stand-On Scales

	400KL	402KL	402LB	402KG SPECIAL ORDER ONLY	450KL
					
FEATURES					
Capacity	390 lb / 180 kg	390 lb / 180 kg	400 lb	180 kg	500 lb / 200 kg
Graduation	¼ lb / 100 g	¼ lb / 100 g	¼ lb / 100 g	¼ lb / 100 g	¼ lb / 100 g
Height Rod Included		•	•	•	•
Rotating Poise Bars					•
Engraved Numerals					•
Warranty (years)	10	10	10	10	10
PRODUCT DIMS (w x d x h)					
Platform Size	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"
Product Footprint	20 ¼" x 10 ⅜" x 57 ½"	20 ¼" x 10 ⅜" x 57 ½"	20 ¼" x 10 ⅜" x 57 ½"	20 ¼" x 14" x 57 ½"	20 ¼" x 14" x 57 ½"
Product Weight	24 lb / 11 kg	28 lb / 13 kg	28 lb / 13 kg	28 lb / 13 kg	31 lb / 14 kg
ACCESSORIES (Optional)					
Wheels	55000	55000	55000	55000	55000
Height Rod	402KLROD				
Counterweights	55090	55090			Included

Pediatric Scales


	522KL	524KL	549KL	550KL	553KL
					
FEATURES					
Capacity	50 lb / 23 kg	50 lb / 23 kg	44 lb / 20 kg	44 lb / 20 kg	44 lb / 20 kg
Resolution	0-20 lb / 0.2 oz; 20-50 lb / 0.5 oz; 0-9 kg / 5 g; 9-23 kg / 10 g	0-20 lb / 0.2 oz; 20-50 lb / 0.5 oz; 0-9 kg / 5 g; 9-23 kg / 10 g	0.5 oz / 10 g	0.5 oz / 10 g	0.5 oz / 10 g
Measuring Tape Included	•				•
Display Size	1" / 25 mm	1" / 25 mm	5⁄8" / 15 mm	1⁄2" / 13 mm	1" / 25 mm
EMR Capable	Optional	Optional			Optional
Batteries Included	•	•	•	•	•
Warranty (years)	2	2	2	2	2
Wireless Connectivity	Optional	Optional			Optional
PRODUCT DIMS (w x d x h)					
Tray/ Seat Size	24 1⁄8" x 14 1⁄2" x 2 5⁄8"	24 1⁄4" x 14 1⁄4" x 18 1⁄2"	19 3⁄8" x 12 1⁄8" x 2"	23 1⁄2" x 8" x 2"	24 7⁄8" x 13" x 1 3⁄8"
Product Footprint	24 1⁄8" x 21" x 23 1⁄8"	24 1⁄8" x 21" x 23 1⁄8"	19 3⁄8" x 13 1⁄2" x 4 3⁄8"	23 1⁄2" x 13 1⁄4" x 3 5⁄8"	24 7⁄8" x 16 1⁄2" x 3 7⁄8"
Product Weight	20 lb / 9 kg	21 lb / 10 kg	4 lb / 2 kg	8 lb / 4 kg	10 lb / 5 kg
FUNCTIONS					
LB/KG Conversion	•	•	•		•
LB/KG Lockout	•	•	•		•
Hold/ Release	•	•			•
Zero	•	•			•
Tare	•	•	•	•	•
Pretare	•	•			•
Auto Zero	•	•	•		•
Auto Off	•	•	•	•	•
Reweigh	•	•			•
Recall	•	•			•
Everlock®	•	•	•		•
Variable Auto Off Time	•	•	•		•
Audible/ Mute Sound	•	•			
Enabling Unit View	•	•	•		•
Motion Detection	•	•	•		•
Auto Weight Hold	•	•	•	•	•
ACCESSORIES (Optional)					
Power Adapter	ADPT31	ADPT31			ADPT31
Carrying Case			64769	553CASE	553CASE
Height Rod	522EHR or 522HR			eBABYHR + 550EHRBRACKET	553HR or 553EHR
Cart			2210CART	2210CART	2210CART
Wireless Kit	•	•			•
Wireless Printer	•	•			•

	565KL	2210KL-AM	2210KL4-AM	8320KL
				
FEATURES				
Capacity	44 lb / 20 kg	45 lb / 20 kg / 20,000 gm	45 lb / 20 kg / 20,000 gm	33 lb / 15 kg
Resolution	0 – 44 lb x 0.1 oz / 0 – 6,000 g x 1 g / 6,001 g – 20,000 g x 2 g	1 g (to 10,000 g); 5 g (10,001 to 20,000 g)	1 g (to 10,000 g); 5 g (10,001 to 20,000 g)	0.5 oz/ 10 g
Measuring Tape Included	•	•	•	
Display Size	¾" / 19 mm	1 ⅝" / 28 mm	1 ⅝" / 28 mm	¾" / 19 mm
EMR Capable	Optional	•	•	
Batteries Included	•	•	•	•
Warranty (years)	2	2	2	2
Antimicrobial		•	•	
Wireless Connectivity	Optional	Optional	Optional	
PRODUCT DIMS (w x d x h)				
Tray/ Seat Size	24 ⅞" x 9 ⅞" x 2 ¾"	24" x 15 ¼" x 2 ¾"	26 ½" x 15 ¼" x 2 ¾"	23 ½" x 11" x ⅞"
Product Footprint	26 ¼" x 13" x 5 ¼"	24 ¾" x 17 ¾" x 6"	26 ½" x 17 ¾" x 6"	29 ⅞" x 12 ⅝" x 2 ⅝"
Product Weight	9 lb / 4.1 kg	28 lb / 13 kg	28 lb / 13 kg	5 lb / 3 kg
FUNCTIONS				
LB/KG Conversion	•	•	•	•
LB/KG Lockout	•	•	•	
Hold/ Release	•	•	•	
Zero	•	•	•	•
Tare	•	•	•	
Pretare	•	•	•	
Auto Zero	•	•	•	•
Auto Off	•	•	•	•
Reweigh	•	•	•	
Recall	•	•	•	
Everlock®	•	•	•	
Variable Auto Off Time	•	•	•	
Audible/ Mute Sound	•	•	•	
Enabling Unit View	•	•	•	
Motion Detection	•	•	•	
Auto Weight Hold	•	•	•	
Breast Milk Intake		•	•	
ACCESSORIES (Optional)				
Power Adapter	69-00001	ADPT30	ADPT30	
Carrying Case				64771
Height Rod	565EHR	2210EHR		
Cart	2210CART	2210CART-AC	2210CART-AC	
Wireless Kit	•	•	•	
Wireless Printer	•	•	•	

	1522KL	1524KL	386S-01	386KGS-01	BABYSLINGSCALE
					
FEATURES					
Capacity	135 lb / 62 kg	130 lb / 65 kg	50 lb	25 kg	25 kg
Graduation	1 oz / 20 g	1 oz / 20 g	¼ lb	100 g	100 g
Measuring Tape Included	•		•	•	
Dial Size			3" / 76 mm	3" / 76 mm	5 ½" / 140 mm
Warranty (years)	2	2	1	1	1
PRODUCT DIMS (w x d x h)					
Tray/ Seat Size	24" x 16 ½" x 2 ⅝"	24 ½" x 14 ¼" x 18 ½"	19 ⅜" x 12 ⅝" x 3 ¾"	19 ⅜" x 12 ⅝" x 3 ¾"	
Product Footprint	24" x 20" x 22 ½"	24 ½" x 21" x 22 ½"	19 ⅜" x 14 ½" x 5 ⅛"	19 ⅜" x 14 ½" x 5 ⅛"	12" x 43" x 1 ⅞"
Product Weight	25 lb / 11 kg	32 lb / 15 kg	10 lb / 5 kg	10 lb / 5 kg	2 lb / 1 kg

Wheelchair & Chair Scales

	2500KL	2600KL	2610KL	2620KL	2622KL
					
FEATURES					
Capacity	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Display Size	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm
Wheels	•	•	•	•	•
EMR Capable	•	•	•	•	•
Power Adapter Included	•	•	•	•	•
Warranty (years)	2	2	2	2	2
Antimicrobial Keypad	•	•	•	•	•
Wireless Connectivity	Optional	Optional	Optional	Optional	Optional
ADA Compliant		Optional	Optional	•	•
PRODUCT DIMS					
Platform Size	29 ¾" x 28 ¼" x 2 ¾"	36" x 32 ¼" x 2 ½"	36" x 32 ¼" x 2 ½"	40" x 32 ¼" x 1 ⅞"	40" x 32 ¼" x 1 ⅞"
Product Footprint	40 ¼" x 34 ½" x 45 ½"	48 ⅞" x 43 ⅞" x 49 ⅞"	56 ¾" x 43 ⅞" x 49 ⅞"	77 ¼" x 43 ⅞" x 48 ¾"	71 ½" x 45 ⅞" x 39 ½"
Product Weight	81 lb / 37 kg	93 lb / 42 kg	98 lb / 44 kg	126 lb / 57 kg	126 lb / 57 kg
FUNCTIONS					
LB/ KG Conversion	•	•	•	•	•
LB/ KG Lockout	•	•	•	•	•
Hold/ Release	•	•	•	•	•
BMI	•	•	•	•	•
Zero	•	•	•	•	•
Tare	•	•	•	•	•
Pretare	•	•	•	•	•
Auto Zero	•	•	•	•	•
Auto Off	•	•	•	•	•
Reweigh	•	•	•	•	•
Recall	•	•	•	•	•
Everlock®	•	•	•	•	•
Time/ Date	•	•	•	•	•
Parent/ Child Weighing	•	•	•	•	•
Variable Auto Off Time	•	•	•	•	•
Audible/ Mute Sound	•	•	•	•	•
Enabling Unit View	•	•	•	•	•
ACCESSORIES (Optional)					
Power Adapter	ADPT30	ADPT30	ADPT30	ADPT30	ADPT30
Wireless Kit	•	•	•	•	•
Wireless Printer	•	•	•	•	•
Ramp	B2500RAMP	B2600RAMP			
ADA Compliance Upgrade Kit		2622ADAKIT, 2600ADAKIT or 2600ADAKIT-WHITE	2622ADAKIT, 2600ADAKIT or 2600ADAKIT-WHITE	2622HB	
Transport Kit					43-0070

	2400KL	2650KL	2700KL	2900KL-AM	2920KL-AM	
FEATURES						
	Capacity	800 lb / 363 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg
	Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
	Display Size	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm
	Wheels	•	•	•		
	EMR Capable	•	•	•	•	•
	Power Adapter Included	•	•	•	•	•
	Warranty (years)	2	2	2	2	2
	Antimicrobial				•	•
	Antimicrobial Keypad	•	•	•	•	•
Wireless Connectivity	Optional	Optional	Optional	Optional	Optional	
ADA Compliant					•	
PRODUCT DIMS (w x d x h)						
Platform Size	25 ¾" - 50 ⅞"	31 ½" x 31 ½" x 2"	43" x 42" x 2 ½"	29 ½" x 32 ½" x 2 ¾"	40" x 32" x 1 ⅝"	
Product Footprint	51 ¼" x 40 ½" x 3 ½"	51 ¾" x 33 ⅝" x 5 ¼"	69 ⅞" x 51 ½" x 3 ⅜"	58" x 37" x 37"	77 ½" x 37 ⅞" x 37"	
Product Weight	40 lb / 18 kg	82 lb / 37 kg	187 lb / 85 kg	91 lb / 41 kg	118 lb / 53.5 kg	
FUNCTIONS						
LB/ KG Conversion	•	•	•	•	•	
LB/ KG Lockout	•	•	•	•	•	
Hold/ Release	•	•	•	•	•	
BMI	•	•	•	•	•	
Zero	•	•	•	•	•	
Tare	•	•	•	•	•	
Pretare	•	•	•	•	•	
Auto Zero	•	•	•	•	•	
Auto Off	•	•	•	•	•	
Reweigh	•	•	•	•	•	
Recall	•	•	•	•	•	
Everlock®	•	•	•	•	•	
Time/ Date	•	•	•	•	•	
Parent/ Child Weighing	•	•	•	•	•	
Variable Auto Off Time	•	•	•	•	•	
Audible/ Mute Sound	•	•	•	•	•	
Enabling Unit View	•	•	•	•	•	
ACCESSORIES (Optional)						
Power Adapter	ADPT30	ADPT30	ADPT30	ADPT30	ADPT30	
Wireless Kit	•	•	•	•	•	
Wireless Printer	•	•	•	•	•	

	2500CKL	2595KL	594KL	445KL
FEATURES				
	Capacity	1000 lb / 454 kg	600 lb / 270 kg	600 lb / 270 kg
	Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
	Display Size	1 ½" / 38 mm	1 ⅝" / 30 mm	1 " / 25 mm
	Wheels	•	•	•
	EMR Capable	•		
	Power Adapter Included	•	•	
	Warranty (years)	2	2	2
	Antimicrobial Keypad	•		
	Wireless Connectivity	Optional	Optional	
PRODUCT DIMS (w x d x h)				
Platform Size	29 ⅝" x 28 ¼" x 2 ⅝"			
Seat Size	18" x 15"	18 ¼" x 15 ¾"	18" x 18 ½"	18" x 17 ¼"
Product Footprint	40 ¼" x 34 ½" x 45 ½"	23 ½" x 37 ½" x 39 ⅞"	23 ¼" x 40 ⅞" x 37 ⅞"	24 ¾" x 31 ⅝" x 39 ⅞"
Product Weight	121 lb / 55 kg	84 lb / 38 kg	39 lb / 18 kg	60 lb / 27 kg
FUNCTIONS	LB/ KG Conversion	•	•	•
	LB/ KG Lockout	•	•	•
	Hold/ Release	•	•	•
	BMI	•		•
	Zero	•	•	•
	Tare	•	•	•
	Pretare	•	•	
	Auto Zero	•		•
	Auto Off	•	•	•
	Reweigh	•	•	
	Recall	•	•	
	Everlock®	•	•	
	Time/ Date	•		
	Parent/ Child Weighing	•		
	Variable Auto Off Time	•	•	
	Audible/ Mute Sound	•		
	Enabling Unit View	•	•	
ACCESSORIES (Optional)				
Power Adapter	ADPT30	ADPT30	ADPT30	
Wireless Kit	•	•		
Wireless Printer	•			
Ramp	B2500RAMP			

Body Composition Scale

BCS-G7

FEATURES

Capacity	1000 lb / 454 kg
Resolution	0.2 lb / 0.1 kg
Body Fat Accuracy	+/- 2%
Measurement Time	60 seconds
Power Adapter Included	•
6 Window Display	•
Limited Test Restrictions	•
OHMS Modulation Logic	•
Integrated Software	•
Warranty (years)	1

DISPLAY MEASUREMENTS

Weight	•
Body Fat %	•
Body Water %	•
Muscle Mass %	•
BMI	•
REE (Resting Energy Expenditure)/ BMR (Basal Metabolic Rate)	•

IMPACT ON ACCURACY/ PRE-TEST RESTRICTIONS

Recent Food/ Drink	OK to Test
Recent Exercise	OK to Test
Full Night Sleep	OK to Test
Recent Alcohol Intake	OK to Test
Full Bladder	OK to Test

PRODUCT DIMS (w x d x h)

Platform Size	14" x 16" x 2 7/8"
Product Footprint	20" x 18" x 41"
Product Weight	28 lb / 12 kg



Professional Home Care Scales

	800KL	810KL	880KL	844KL	822KL	894KLT/ 895KLT
						
FEATURES						
Capacity	397 lb / 180 kg	440 lb / 200 kg	500 lb / 220 kg	440 lb / 200 kg	400 lb / 180 kg	894KLT: 400 lb / 180 kg 895KLT: 550 lb / 250 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.1 lb / 0.05 kg	0.1 lb / 0.05 kg	0.2 lb / 0.1 kg
Display Size	1 ½" / 38 mm	1 ⅞" / 48 mm	1 ⅞" / 48 mm	1 ¼" / 32 mm	1 ½" / 38 mm	1 ½" / 38 mm
Batteries Included	•	•	•	•	•	•
Case Pack Quantity	4	1	4	2	2	1
Warranty (years)	1	1	1	1	1	1
PRODUCT DIMS (w x d x h)						
Platform Size	12 ¼" x 12 ¼" x 1 ¾"	13" x 13" x ⅞"	17" x 13" x 1 ¼"	12 ⅝" x 12 ⅝" x 1 ⅞"	14 ¼" x 11 ¾" x 1 ¾"	14 ¼" x 11 ¾" x 1 ¾"
Product Footprint	12 ¼" x 12 ¼" x 1 ¾"	13" x 13" x ⅞"	17" x 13" x 1 ¼"	12 ⅝" x 12 ⅝" x 1 ⅞"	14 ¼" x 11 ¾" x 1 ¾"	14 ¼" x 11 ¾" x 1 ¾"
Product Weight	4 lb / 2 kg	4 lb / 2 kg	8 lb / 4 kg	4 lb / 2 kg	4 lb / 2 kg	4 lb / 2 kg
FUNCTIONS						
LB/ KG Conversion		•				•
Auto Zero	•	•	•	•	•	•
Auto Off	•	•	•	•	•	•
Weight Hold		•				•
Bilingual Talking Scale						•
Voice Disable						•
ACCESSORIES (Optional)						
Carrying Case	64771		64769	64771	64771	64771

	142KL	160KL	160LB	160KG
				
FEATURES				
Capacity	330 lb / 150 kg	400 lb / 180 kg	400 lb	180 kg
Resolution	1 lb / 0.5 kg	1 lb / 1 kg	1 lb	1 kg
Dial Size	6 ½" / 165 mm	8" / 203 mm	8" / 203 mm	8" / 203 mm
Case Pack Quantity	2	2	2	2
Warranty (years)	1	1	1	1
PRODUCT DIMS (w x d x h)				
Platform Size	11 ½" x 10 ¾" x 2 ¼"	11" x 12 ½" x 3"	11" x 12 ½" x 3"	11" x 12 ½" x 3"
Product Footprint	11 ¾" x 17 ¾" x 3 ½"	11 ¾" x 18 ¼" x 7 ¼"	11 ¾" x 18 ¼" x 7 ¼"	11 ¾" x 18 ¼" x 7 ¼"
Product Weight	6 lb / 3 kg	9 lb / 4 kg	9 lb / 4 kg	9 lb / 4 kg

	100LB	100KG
		
FEATURES		
Capacity	270 lb	125 kg
Resolution	1 lb	1 lb
Dial Size	2" / 51 mm	2" / 51 mm
Case Pack Quantity	3	3
Warranty (years)	1	1
PRODUCT DIMS (w x d x h)		
Platform Size	10 ¼" x 10" x 1 ⅞"	10 ¼" x 10" x 1 ⅞"
Product Footprint	10 ¼" x 10" x 2"	10 ¼" x 10" x 2"
Product Weight	2 lb / .9 kg	2 lb / .9 kg
ACCESSORIES (Optional)		
Carrying Case	64771	64771

Specialty Scales

	PTS-1000KL	2000KL	3400KL	3401KL	YG500R
					
FEATURES					
Capacity	550 lb / 250 kg	400 lb / 181 kg	176 oz / 5000 g	176 oz / 5000 g	500 g
Resolution/ Graduation	1 lb / 0.5 kg	0.2 lb / 0.1 kg	0.05 oz / 1 g	0.05 oz / 1 g	2 g
Display/ Dial Size	1" / 25 mm	1 1/2" / 38 mm	5/8" / 16 mm	5/8" / 16 mm	6" / 152 mm
Wheels		•			
EMR Capable		•			
Power Adapter Included	•	•	•	•	
Warranty (years)	2	2	2	2	1
Wireless Connectivity		Optional			
PRODUCT DIMS (w x d x h)					
Platform Size	71" x 21"	70 1/2" x 30"	9 1/2" x 11 3/4" x 2 1/2"	7 3/4" x 9 3/4" x 5/8"	6 1/2" x 4 2/3"
Product Footprint	71" x 27 1/2" x 1"	70 1/2" x 46 1/2" x 60 1/2"	10 7/8" x 13" x 8"	10" x 11 3/4" x 5 3/8"	6 1/2" x 6 7/8" x 8 5/8"
Product Weight	25 lb / 11.4 kg	104 lb / 47 kg	11 lb / 5 kg	10 lb / 4.5 kg	3 lb / 1 kg
FUNCTIONS					
LB/ KG Conversion	•	•	•	•	
LB/ KG Lockout		•			
Hold/ Release	•	•			
BMI		•			
Zero		•	•	•	
Tare		•	•	•	
Pretare		•			
Auto Zero		•	•	•	
Auto Off	•	•	•	•	
Reweigh		•			
Recall		•			
Everlock®		•			
Time/ Date		•			
Variable Auto Off Time		•			
Audible/ Mute Sound		•			
Enabling Unit View		•			
ACCESSORIES (Optional)					
Power Adapter	69-00012	ADPT30	3400ADPT	3400ADPT	
Wireless Kit		•			
Wireless Printer		•			
Other			3400PAN (3 pack pans) 3400TRAY (3 pack trays)	3400PAN (3 pack pans) 3400TRAY (3 pack trays)	

Height Measuring Instruments

	201HR	205HR	PORTROD	WMROD
				
FEATURES				
Measurement Range	2° - 90° / 6 - 230 cm	0 - 90° / 0 - 230 cm	24° - 83° / 60 - 213 cm	11 7/8° - 84° / 30 - 213 cm
Graduation	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/8" / 1 mm
Warranty (years)	2	2	2	2
PRODUCT DIMS (w x d x h)				
Product Footprint	2 1/4" x 11" x 54"	2 5/8" x 8" x 60"	2 3/4" x 8 1/4" x 60"	1 1/2" x 13 1/2" x 80 1/4"
Product Weight	2 lb / 1 kg	2 lb / 1 kg	1 lb / .5 kg	4.2 lb / 2 kg

	221HR	280HR	BABYHR	EBABYHR
				
FEATURES				
Measurement Range	7° - 82.6° / 17.7 - 210 cm	3.9 - 43.3° / 10 - 110 cm	3° - 39° / 7 - 100 cm	13 3/4° - 31 1/2° / 35 - 80 cm
Graduation	1/16" / 1 mm	1/4" / 5 mm	1/16" / 1 mm	1/16" / 1 mm
Warranty (years)	2	2	2	2
PRODUCT DIMS (w x d x h)				
Product Footprint	13 3/4" x 23" x 89"	55" x 12"	41" x 11 1/4" x 1 1/2"	24 3/4" x 11 1/2" x 2 3/4"
Product Weight	7.9 lb / 3.6 kg	.88 lb / .4 kg	1 lb / .5 kg	1 lb / .5 kg

Height Measuring Accessories for Platform Scales

	Height Rods for 1100KL Series	Height Rods for 1110KL Series*	Height Rods for 2101KL Series	Height Rods for 3001KL-AM Series	Height Rods for 3105KL-AM Series
	PROPLUSROD	201HR-1110	STROD	3001HR	201HR-3105
					
FEATURES					
Mechanical Measurement Range	2" - 90" / 6 - 230 cm	2" - 90" / 6 - 230 cm	2" - 90" / 6 - 230 cm	2" - 90" / 6 - 230 cm	2" - 90" / 6 - 230 cm
Mechanical Measurement Graduation	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm
Warranty (years)	2	2	2	2	2
PRODUCT DIMS (w x d x h)					
Product Footprint	2 ¼" x ¾" x 54"	2 ¼" x ¾" x 54"	2 ¼" x ¾" x 54"	2 ¼" x ¾" x 54"	2 ¼" x ¾" x 54"
Product Weight	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg
	245EHR-1100	245EHR-1110	245EHR-2101	245EHR-3001	245EHR-3105
					
FEATURES					
Mechanical Measurement Range	11 ¾" - 79 ½" / 30 - 202 cm	11 ¾" - 79 ½" / 30 - 202 cm	11 ¾" - 79 ½" / 30 - 202 cm	11 ¾" - 79 ½" / 30 - 202 cm	11 ¾" - 79 ½" / 30 - 202 cm
Mechanical Measurement Graduation	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm
Digital Measurement Range	3 ft. 11.2" - 6 ft. 7.5" / 120 cm - 202 cm	3 ft. 11.2" - 6 ft. 7.5" / 120 cm - 202 cm	3 ft. 11.2" - 6 ft. 7.5" / 120 cm - 202 cm	3 ft. 11.2" - 6 ft. 7.5" / 120 cm - 202 cm	3 ft. 11.2" - 6 ft. 7.5" / 120 cm - 202 cm
Digital Measurement Resolution	0.1" / 1 mm	0.1" / 1 mm	0.1" / 1 mm	0.1" / 1 mm	0.1" / 1 mm
Warranty (years)	2	2	2	2	2
PRODUCT DIMS (w x d x h)					
Product Footprint	2 ¼" x 2" x 48"	2 ¼" x 2" x 48"	2 ¼" x 2" x 48"	2 ¼" x 2" x 48"	2 ¼" x 2" x 48"
Product Weight	4 lb / 2 kg	4 lb / 2 kg	4 lb / 2 kg	4 lb / 2 kg	4 lb / 2 kg

	Height Rod for 400KL Series	Height Rod for 499KL & 599KL Series
	402KLROD	499KLROD
		
FEATURES		
Mechanical Measurement Range	11 ⅞" - 84" / 30 - 213 cm	11 ⅞" - 84" / 30 - 213 cm
Mechanical Measurement Graduation	1/8" / 1 mm	1/8" / 1 mm
Warranty (years)	2	2
PRODUCT DIMS (w x d x h)		
Product Footprint	1 ½" x 13 ½" x 80 ¼"	1 ½" x 13 ½" x 80 ¼"
Product Weight	4.2 lb / 2 kg	4.2 lb / 2 kg

*Must have an 1110KL and 1110HB handlebars

Height Measuring Accessories for Pediatric Scales

	Height Rods for 522KL Series	Height Rods for 553KL Series	Height Rod for 550KL	Height Rod for 565KL Series	Height Rod for 2210KL-AM Series
	522EHR	553EHR	eBABYHR + 550EHRBRACKET	565EHR	2210EHR
					
FEATURES					
Digital Measurement Range	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm
Digital Measurement Resolution	1⁄16" / 0.1 cm	1⁄16" / 0.1 cm	1⁄16" / 0.1 cm	1⁄16" / 0.1 cm	1⁄16" / 0.1 cm
Warranty (years)	2	2	2	2	2
PRODUCT DIMS (w x d x h)					
Product Footprint	2 3⁄4" x 11 1⁄2" x 24 3⁄4"	2 3⁄4" x 11 1⁄2" x 24 3⁄4"	2 3⁄4" x 11 1⁄2" x 24 3⁄4"	2 3⁄4" x 11 1⁄2" x 24 3⁄4"	2 3⁄4" x 11 1⁄2" x 24 3⁄4"
Product Weight	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg
	522HR	553HR			
					
FEATURES					
Mechanical Measurement Range	3" - 39" / 7 - 100 cm	3" - 39" / 7 - 100 cm			
Mechanical Measurement Graduation	1⁄16" / 1 mm	1⁄16" / 1 mm			
Warranty (years)	2	2			
PRODUCT DIMS (w x d x h)					
Product Footprint	1 1⁄2" x 11 1⁄4" x 41"	1 1⁄2" x 11 1⁄4" x 41"			
Product Weight	1 lb / .5 kg	1 lb / .5 kg			