

Designed for Accessibility. Certified for Compliance.

Health o meter® Professional's ADA compliant scales are thoughtfully engineered to support facilities in meeting accessibility standards. These scales comply with the latest Medical Diagnostic Equipment (MDE) standards and are reviewed for ADA Compliance and Platinum Certified by Corada.


Wheelchair Scales

MDE Requirements:

- 1 32" wide by 40" deep platform
- 2 2" edge protection on platform
- 3 2" edge protection on ramps
- 4 1:8 running ramp slope


2620KL

Digital Wheelchair Scale



2920KL-AM

Antimicrobial Digital Wall-Mounted Fold-Up Wheelchair Scale



Stand-On Scales

MDE Requirements:

- 1 Slip resistant platform
- 2 Handrail supports on both sides of platform (Handrail along length of platform for wheelchair scales)
- 3 34" – 38" max. handrail height from platform

2101KL

Digital Platform Scale



1100KL

Digital Platform Scale with Extra-Wide Handrails



Meets Wheelchair & Stand-On MDE Requirements

2622KL

Digital Wheelchair and Stand-On Scale



Scan to learn more about our certified ADA Compliant scales and the new MDE standards.
or visit www.homscales.com/ada-mde

