Health o meter®
Professional

The Acute Care Line

Engineered with a Purpose®

Designed to meet the demands of acute care environments, Health o meter® Professional Acute Care scales deliver reliable performance, efficient workflows, and safer weighing experiences for patients and caregivers.



Reliable & Durable

Robust construction designed to withstand frequent use, movement, and routine cleaning in demanding hospital environments.

Fast, Efficient Workflows

Quick, accurate weight measurements and intuitive features help streamline daily clinical workflows.

Patient & Caregiver Safety

Designed to support patient stability, help reduce fall risks, and promote safer weighing experiences.

Built for Every Patient

From neonatal to bariatric and wheelchair accessible solutions, the Acute Care scale line supports a wide range of patient needs.



Heavy-Duty Mobile Antimicrobial Platform Scale

3105KL-AM



- Capacity of 1000 lb / 454 kg
- EMR connectivity, BMI, hold / release, reweigh, recall, zero, KG / LB, tare, pre-tare, Everlock®
- Optional electronic and mechanical height rod
- Safety features include antimicrobial protection on the head, keypad, transport handle, handrails, body, and kickplate
- Options such as KG only and built-in wireless models are available

High-Resolution Antimicrobial Digital Neonatal / Pediatric Tray Scale

2210KL-AM (2-sided tray)
2210KL4-AM (4-sided tray)
2210CART-AC (acute care cart)





- High 1 gram resolution with a capacity of 20,000 g / 20 kg / 45 lb
- EMR connectivity, breast milk intake, auto-weight hold, hold/release, zero, tare, pre-tare, reweigh, recall, GM / KG / LB conversion and lockout, auto-zero, and auto-off
- Round-edged removable tray and backlit display
- Optional electronic height rod
- Options such as KG only and built-in wireless models are available

Antimicrobial Wall-Mounted Fold Down Wheelchair Scale

2900KL-AM
2920KL-AM (ADA Compliant)





- Capacity of 1000 lb / 454 kg
- Large easy-to-lift platform conserves space
- EMR connectivity, BMI, hold / release, reweigh, recall, zero, KG / LB, tare, pre-tare, Everlock®
- Safety features include anti-slip caution tape, raised edges, recessed handle, and no pinch points
- Options such as KG only and built-in wireless models are available
- ADA Compliant model available to support facilities in meeting accessibility standards

Digital Wet Diaper / Lap Sponge / Organ Scale

3400KL / 3400KG (Pan)
3401KL (Tray)





- Capacity of 176 oz / 5000 g
- Stable stainless steel base and removable stainless steel weighing pan or tray that can be autoclaved for sterilization
- Powered by a long-life rechargeable battery
- Backlit display that's easy to read from any angle and light condition

Patient Transfer Scale

PTS-1000KL / PTS-1000KG



- Capacity of 550 lb / 250 kg
- A transfer board with a built-in medical scale that provides a weight measurement during patient transfers
- Provides quick, accurate weight measurements for immobile patients
- Helps shorten door-to-needle times for time-critical patients
- The PTS supports workflows across multiple departments, including Stroke Units, Emergency Departments, ICUs, and Radiology
- KG only model is a safe solution for facilities standardizing on metric-only measurements for critical medication dosages based on kilograms

Experience the benefits of our Acute Care Scales with a no-obligation free product trial and live virtual demo!

Request your trial today or join a virtual demo led by our expert Product Specialists to see the features in action and get all your questions answered.

Learn more about our Acute Care scale line or request a free trial: www.homscales.com/acutecarescales

Scan to learn more.